

#### Description

The EHR system receives and incorporates an electronic lab result message. This test plan and the previous test plan are designed to test the EHR-S' capability to process the progression from a "preliminary" to "final" to "corrected" microbiology Culture and Sensitivity reports.

#### PreCondition

Patient information is pre-loaded in the EHR-S.  
No other Pre-Condition.

#### PostCondition

The test message information received by the EHR-S has been incorporated with the patient's record.

#### TestObjectives

- Demonstrate ability to import and incorporate a final and corrected result microbiology message when the filler order number does not uniquely identifies the test ordered. (FRN profiles).