## - Test Case Information-

## EDOS\_2.3\_1.1-M08\_GU - M08\_Update\_reactivate

Test Case ID EDOS\_2.3\_1.1-M08\_GU

-MSH-

		MSU	
Location	Data Element	Data	Categorization
MSH.1	Field Separator		IG Fixed Data
MSH.2	<b>Encoding Characters</b>	^\&#</td><td>IG Fixed Data</td></tr><tr><td>MSH.3</td><td>Sending Application</td><td></td><td></td></tr><tr><td>MSH.3.1</td><td>Namespace ID</td><td>NIST Test Lab APP</td><td>Configurable Data</td></tr><tr><td>MSH.3.2</td><td>Universal ID</td><td>2.16.840.1.113883.3.72.5.20</td><td>Configurable Data</td></tr><tr><td>MSH.3.3</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr><tr><td>MSH.4</td><td>Sending Facility</td><td></td><td></td></tr><tr><td>MSH.4.1</td><td>Namespace ID</td><td>NIST Lab Facility</td><td>Configurable Data</td></tr><tr><td>MSH.4.2</td><td>Universal ID</td><td>2.16.840.1.113883.3.72.5.21</td><td>Configurable Data</td></tr><tr><td>MSH.4.3</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr><tr><td>MSH.6</td><td>Receiving Facility</td><td></td><td></td></tr><tr><td>MSH.6.1</td><td>Namespace ID</td><td>NIST EHR Facility</td><td>Configurable Data</td></tr><tr><td>MSH.6.2</td><td>Universal ID</td><td>2.16.840.1.113883.3.72.5.23</td><td>Configurable Data</td></tr><tr><td>MSH.6.3</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr><tr><td>MSH.7</td><td>Date/Time Of Message</td><td></td><td></td></tr><tr><td>MSH.7.1</td><td>Time</td><td>20130421113601-0700</td><td>System Generated</td></tr><tr><td>MSH.9</td><td>Message Type</td><td></td><td></td></tr><tr><td>MSH.9.1</td><td>Message Code</td><td>MFN</td><td>IG Fixed Data</td></tr><tr><td>MSH.9.2</td><td>Trigger Event</td><td>M08</td><td>IG Fixed Data</td></tr><tr><td>MSH.9.3</td><td>Message Structure</td><td>MFN_M08</td><td>IG Fixed Data</td></tr><tr><td>MSH.10</td><td>Message Control ID</td><td>EDOS_2.3_1.1-M08_GU</td><td>System Generated</td></tr><tr><td>MSH.11</td><td>Processing ID</td><td></td><td></td></tr><tr><td>MSH.11.1</td><td>Processing ID</td><td>D</td><td>Changeable Data</td></tr><tr><td>MSH.12</td><td>VersionID</td><td></td><td></td></tr><tr><td>MSH.12.1</td><td>Version ID</td><td>2.5.1</td><td>IG Fixed Data</td></tr><tr><td>MSH.21</td><td>Message Profile Identifier</td><td></td><td></td></tr><tr><td>MSH.21.1</td><td>Entity Identifier</td><td>EDOS_Common_Component</td><td>Test Case Fixed Data</td></tr><tr><td>MSH.21.2</td><td>Namespace ID</td><td>EDOS Base Profile</td><td>Changeable Data</td></tr><tr><td>MSH.21.3</td><td>Universal ID</td><td>2.16.840.1.113883.9.67</td><td>Test Case Fixed Data</td></tr><tr><td>MSH.21.4</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr><tr><td>MSH.21[2]</td><td>Message Profile Identifier</td><td></td><td></td></tr><tr><td>MSH.21[2].1</td><td>Entity Identifier</td><td>EDOS_GU_Component</td><td>Test Case Fixed Data</td></tr><tr><td>MSH.21[2].2</td><td>Namespace ID</td><td>EDOS GU Profile</td><td>Changeable Data</td></tr><tr><td>MSH.21[2].3</td><td>Universal ID</td><td>2.16.840.1.113883.9.68</td><td>Test Case Fixed Data</td></tr><tr><td>MSH.21[2].4</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr></tbody></table>	

-MFI-

Location	Data Element	Data	Categorization	
MFI.1	Master File Identifier			
MFI.1.1	Identifier	OMM	IG Fixed Data	
MFI.1.2	Text	Mixed type observation master file	Changeable Data	
MFI.1.3	Name of Coding System	HL70175	IG Fixed Data	
MFI.1.7	Coding System Version ID	2.5.1	Test Case Fixed Data	
MFI.3	File-Level Event Code	UPD	Test Case Fixed Data	
MFI.6	Response Level Code	NE	IG Fixed Data	

-MFE-

Location	Data Element	Data	Categorization
MFE.1	Record-Level Event Code	MAC	Test Case Fixed Data

MFE <sub>Bocation</sub>	Effective Datp/AinElement Data		Categorization	
MFE.3.1	Time	20131219145310	System Generated	
MFE.4	Primary Key Value - MFE			
MFE.4.1	Identifier	500	Changeable Data	
MFE.4.2	Text	Erythrocyte sedimentation rate	Changeable Data	
MFE.4.3	Name of Coding System	99USL	Changeable Data	
MFE.4.7	Coding System Version ID	20130421	Changeable Data	
MFE.5	Primary Key Value Type	CWE	IG Fixed Data	

-OM1-

Location	Data Element	Data	Categorization	
OM1.1	Sequence Number - Test/Observation Master File		IG Fixed Data	
OM1.2	Producer's Service/Test/Observation ID			
OM1.2.1	Identifier	500	Changeable Data	
OM1.2.2	Text	Erythrocyte sedimentation rate	Changeable Data	
OM1.2.3	Name of Coding System	99USL	Changeable Data	
OM1.2.7	Coding System Version ID	20130421	Changeable Data	
OM1.4	Specimen Required	Y	Changeable Data	
OM1.5	Producer ID			
OM1.5.1	Identifier	05D0669071	Changeable Data	
OM1.5.2	Text	Century Hospital Clinical Laboratory	Changeable Data	
OM1.5.3	Name of Coding System	99USL	Changeable Data	
OM1.5.7	Coding System Version ID	2013	Changeable Data	
OM1.7	Other Service/Test/Observation IDs for the Observation			
OM1.7.1	Identifier	30341-2	Test Case Fixed Data	
OM1.7.1	Text	Erythrocyte sedimentation rate	Changeable Data	
OM1.7.3	Name of Coding System	LN	Test Case Fixed Data	
OM1.7.7	Coding System Version ID	2.42	Changeable Data	
OM1.7[2]	Other Service/Test/Observation IDs	2.42	Changeaok Data	
	for the Observation  Identifier	416838001		
OM1.7[2].1	Identifier	Erythrocyte sedimentation rate		
OM1.7[2].2	Text	measurement		
OM1.7[2].3	Name of Coding System	SCT		
OM1.7[2].7	Coding System Version ID	201509-US Ed		
OM1.9	Preferred Report Name for the Observation	Erythrocyte sedimentation rate	Changeable Data	
OM1.12	Orde rability	Y	Test Case Fixed Data	
OM1.12 OM1.18	Nature of Service/Test/Observation	A	Test Case Fixed Data	
OM1.32	Interpretation of Observations	The erythrocyte sedimentation rate is a nonspecific measure of inflammatory disease.	Changeable Data	
OM1.39	Factors that may Affect the Observation	Insufficient blood, Clotting, Hemolysis, Blood specimen received > 12 hours after collection.	Changeable Data	
OM1.40	Service/Test/Observation Performance Schedule	Daily	Changeable Data	
OM1.40[2]	Service/Test/Observation Performance Schedule	Continuously	Changeable Data	
OM1.48	<b>Exclusive Test</b>	N	Changeable Data	
OM1.49	Diagnostic Service Sector ID	LAB	Changeable Data	
OM1.51	Other Names	Westergren	Changeable Data	
OM1.51[2]	Other Names	Sed Rate	Changeable Data	
OM1.54	Special Instructions	Please include tentative diagnosis/treatment on the request form.	Changeable Data	
OM1.54[2]	Special Instructions	Please direct any questions regarding this test to the hematology division.	Changeable Data	

OM155cation			Categorization
OM1.55.1	Identifier	CP	Changeable Data
OM1.55.2	Text	Clinical Pathology	Changeable Data
OM1.55.3	Name of Coding System	99USL	Changeable Data
OM1.55.7	Coding System Version ID	2015	Changeable Data
OM1.55[2]	Test Relationship Category		
OM1.55[2].1	Identifier	HM	Changeable Data
OM1.55[2].2	Text	Hematology	Changeable Data
OM1.55[2].3	Name of Coding System	99USL	Changeable Data
OM1.55[2].7	Coding System Version ID	2015	Changeable Data
OM1.56	Observation Identifier associated with Producer's Service/Test/Observation ID		
OM1.56.1	Identifier	30341-2	Changeable Data
OM1.56.2	Text	Erythrocyte sedimentation rate	Changeable Data
OM1.56.3	Name of Coding System	LN	Test Case Fixed Data
OM1.56.4	Alternate Identifier	500	Changeable Data
OM1.56.5	Alternate Text	Erythrocyte sedimentation rate	Changeable Data
OM1.56.6	Name of Alternate Coding System	99USL	Changeable Data
OM1.56.7	Coding System Version ID	2.48	Changeable Data
OM1.56.8	Alternate Coding System Version ID	20130421	Changeable Data
OM1.56.9	Original Text	Erythrocyte sedimentation rate	Changeable Data
OM1.57	<b>Expected Turn-Around Time</b>		
OM1.57.1	Quantity	1	Changeable Data
OM1.57.2	Units		
OM1.57.2.1	Identifier	d	Changeable Data
OM1.57.2.2	Text	day	Changeable Data
OM1.57.2.3	Name of Coding System	UCUM	Changeable Data
OM1.57.2.7	Coding System Version ID	1.8	Changeable Data

-OM2-

Location	Data Element	Data	Categorization
OM2.1	Sequence Number - Test/Observation Master File		
OM2.2	Units of Measure		
OM2.2.1	Identifier	mm/h	Changeable Data
OM2.2.2	Text	millimeter per hour	Changeable Data
OM2.2.3	Name of Coding System	UCUM	Changeable Data
OM2.2.4	Alternate Identifier	mm/hour	Changeable Data
OM2.2.5	Alternate Text	mm/hour	Changeable Data
OM2.2.6	Name of Alternate Coding System	99USL	Changeable Data
OM2.2.7	Coding System Version ID	1.8	Changeable Data
OM2.2.8	Alternate Coding System Version ID	2014	Changeable Data
OM2.2.9	Original Text	mm/hour	Changeable Data
OM2.6	Reference (Normal) Range for Ordinal and Continuous Observations		
OM2.6.1	Numeric Range		
OM2.6.1.1	Low Value	0	Changeable Data
OM2.6.1.2	High Value	15	Changeable Data
OM2.6.2	Administrative Sex		
OM2.6.2.1	Identifier	M	Changeable Data
OM2.6.2.2	Text	Male	Changeable Data
OM2.6.2.3	Name of Coding System	HL70001	IG Fixed Data
OM2.6.2.7	Coding System Version ID	2.5.1	IG Fixed Data
OM2.6[2]	Reference (Normal) Range for Ordinal and Continuous Observations		
OM2.6[2].1	Numeric Range		
OM2.6[2].1.1	Low Value	0	Changeable Data

OMacalion.2	High Valu <b>Data Element</b>	25	Data	Changealtego attation
OM2.6[2].2	Administrative Sex			
OM2.6[2].2.1	Identifier	F		Changeable Data
OM2.6[2].2.2	Text	Female		Changeable Data
OM2.6[2].2.3	Name of Coding System	HL70001		IG Fixed Data

-OM4-

Location	Data Element	Data	Categorization	
OM4.1	Sequence Number - Test/Observation Master File	1.1	IG Fixed Data	
OM4.3	Container Description	Black Top Tube (Vac-Tec)	Changeable Data	
OM4.4	Container Volume	3.0	Changeable Data	
OM4.5	Container Units			
OM4.5.1	Identifier	mL	Changeable Data	
OM4.5.2	Text	milliliters	Changeable Data	
OM4.5.3	Name of Coding System	UCUM	Changeable Data	
OM4.5.7	Coding System Version ID	1.8	Changeable Data	
OM4.6	Specimen			
OM4.6.1	Identifier	119297000	Test Case Fixed Data	
OM4.6.2	Text	Blood sample	Changeable Data	
OM4.6.3	Name of Coding System	SCT	Test Case Fixed Data	
OM4.6.4	Alternate Identifier	WBLD	Changeable Data	
OM4.6.5	Alternate Text	Whole blood	Changeable Data	
OM4.6.6	Name of Alternate Coding System	99USL	Changeable Data	
OM4.6.7	Coding System Version ID	201509-US Ed	Changeable Data	
OM4.6.8	Alternate Coding System Version ID	2014	Changeable Data	
OM4.6.9	Original Text	Whole blood	Changeable Data	
OM4.7	Additive			
OM4.7.1	Identifier	WEST	Changeable Data	
OM4.7.2	Text	Buffered Citrate (Westergren Sedimentation Rate)	Changeable Data	
OM4.7.3	Name of Coding System	HL70371	Changeable Data	
OM4.7.7	Coding System Version ID	2.5.1	Changeable Data	
OM4.10	Normal Collection Volume			
OM4.10.1	Quantity	2.4	Changeable Data	
OM4.10.2	Units			
OM4.10.2.1	Identifier	mL	Changeable Data	
OM4.10.2.2	Text	milliliters	Changeable Data	
OM4.10.2.3	Name of Coding System	UCUM	IG Fixed Data	
OM4.10.2.7	Coding System Version ID	1.8	Changeable Data	
OM4.11	Minimum Collection Volume			
OM4.11.1	Quantity	2.4	Changeable Data	
OM4.11.2	Units			
OM4.11.2.1	Identifier	mL	Changeable Data	
OM4.11.2.2	Text	milliliters	Changeable Data	
OM4.11.2.3	Name of Coding System	UCUM	IG Fixed Data	
OM4.11.2.7	Coding System Version ID	1.8	Changeable Data	
OM4.12	Specimen Requirements	Specimen must be received and tested within 4-6 hrs of collection	Changeable Data	
OM4.15	Specimen Handling Code			
OM4.15.1	Identifier	CREF	Changeable Data	
OM4.15.2	Text	Critical refrigerated	Changeable Data	
OM4.15.3	Name of Coding System	HL70376	IG Fixed Data	
OM4.15.7	Coding System Version ID	2.5.1	Changeable Data	
OM4.15[2]	Specimen Handling Code		-	
OM4.15[2].1	Identifier	MTLF	Changeable Data	
OM4.15[2].2	Text	Metal Free	Changeable Data	

OMHolestion	Name of Coping Spreement	HL70376	Data	IG Fixed Regorization
OM4.15[2].7	Coding System Version ID	2.5.1		Changeable Data
OM4.16	Specimen Preference	P		Test Case Fixed Data

-OM4-

Location	Data Element	Data	Categorization	
OM4.1	4.1 Sequence Number - Test/Observation Master File 1.2		IG Fixed Data	
OM4.3	Container Description	Lavender Top (EDTA) tube	Changeable Data	
OM4.4	Container Volume	3.0	Changeable Data	
OM4.5	Container Units			
OM4.5.1	Identifier	mL	Changeable Data	
OM4.5.2	Text	milliliters	Changeable Data	
OM4.5.3	Name of Coding System	UCUM	Changeable Data	
OM4.5.7	Coding System Version ID	1.8	Changeable Data	
OM4.6	Specimen			
OM4.6.1	Identifier	119297000	Test Case Fixed Data	
OM4.6.2	Text	Blood sample	Changeable Data	
OM4.6.3	Name of Coding System	SCT	Test Case Fixed Data	
OM4.6.4	Alternate Identifier	WBLD	Changeable Data	
OM4.6.5	Alternate Text	Whole blood	Changeable Data	
OM4.6.6	Name of Alternate Coding System	99USL	Changeable Data	
OM4.6.7	Coding System Version ID	201509-US Ed	Changeable Data	
OM4.6.8	Alternate Coding System Version ID	2014	Changeable Data	
OM4.6.9	Original Text	Whole blood	Changeable Data	
OM4.7	Additive			
OM4.7.1	Identifier	EDTK	Changeable Data	
OM4.7.2	Text	Potassium/K EDTA	Changeable Data	
OM4.7.3	Name of Coding System	HL70371	Changeable Data	
OM4.7.7	Coding System Version ID	2.5.1	Changeable Data	
OM4.10	Normal Collection Volume			
OM4.10.1	Quantity	2.4	Changeable Data	
OM4.10.2	Units			
OM4.10.2.1	Identifier	mL	Changeable Data	
OM4.10,2.2	Text	milliliters	Changeable Data	
OM4.10.2.3	Name of Coding System	UCUM	IG Fixed Data	
OM4.10.2.7	Coding System Version ID	1.8	Changeable Data	
OM4.11	Minimum Collection Volume		Classification   Clas	
OM4.11.1	Quantity	2.4	Changeable Data	
OM4.11.2	Units		Commission of the commission	
OM4.11.2.1	Identifier	mL	Changeable Data	
OM4.11.2.2	Text	milliliters	Changeable Data	
OM4.11.2.3	Name of Coding System	UCUM	IG Fixed Data	
OM4.11.2.7	Coding System Version ID	1.8	Changeable Data	
OM4.12	Specimen Requirements	Specimen must be received and tested within 4-6 hrs of collection.	Changeable Data	
OM4.15	Specimen Handling Code			
OM4.15.1	Identifier	CREF	Changeable Data	
OM4.15.2	Text	Critical refrigerated	Changeable Data	
OM4.15.3	Name of Coding System	HL70376	IG Fixed Data	
OM4.15.7	Coding System Version ID	2.5.1	Changeable Data	
OM4.15[2]	Specimen Handling Code		J	
OM4.15[2].1	Identifier	MTLF	Changeable Data	
OM4.15[2].2	Text	Metal Free	Changeable Data	
OM4.15[2].3	Name of Coding System	HL70376	IG Fixed Data	
OM4.15[2].7	Coding System Version ID	2.5.1	Changeable Data	
	Specimen Preference	A	Test Case Fixed Data	

OM4 <u>Ll</u> dcation	Preferred Sprzingen Attribute	1.1	Data	Test Coategordzation
	Sequence			