| Element | Data |
|------------|------|
| Identifier | OMC |
| Event code | UPD |

Panel/Test battery Information - CMP Record information

| Element | Data |
|-------------------------------------|----------------|
| Event code | MAD |
| Effective date and time | 20131219145310 |
| Test/Panel Name | CMP |
| Test/Panel identifier (code system) | 100 (99USL) |

General information-

| Element | Data |
|---|--|
| Specimen required ? | Y |
| Test/Panel other name | Comprehensive metabolic 2000 panel - Serum or Plasma |
| Test/Panel other identifier (code system) | 24323-8 (LN) |
| Preferred Short Name | CMP |
| Is the test/panel orderable ? | Y |
| Nature | P |
| Interpretation | Test used to measure blood sugar, electrolytes and fluid balance, kidney and liver function. |
| Patient preparation | Patient fasting required for 12 hours. |
| Factors that may affect the test | Insufficient specimen, Gross hemolysis, Improper labeling. |
| Schedule | Daily |
| Is the test exclusive? | N |
| Diagnostic Service Sector | LAB |
| Test alias | CMP |
| Expected turn around time | 1 day |

- Panel/Test battery Details-

| Element | Data |
|--|---------------------------|
| Test included in panel/test battery | Serum Glucose |
| Test included in panel/test battery (code & code system) | 104 (99USL) |
| Test included in panel/test battery | Blood Urea Nitrogen (BUN) |
| Test included in panel/test battery (code & code system) | 106 (99USL) |
| Test included in panel/test battery | Creatinine |
| Test included in panel/test battery (code & code system) | 102 (99USL) |
| Test included in panel/test battery | BUN/Creatinine Ratio |
| Test included in panel/test battery (code & code system) | 108 (99USL) |
| Test included in panel/test battery | GFR, calculated |
| Test included in panel/test battery (code & code system) | 110 (99USL) |
| Test included in panel/test battery | Calcium |
| Test included in panel/test battery (code & code system) | 112 (99USL) |
| Test included in panel/test battery | Total protein, serum |
| Test included in panel/test battery (code & code system) | 114 (99USL) |
| Test included in panel/test battery | Albumin |
| Test included in panel/test battery (code & code system) | 116 (99USL) |
| Test included in panel/test battery | Globulin |
| Test included in panel/test battery (code & code system) | 118 (99USL) |
| Test included in panel/test battery | Albumin/globulin ratio |
| Test included in panel/test battery (code & code system) | 120 (99USL) |
| Test included in panel/test battery | Total bilirubin, serum |
| Test included in panel/test battery (code & code system) | 122 (99USL) |

| Test included in panel/test battery | Alkaline phosphatase (ALP) |
|--|----------------------------------|
| Test included in panel/test battery (code & code system) | 124 (99USL) |
| Test included in panel/test battery | Alanine aminotransferase (ALT) |
| Test included in panel/test battery (code & code system) | 126 (99USL) |
| Test included in panel/test battery | Aspartate aminotransferase (ASP) |
| Test included in panel/test battery (code & code system) | 128 (99USL) |
| Test included in panel/test battery | Sodium, serum |
| Test included in panel/test battery (code & code system) | 130 (99USL) |
| Test included in panel/test battery | Potassium, serum |
| Test included in panel/test battery (code & code system) | 132 (99USL) |
| Test included in panel/test battery | Chloride, serum |
| Test included in panel/test battery (code & code system) | 134 (99USL) |
| Test included in panel/test battery | Carbon dioxide, serum |
| Test included in panel/test battery (code & code system) | 136 (99USL) |
| Test included in panel/test battery | Anion gap |
| Test included in panel/test battery (code & code system) | 138 (99USL) |

- Specimen information details -

| Specified information details | |
|--|---|
| Element | Data |
| Sequence ID | 1.1 |
| Container description | Gold Serum Separator tube |
| Container description | Red, No Additive tube |
| Container volume | 5.0 |
| Container volume | 5.0 |
| Container units | milliliter |
| Container units (code & code system) | mL (UCUM) |
| Container units | milliliter |
| Container units (code & code system) | mL (UCUM) |
| Specimen (code & code system) | 119364003 (SCT) |
| Normal collection volume | 1 |
| Normal collection volume unit | milliliter |
| Normal collection volume unit (code & code system) | mL (UCUM) |
| Normal collection volume | 0.5 |
| Normal collection volume unit | milliliter |
| Normal collection volume unit (code & code system) | mL (UCUM) |
| Specimen Requirements | Protect from light. Allow serum tube to clot completely at room temperature. Separate serum or plasma from cells ASAP or within 30 minutes of collection. |
| Specimen Handling Code | Refrigerated temperature |
| Specimen Handling Code (code & code system) | REF (HL70376) |

- Specimen information details-

| Element | Data |
|--|----------------------------|
| Sequence ID | 1.2 |
| Container description | Green Lithium Heparin tube |
| Container volume | 3.0 |
| Container units | milliliter |
| Container units (code & code system) | mL (UCUM) |
| Specimen (code & code system) | 119361006 (SCT) |
| Specimen (code & code system) | HEPL (HL70371) |
| Normal collection volume | 2 |
| Normal collection volume unit | milliliter |
| Normal collection volume unit (code & code system) | mL (UCUM) |

| 0.7 |
|--|
| milliliter |
| mL (UCUM) |
| Centrifuge at 1100-2000 g for a minimum of 10 minutes. |
| Refrigerated temperature |
| REF (HL70376) |
| |