-Description

Based on LRI 1.0 Sed Rate and LRI 1.1 Sed Rate:

Mr. William A. Jones is a 51 year old white male who presented with pain, swelling, and redness in several joints including his elbows, wrists, knees, and hips. In order to rule out rheumatoid arthritis, Dr. Nicholas Radon ordered an erythrocyte sedimentation rate (ESR or Sed Rate) blood test to be performed, specifying that the laboratory was to draw the blood specimen from the patient the same day. The final result of this test was 10 mm/hr, a normal value with the normal reference range being 0-17 mm/hr. The final result report for this lab test was generated by the LIS and transmitted to the patient's record in the ambulatory EHR used in Dr. Radon's office practice later that same day.

The ESR test instrument was recalibrated and the ESR specimen was re-analyzed. After the recalibration and re-analysis, the result value for Mr. Jones' ESR test changed from 10 mm/hr to 20 mm/hr. A Corrected result report with a result value of 20 mm/hr and an abnormal flag of "H" was generated by the LIS and transmitted to the patient's record in the ambulatory EHR used in Dr. Radon's office practice.

Comments

This test case is evaluating the correct use of the status codes in OBR-25 and OBX-11.

-PreCondition-

LRI_1.0_1.1-NG/LRI_1.0_1.1-NG have been previously sent which includes specific content in which the data type for OBX.2 is NM, the result is "10" (OBX.5), the abnormal flag is "N" (OBX.8), and the observation results status is "F" (OBX.11).

PostCondition-

No Post Condition.

TestObjectives

- Demonstrate the ability to message corrected results based on previous results.
- Specific test points include the use of the NM data type for OBX.2, the value of "20" for results (OBX.5), the value of "H" for abnormal flags (OBX.8), and the observation results status is "C" (OBX.11).

Notes to Testers

Run LRI_1.0_1.1-NG_Sed_Rate_Final test case prior to this test case. Note that this test case is explicitly requiring the use of the NM data type for OBX-2 (even though the SN data type is equally correct).

Note that data element SPM-24 is provided in the test data. This element is optional for the sender and required-but maybe empty for the receiver. This element is not subject to conformance testing for this test case (sender).

LRI_1.0_2.1-NG - Maximally Populated SED Rate message - Corrected Results

Test Case ID LRI_1.0_2.1-NG

-MSH-

Location	Data Element	Data	Categorization
MSH.1	Field Separator		IG Fixed Data
MSH.2	Encoding Characters	~\&#</td><td>IG Fixed Data</td></tr><tr><td>MSH.3</td><td>Sending Application</td><td></td><td></td></tr><tr><td>MSH.3.1</td><td>Namespace ID</td><td>NIST Test Lab APP</td><td>Configurable Data</td></tr><tr><td>MSH.3.2</td><td>Universal ID</td><td>372520</td><td>Configurable Data</td></tr><tr><td>MSH.3.3</td><td>Universal ID Type</td><td>L</td><td>Configurable Data</td></tr><tr><td>MSH.4</td><td>Sending Facility</td><td></td><td></td></tr><tr><td>MSH.4.1</td><td>Namespace ID</td><td>NIST Lab Facility</td><td>Configurable Data</td></tr><tr><td>MSH.4.2</td><td>Universal ID</td><td>372521</td><td>Configurable Data</td></tr><tr><td>MSH.4.3</td><td>Universal ID Type</td><td>L</td><td>Configurable Data</td></tr><tr><td>MSH.6</td><td>Receiving Facility</td><td></td><td></td></tr><tr><td>MSH.6.1</td><td>Namespace ID</td><td>NIST EHR Facility</td><td>Configurable Data</td></tr><tr><td>MSH.6.2</td><td>Universal ID</td><td>372523</td><td>Configurable Data</td></tr><tr><td>MSH.6.3</td><td>Universal ID Type</td><td>L</td><td>Configurable Data</td></tr><tr><td>MSH.7</td><td>Date/Time Of Message</td><td></td><td></td></tr><tr><td>MSH.7.1</td><td>Time</td><td>20150926140551</td><td>System Generated</td></tr><tr><td>MSH.9</td><td>Message Type</td><td></td><td></td></tr><tr><td>MSH.9.1</td><td>Message Code</td><td>ORU</td><td>IG Fixed Data</td></tr><tr><td>MSH.9.2</td><td>Event Type</td><td>R01</td><td>IG Fixed Data</td></tr><tr><td>MSH.9.3</td><td>Message Structure</td><td>ORU_R01</td><td>IG Fixed Data</td></tr><tr><td>MSH.10</td><td>Message Control ID</td><td>LRI_1.0_2.1-NG</td><td>System Generated</td></tr><tr><td>MSH.11</td><td>Processing ID</td><td></td><td></td></tr><tr><td>MSH.11.1</td><td>Processing ID</td><td>D</td><td>Changeable Data</td></tr><tr><td>MSH.12</td><td>VersionID</td><td></td><td></td></tr><tr><td>MSH.12.1</td><td>Version ID</td><td>2.5.1</td><td>IG Fixed Data</td></tr><tr><td>MSH.15</td><td>Accept Acknowledgment Type</td><td>AL</td><td>IG Fixed Data</td></tr><tr><td>MSH.16</td><td>Application Acknowledgment Type</td><td>AL</td><td>IG Fixed Data</td></tr><tr><td>MSH.21</td><td>Message Profile Identifier</td><td></td><td></td></tr><tr><td>MSH.21.1</td><td>Entity Identifier</td><td>LRI Common Component</td><td>IG Fixed Data</td></tr><tr><td>MSH.21.2</td><td>Namespace ID</td><td>LRI Base Profile</td><td>Changeable Data</td></tr><tr><td>MSH.21.3</td><td>Universal ID</td><td>2.16.840.1.113883.9.16</td><td>IG Fixed Data</td></tr><tr><td>MSH.21.4</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr><tr><td>MSH.21[2]</td><td>Message Profile Identifier</td><td></td><td></td></tr><tr><td>MSH.21[2].1</td><td>Entity Identifier</td><td>LRI_NG_Component</td><td>IG Fixed Data</td></tr><tr><td>MSH.21[2].2</td><td>Namespace ID</td><td>LRI Base Profile</td><td>Changeable Data</td></tr><tr><td>MSH.21[2].3</td><td>Universal ID</td><td>2.16.840.1.113883.9.13</td><td>IG Fixed Data</td></tr><tr><td>MSH.21[2].4</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr><tr><td>MSH.21[3]</td><td>Message Profile Identifier</td><td></td><td></td></tr><tr><td>MSH.21[3].1</td><td>Entity Identifier</td><td>LRI_FRU_Component</td><td>IG Fixed Data</td></tr><tr><td>MSH.21[3].2</td><td>Namespace ID</td><td>LRI NG Profile</td><td>Changeable Data</td></tr><tr><td>MSH.21[3].3</td><td>Universal ID</td><td>2.16.840.1.113883.9.83</td><td>IG Fixed Data</td></tr><tr><td>MSH.21[3].4</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr></tbody></table>	

-PID-

Location	Data Element	Data	Categorization
PID.1	Set ID - PID	1	IG Fixed Data
PID.3	Patient Identifier List		
PID.3.1	ID Number	PATID1234	Configurable Data
PID.3.4	Assigning Authority		

P Doğatibn	Namespace Element	NIST MPI Data	Changeallegorization	
PID.3.5	Identifier Type Code	MR	Configurable Data	
PID.5	Patient Name			
PID.5.1	Family Name			
PID.5.1.1	Surname	Jones	Changeable Data	
PID.5.2	Given Name	William	Changeable Data	
PID.5.3	Second and Further Given Names or Initials Thereof	A	Changeable Data	
PID.5.7	Name Type Code	L	Changeable Data	
PID.7	Date/Time of Birth			
PID.7.1	Time	19610615	Changeable Data	
PID.8	Administrative Sex	M	Changeable Data	
PID.10	Race			
PID.10.1	Identifier	2106-3	Changeable Data	
PID.10.2	Text	White	Changeable Data	
PID.10.3	Name of the Coding System	HL70005	Changeable Data	
PID.10.9	Original Text	White	Changeable Data	
PID.10[2]	Race			
PID.10[2].1	Identifier	1002-5	Changeable Data	
PID.10[2].2	Text	American Indian or Alaska Native	Changeable Data	
PID.10[2].3	Name of the Coding System	HL70005	Changeable Data	
PID.10[2].9	Original Text	American Indian	Changeable Data	
PID.18	Patient Account Number			
PID.18.1	ID Number	PATID1234	Configurable Data	
PID.18.4	Assigning Authority			
PID.18.4.1	Namespace ID	NIST MPI	Changeable Data	
PID.18.5	Identifier Type Code	AN	Configurable Data	

-ORC-

- ORC -				
Location	Data Element	Data	Categorization	
ORC.1	Order Control	RE	Test Case Fixed Data	
ORC.2	Placer Order Number			
ORC.2.1	Entity Identifier	ORD723222	Changeable Data	
ORC.2.2	Namespace ID	NIST EHR	Changeable Data	
ORC.3	Filler Order Number			
ORC.3.1	Entity Identifier	R-783274	Changeable Data	
ORC.3.2	Namespace ID	NIST Lab Filler		
ORC.4	Placer Group Number			
ORC.4.1	Entity Identifier	GORD874211	Changeable Data	
ORC.4.2	Namespace ID	NIST EHR		
ORC.12	Ordering Provider			
ORC.12.1	ID Number	5742200012	Changeable Data	
ORC.12.2	Family Name			
ORC.12.2.1	Surname	Radon	Changeable Data	
ORC.12.3	Given Name	Nicholas	Changeable Data	
ORC.12.4	Second and Further Given Names or Initials Thereof	M	Changeable Data	
ORC.12.5	Suffix (e.g., JR or III)	JR	Changeable Data	
ORC.12.6	Prefix (e.g., DR)	DR	Changeable Data	
ORC.12.9	Assigning Authority			
ORC.12.9.1	Namespace ID	NPI	Changeable Data	
ORC.12.10	Name Type Code	L	Changeable Data	
ORC.12.13	Identifier Type Code	NPI Changeable Data		

-OBR-

Location	Data Element	Data	Categorization
OBR.1	Set ID - OBR	1	IG Fixed Data

OBR _{Location}	Placer Orden National Placer Orden	Data	Categorization	
OBR.2.1	Entity Identifier	ORD723222	Changeable Data	
OBR.2.2	Namespace ID	NIST EHR	Changeable Data	
OBR.3	Filler Order Number			
OBR.3.1	Entity Identifier	R-783274	Changeable Data	
OBR.3.2	Namespace ID	NIST Lab Filler	Changeable Data	
OBR.4	Universal Service Identifier			
OBR.4.1	Identifier	30341-2	Test Case Fixed Data	
OBR.4.2	Text	Erythrocyte sedimentation rate	Test Case Fixed Data	
OBR.4.3	Name of Coding System	LN	Test Case Fixed Data	
OBR.4.4	Alternate Identifier	500	Changeable Data	
OBR.4.5	Alternate Text	Erythrocyte sedimentation rate	Changeable Data	
OBR.4.6	Name of Alternate Coding System	99USL	Changeable Data	
OBR.4.7	Coding System Version	2.52	Changeable Data	
OBR.4.9	Original Text	Erythrocyte sedimentation rate	Changeable Data	
OBR.7	Observation Date/Time			
OBR.7.1	Time	201509251400	Changeable Data	
OBR.11	Specimen Action Code	L	Changeable Data	
OBR.16	Ordering Provider			
OBR.16.1	ID Number	5742200012	Changeable Data	
OBR.16.2	Family Name			
OBR.16.2.1	Surname	Radon	Changeable Data	
OBR.16.3	Given Name	Nicholas	Changeable Data	
OBR.16.4	Second and Further Given Names or Initials Thereof	M	Changeable Data	
OBR.16.5	Suffix (e.g., JR or III)	JR	Changeable Data	
OBR.16.6	Prefix (e.g., DR)	DR	Changeable Data	
OBR.16.9	Assigning Authority			
OBR.16.9.1	Namespace ID	NPI	Changeable Data	
OBR.16.10	Name Type Code	L	Changeable Data	
OBR.16.13	Identifier Type Code	NPI	Changeable Data	
OBR.22	Results Rpt/Status Chng - Date/Time			
OBR.22.1	Time	20150926140551	Changeable Data	
OBR.25	Result Status	С	Test Case Fixed Data	
OBR.28	Result Copies To			
OBR.28.1	ID Number	10092000194	Changeable Data	
OBR.28.2	Family Name			
OBR.28.2.1	Surname	Hamlin	Changeable Data	
OBR.28.3	Given Name	Pafford	Changeable Data	
OBR.28.4	Second and Further Given Names or Initials Thereof	M	Changeable Data	
OBR.28.5	Suffix (e.g., JR or III)	Sr.	Changeable Data	
OBR.28.6	Prefix (e.g., DR)	Dr.	Changeable Data	
OBR.28.9	Assigning Authority		<i>G</i>	
OBR.28.9.1	Namespace ID	NPI	Changeable Data	
OBR.28.10	Name Type Code	L	Changeable Data	
OBR.28.13	Identifier Type Code	NPI	Changeable Data	
OBR.28.13	Results Handling	114.1	Changeaoic Data	
OBR.49.1	Identifier	CC	Test Case Fixed Data	
ODK.47.1	тиспинст	CC Test Case Fixed Date		
OBR.49.2	Text	Copies Requested Test Case Fixed Data HL70507 IG Fixed Data		

-NTE-

Location	Data Element	Data	Categorization
NTE.1	Set ID - NTE	1	IG Fixed Data
NTE.3	Comment	Patient is extremely anxious about needles used for drawing blood.\.br\If patient is overly frightened, nervous, or anxious	Changeable Data

Location	Data Element	please reschedule beggdadraw.	Categorization	
		NTE-		
Location	Data Element	Data	Categorization	
NTE.1	Set ID - NTE	2	IG Fixed Data	
NTE.3	Comment	Patient is allergic to latex	Changeable Data	
		TQ1		
Location	Data Element	Data	Categorization	
TQ1.1	Set ID - TQ1	1	IG Fixed Data	
TQ1.7	Start date/time			
TQ1.7.1	Time	201509251400	Changeable Data	
TQ1.8	End date/time			
TQ1.8.1	Time	201509261400	Changeable Data	
TQ1.9	Priority			
TQ1.9.1	Identifier	R	Changeable Data	
TQ1.9.2	Text	Routine	Changeable Data	
TQ1.9.3	Name of the Coding System	HL70485	Changeable Data	
TQ1.9.9	Original Text	Routine	Changeable Data	
		—OBX—		
Location	Data Element	Data	Categorization	
OBX.1	Set ID - OBX	1	IG Fixed Data	
OBX.2	Value Type	NM	Test Case Fixed Data	
OBX.3	Observation Identifier			
OBX.3.1	Identifier	30341-2	Test Case Fixed Data	
OBX.3.2	Text	Erythrocyte sedimentation rate	Test Case Fixed Data	
OBX.3.3	Name of the Coding System	LN	Test Case Fixed Data	
OBX.3.4	Alternate Identifier	815117	Changeable Data	
OBX.3.5	Alternate Text	ESR	Changeable Data	
OBX.3.6	Name of Alternate Coding System	99USL	Changeable Data	
OBX.3.7	Coding System Version	2.52	Changeable Data	
OBX.3.9	Original Text	Erythrocyte sedimentation rate	Changeable Data	
OBX.4	Observation Sub-ID			
OBX.4.2	Group	1	Test Case Fixed Data	
OBX.4.3	Sequence	1	Test Case Fixed Data	
OBX.4.4	Identifier	1	Test Case Fixed Data	
OBX.5	Observation Value	20	Test Case Fixed Data	
OBX.6	Units			
OBX.6.1	Identifier	mm/h	Changeable Data	
		mm/h Changeable Data millimeter per hour Changeable Data		
OBX.6.2	Text	millimeter per nour	Changeable Data	

UCUM

0 to 17

201509251400

20150926130550

Century Hospital

CLIA

24D9871327

XX

1.9

Η

С

Changeable Data

Test Case Fixed Data

Test Case Fixed Data

OBX.6.3

OBX.6.7

OBX.7

OBX.8

OBX.11

OBX.14

OBX.19

OBX.23

OBX.14.1

OBX.19.1

OBX.23.1

OBX.23.6

OBX.23.7

OBX.23.10

OBX.24

OBX.23.6.1

Name of the Coding System

Coding System Version

Observation Result Status

Date/Time of the Analysis

Organization Name

Assigning Authority

Namespace ID

Identifier Type Code

Organization Identifier

Performing Organization Address

Date/Time of the Observation

Performing Organization Name

Reference Range

Abnormal Flags

Time

OB X ocation	Street Addr Data Element	Data	Categorization
OBX.24.1.1	Street or Mailing Address	2070 Test Park	Changeable Data
OBX.24.3	City	Los Angeles	Changeable Data
OBX.24.4	State or Province	CA	Changeable Data
OBX.24.5	Zip or Postal Code	90067	Changeable Data
OBX.24.6	Country	USA	Changeable Data
OBX.24.7	Address Type	В	Changeable Data
OBX.24.9	County/Parish Code	06037	Changeable Data
OBX.25	Performing Organization Medical Director		
OBX.25.1	ID Number	5432178916	Changeable Data
OBX.25.2	Family Name		
OBX.25.2.1	Surname	Knowsalot	Changeable Data
OBX.25.3	Given Name	Phil	Changeable Data
OBX.25.4	Second and Further Given Names or Initials Thereof	J.	Changeable Data
OBX.25.5	Suffix (e.g., JR or III)	III	Changeable Data
OBX.25.6	Prefix (e.g., DR)	Dr.	Changeable Data
OBX.25.9	Assigning Authority		
OBX.25.9.1	Namespace ID	NPI	Changeable Data
OBX.25.10	Name Type Code	L	Changeable Data
OBX.25.13	Identifier Type Code	NPI	Changeable Data
OBX.29	Observation Type	RSLT	Test Case Fixed Data
OBX.30	Observation SubType	UNSP	Test Case Fixed Data

-NTE-

Location	Data Element	Data	Categorization
NTE.1	Set ID - NTE	1	IG Fixed Data
NTE.3	('ammant	Specimen re-analyzed per request of ordering provider.	Changeable Data

-SPM-

Location	Data Element	Data	Categorization	
SPM.1	Set ID - SPM	1	IG Fixed Data	
SPM.2	Specimen ID			
SPM.2.1				
SPM.2.1.1		S-2015-66	Configurable Data	
SPM.2.1.2		GoodHealthC_EHR	Configurable Data	
SPM.2.2				
SPM.2.2.1		S-9911-33	Changeable Data	
SPM.2.2.2		NIST Lab Filler	Changeable Data	
SPM.4	Specimen Type			
SPM.4.1	Identifier	119297000	Changeable Data	
SPM.4.2	Text	BLD	Changeable Data	
SPM.4.3	Name of the Coding System	SCT	Changeable Data	
SPM.4.4	Alternate Identifier	BldSpc	Changeable Data	
SPM.4.5	Alternate Text	Blood	Changeable Data	
SPM.4.6	Name of Alternate Coding System	99USL	Changeable Data	
SPM.4.7	Coding System Version	201509USEd	Changeable Data	
SPM.4.9	Original Text	Blood Specimen	Changeable Data	
SPM.17	Specimen Collection Date/Time			
SPM.17.1	Range Start Date/Time			
SPM.17.1.1	Time	201509251400	Changeable Data	
SPM.24	Specimen Condition			
SPM.24.1	Identifier	COOL	Changeable Data	
SPM.24.2	Text	Cool	Changeable Data	
SPM.24.3	Name of the Coding System	HL70493	Changeable Data	
SPM.24.4	Alternate Identifier	CL Changeable Dat		

SPMocation	Alternate T Data Element	Cool	Data	Chang Caltle dorization
SPM.24.6	Name of Alternate Coding System	99USL		Changeable Data
SPM.24.9	Original Text	Cool		Changeable Data

Patient Information			
Data			
William A Jones			
06/15/1961			
Male			
White			
American Indian or Alaska Native			

Order Observation

Ordering Provider

Element	Data
Name	DR Nicholas M Radon JR
Identifier number	5742200012

Observation Details

Element	Data			
Observation General Information				
Placer Order Number	ORD723222			
Filler Order Number	R-783274			
Placer Group Number	GORD874211			
Parent Universal Service Identifier				
Identifier				
Text				
Alt Identifier				
Alt Text				
Original Text				
Observation Details				
Universal Service Identifier	Erythrocyte sedimentation rate			
Observation Date/Time	09/25/2015 2:00 PM			
Observation end Date/Time	07/25/2015 2.00 1 141			
Specimen Action Code	L			
Relevant Clinical Information				
Relevant Clinical Information Original Text				
Observation Result Information				
Result Status	C			
Results Report/Status Change - Date/Time	09/26/2015 2:05 PM			
Results Copy To				
Name	Dr. Pafford M Sr.			
Identifier	10092000194			
Results Handling				
Standard				
Observation Notes				

Notes and comments	Patient is extremely anxious about needles used for drawing blood.\br\If patient is overly frightened, nervous, or anxious please reschedule blood draw.
Notes and comments	Patient is allergic to latex

Element	Data
Priority	Routine
Start Date/time	09/25/2015 2:00 PM
End Date/time	09/26/2015 2:00 PM

Results Performing Laboratory

Element	Data
Laboratory Name	Century Hospital
Organization identifier	24D9871327
Address	2070 Test Park Los Angeles CA 90067 USA
Director Name	Dr. Phil J. Knowsalot III
Director identifier	5432178916

-Specimen Information-

Element	Data	
Specimen Type	BLD	
Alt Specimen Type	Blood	
Specimen Original Text	Blood Specimen	
Start date/time	201509251400	
Specimen Condition	Cool	
Alt Specimen Condition	Cool	
Condition Original Text	Cool	

-Lab results-

Element				Data				
Test performed			Erythrocyte sedimentation rate					
Test Report date			09/26/2015 14:05					
Result Observation Name	Result	UOM	Range	Abnormal Flag	Status	Date/Time of Observation		Notes
Erythrocyte sedimentation rate	20	millimeter per hour	0 to 17	Н	С	09/25/2015 14:00	09/26/2015 13:05	Specimen re- analyzed per request of ordering provider.

MSH|^~\&#|NIST Test Lab APP^372520^L|NIST Lab Facility^372521^L||NIST EHR Facility^372523^L|20150926140551||ORU^R01^ORU

_R01|LRI_1.0_2.1-NG|D|2.5.1|||AL|AL|||||LRI_Common_Component^LRI Base Profile^2.16.840.1.113883.9.16^ISO~LRI_NG_Component t^LRI Base Profile^2.16.840.1.113883.9.13^ISO~LRI_FRU_Component^LRI NG Profile^2.16.840.1.113883.9.83^ISO

 $PID|1||PATID1234^{\wedge\wedge}NIST\ MPI^MR||Jones^William^A^{\wedge\wedge\wedge}L||19610615|M||2106-3^White^HL70005^{\wedge\wedge\wedge\wedge}White\sim 1002-5^American\ India$

n or Alaska Native^HL70005^^^^^American Indian|||||||PATID1234^^^NIST MPI^AN

ORC|RE|ORD723222^NIST EHR|R-783274^NIST Lab Filler|GORD874211^NIST EHR|||||||||5742200012^Radon^Nicholas^M^JR^DR^^^NPI^L ^^^NPI

 $OBR|1|ORD723222^{NIST}EHR|R-783274^{NIST}Lab\ Filler|30341-2^{Erythrocyte}\ sedimentation\ rate^{LN^500^{Erythrocyte}\ sedimentation\ rate^{99USL^2.52^{Erythrocyte}\ sedimentation\ rate}||201509251400||||L||||5742200012^{Radon^{Nicholas^{M^{PI^{N-N}PI^$

NTE|1|Patient is extremely anxious about needles used for drawing blood.\br\If patient is overly frightened, nervous, or anxious please reschedule blood draw.

NTE|2||Patient is allergic to latex

TQ1|1||||||201509251400|201509261400|R^Routine^HL70485^^^^Routine

 $OBX|1|NM|30341-2^Erythrocyte \ sedimentation \ rate^LN^815117^ESR^99USL^2.52^Erythrocyte \ sedimentation \ rate|^1^1^1|20|mm'h^millimeter \ per \ hour^UCUM^^^1.9|0 \ to \ 17|H|||C|||201509251400|||||20150926130550||||Century\ Hospital^^^^CLIA^XX^^^224D98 \ 71327|2070\ Test\ Park^^Los$

Angeles^CA^90067^USA^B^06037|5432178916^Knowsalot^Phil^J.^IIIA^Dr.^^^NPI^L^^^NPI||||RSLT|UNSP

NTE|1||Specimen re-analyzed per request of ordering provider.

SPM|1|S-2015-66&GoodHealthC_EHR^S-9911-33&NIST Lab Filler||119297000^BLD^SCT^BldSpc^Blood^99USL^201509USEd^^Blood Specimen||||||||||||201509251400|||||||COOL^Cool^HL70493^CL^Cool^99USL^^^Cool