

NIST HL7 V2 Validation Tool Test Data Electronic Directory Of Service (eDOS)

Version 1.0

13th March, 2015

eDOS Test Case EDOS_2.3_3.1-M04_NG

M04_Update_reactivate

Test Story

Description

The originator of the eDOS is the Century Hospital Clinical Laboratory located at 2070 Test Park, Los Angeles, CA, 90067, CLIA: 24D9871327 , phone number: (310) 461-3666.
The Sender is the Century Hospital Clinical Laboratory LIS.
The Receiver is Dr. Radon's Office EHR.

The Laboratory Test Compendium from Century Hospital Clinical Laboratory is a subset of their Directory of Services. It based upon the most commonly ordered lab tests by Dr. Radon. This compendium is being used by Dr. Nicholas Radon's practice EHR to place orders electronically to Century Hospital's Clinical Laboratory. Century Hospital has made changes to their test compendium:
- Reactivated the CPT code for the previously de-activated Sed Rate.
Century Hospital's LIS transmits the appropriate eDOS update message with these changes to Dr. Radon's EHR.

Comments

Update changes for reactivation of charge codes (CPT codes).

PreCondition

EDOS_1.0_1.1-M08_New
and
EDOS_1.0_2.1-M10_New
and
EDOS_1.0_3.1-M04_New
and
EDOS_2.3_1.1-M08_Update_Reactivate
and
EDOS_2.3_2.1-M10_Update_Reactivate
have been completed.

PostCondition

No Post-Conditions.

TestObjectives

- Demonstrate capability to support a M04 update message - action: reactivate.

Notes to Testers

No Notes.

Message Content Data Sheet

Test Case Information

| EDOS_2.3_3.1-M04-NG - M04_Update_reactivate | |
|---|---------------------|
| Test Case ID | EDOS_2.3_3.1-M04_NG |

MSH :

| Location | Data Element | Data | Categorization |
|-------------|---------------------------------|------------------------|----------------------|
| MSH.1 | Field Separator | | IG Fixed Data |
| MSH.2 | Encoding Characters | ^~\&# | IG Fixed Data |
| MSH.3 | Sending Application | | |
| MSH.3.1 | Namespace ID | NIST Test Lab APP | Configurable Data |
| MSH.4 | Sending Facility | | |
| MSH.4.1 | Namespace ID | NIST Lab Facility | Configurable Data |
| MSH.6 | Receiving Facility | | |
| MSH.6.1 | Namespace ID | NIST EHR Facility | Configurable Data |
| MSH.7 | Date/Time Of Message | | |
| MSH.7.1 | Time | 20130421113601-0700 | System Generated |
| MSH.9 | Message Type | | |
| MSH.9.1 | Message Code | MFN | IG Fixed Data |
| MSH.9.2 | Trigger Event | M04 | IG Fixed Data |
| MSH.9.3 | Message Structure | MFN_M04 | IG Fixed Data |
| MSH.10 | Message Control ID | EDOS_2.3_3.1-M04-NG | System Generated |
| MSH.11 | Processing ID | | |
| MSH.11.1 | Processing ID | T | Changeable Data |
| MSH.12 | VersionID | | |
| MSH.12.1 | Version ID | 2.5.1 | IG Fixed Data |
| MSH.15 | Accept Acknowledgment Type | AL | Changeable Data |
| MSH.16 | Application Acknowledgment Type | NE | Changeable Data |
| MSH.21[1] | Message Profile Identifier | | |
| MSH.21[1].1 | Entity Identifier | EDOS_Common_Component | Test Case Fixed Data |
| MSH.21[1].2 | Namespace ID | EDOS Base Profile | Changeable Data |
| MSH.21[1].3 | Universal ID | 2.16.840.1.113883.9.67 | Test Case Fixed Data |
| MSH.21[1].4 | Universal ID Type | ISO | IG Fixed Data |

| | | | |
|------------------|-----------------------------------|------------------------|----------------------|
| MSH.21[2] | Message Profile Identifier | | |
| MSH.21[2].1 | Entity Identifier | EDOS_NG_Component | Test Case Fixed Data |
| MSH.21[2].2 | Namespace ID | EDOS NG Profile | Changeable Data |
| MSH.21[2].3 | Universal ID | 2.16.840.1.113883.9.69 | Test Case Fixed Data |
| MSH.21[2].4 | Universal ID Type | ISO | IG Fixed Data |

MFI :

| Location | Data Element | Data | Categorization |
|--------------|-------------------------------|--------------------------------|----------------------|
| MFI.1 | Master File Identifier | | |
| MFI.1.1 | Identifier | CDM | IG Fixed Data |
| MFI.1.2 | Text | Charge description master file | Changeable Data |
| MFI.1.3 | Name of Coding System | HL70175 | IG Fixed Data |
| MFI.1.7 | Coding System Version ID | 2.5.1 | Test Case Fixed Data |
| MFI.3 | File-Level Event Code | UPD | Test Case Fixed Data |
| MFI.6 | Response Level Code | NE | IG Fixed Data |

MFE :

| Location | Data Element | Data | Categorization |
|--------------|--------------------------------|--------------------------------|----------------------|
| MFE.1 | Record-Level Event Code | MAC | Test Case Fixed Data |
| MFE.3 | Effective Date/Time | | |
| MFE.3.1 | Time | 20131219145310 | System Generated |
| MFE.4 | Primary Key Value - MFE | | |
| MFE.4.1 | Identifier | 500 | Changeable Data |
| MFE.4.2 | Text | Erythrocyte sedimentation rate | Changeable Data |
| MFE.4.3 | Name of Coding System | 99USI | Changeable Data |
| MFE.4.7 | Coding System Version ID | 20130421 | Changeable Data |
| MFE.5 | Primary Key Value Type | CWE | IG Fixed Data |

CDM :

| Location | Data Element | Data | Categorization |
|--------------|--------------------------------|--------------------------------|----------------------|
| CDM.1 | Primary Key Value - CDM | | |
| CDM.1.1 | Identifier | 500 | Test Case Fixed Data |
| CDM.1.2 | Text | Erythrocyte sedimentation rate | Changeable Data |
| CDM.1.3 | Name of Coding System | 99USI | Changeable Data |

NIST HL7 V2 eDOS Validation Tool Test Data

| | | | |
|-----------------|---------------------------------|---|-----------------|
| CDM.1.7 | Coding System Version ID | 20130421 | Changeable Data |
| CDM.3 | Charge Description Short | Erythrocyte sedimentation rate | Changeable Data |
| CDM.7[1] | Procedure Code | | |
| CDM.7[1].1 | Identifier | 85652 | Changeable Data |
| CDM.7[1].2 | Text | Sedimentation rate, erythrocyte; automated | Changeable Data |
| CDM.7[1].3 | Name of Coding System | C4 | IG Fixed Data |
| CDM.7[1].7 | Coding System Version ID | 19.2 | Changeable Data |

Test Data Specification

Master File Identification

| Element name | Data |
|-----------------------|---|
| File-Level Event Code | Change file records as defined in the record-level event codes for each record that follows |

Charge description information [1]

Master File Entry

| Element name | Data |
|-------------------------|--------------------------------|
| Record-Level Event Code | Reactivate deactivated record |
| Effective Date/Time | 2013-12-19 2:53:10pm |
| Primary Key Value - MFE | Erythrocyte sedimentation rate |

Charge description

| Element name | Data |
|--------------------------|--|
| Charge Description Short | Erythrocyte sedimentation rate |
| Procedure Code [1] | Sedimentation rate, erythrocyte; automated |

HL7 Message

MSH|^~\&#|NIST Test Lab APP|NIST Lab Facility||NIST EHR Facility|20130421113601-0700||MFN^M04^MFN_M04|EDOS_2.3_3.1-M04-NG|T|2.5.1|||AL|NE||||EDOS_Common_Component^EDOS Base Profile^2.16.840.1.113883.9.67^ISO~EDOS_NG_Component^EDOS NG Profile^2.16.840.1.113883.9.69^ISO

MFI|CDM^Charge description master file^HL70175^^^^2.5.1||UPD|||NE

MFE|MAC||20131219145310|500^Erythrocyte sedimentation rate^99USI^^^^20130421|CWE

CDM|500^Erythrocyte sedimentation rate^99USI^^^^20130421||Erythrocyte sedimentation rate||||85652^Sedimentation rate, erythrocyte; automated^C4^^^^19.2