- Test Case Information-

NIST-LOI_1.0_2.1-GU_CP - Sed_Rate

Test Case ID

LOI_1.0_2.1-GU_CP

-MSH-

Location	Data Element	Data	Categorization
MSH.1	Field Separator		IG Fixed Data
MSH.2	Encoding Characters	^\&#</td><td>IG Fixed Data</td></tr><tr><td>MSH.3</td><td>Sending Application</td><td></td><td></td></tr><tr><td>MSH.3.1</td><td>Namespace ID</td><td>NIST EHR</td><td>Configurable Data</td></tr><tr><td>MSH.3.2</td><td>Universal ID</td><td>2.16.840.1.113883.3.72.5.22</td><td>Configurable Data</td></tr><tr><td>MSH.3.3</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr><tr><td>MSH.4</td><td>Sending Facility</td><td></td><td> </td></tr><tr><td>MSH.4.1</td><td>Namespace ID</td><td>NIST EHR Facility</td><td>Configurable Data</td></tr><tr><td>MSH.4.2</td><td>Universal ID</td><td>2.16.840.1.113883.3.72.5.23</td><td>Configurable Data</td></tr><tr><td>MSH.4.3</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr><tr><td>MSH.5</td><td>Receiving Application</td><td>iso iso</td><td>10 1 Act Data</td></tr><tr><td>MSH.5.1</td><td>Namespace ID</td><td>NIST Test Lab APP</td><td>Configurable Data</td></tr><tr><td>MSH.5.2</td><td>Universal ID</td><td>2.16.840.1.113883.3.72.5.20</td><td>Configurable Data</td></tr><tr><td></td><td></td><td></td><td>IG Fixed Data</td></tr><tr><td>MSH.5.3</td><td>Universal ID Type</td><td>ISO</td><td>IO FIXED Data</td></tr><tr><td>MSH.6</td><td>Receiving Facility</td><td>NICE L. L. F TV</td><td>CC. 11 D :</td></tr><tr><td>MSH.6.1</td><td>Namespace ID</td><td>NIST Lab Facility</td><td>Configurable Data</td></tr><tr><td>MSH.6.2</td><td>Universal ID</td><td>2.16.840.1.113883.3.72.5.21</td><td>Configurable Data</td></tr><tr><td>MSH.6.3</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr><tr><td>MSH.7</td><td>Date/Time Of Message</td><td></td><td></td></tr><tr><td>MSH.7.1</td><td>Time</td><td>20130211184101-0500</td><td>System Generated</td></tr><tr><td>MSH.9</td><td>Message Type</td><td></td><td></td></tr><tr><td>MSH.9.1</td><td>Message Code</td><td>OML</td><td>IG Fixed Data</td></tr><tr><td>MSH.9.2</td><td>Trigger Event</td><td>O21</td><td>IG Fixed Data</td></tr><tr><td>MSH.9.3</td><td>Message Structure</td><td>OML_O21</td><td>IG Fixed Data</td></tr><tr><td>MSH.10</td><td>Message Control ID</td><td>NIST-LOI_1.0_2.1-GU_CP</td><td>System Generated</td></tr><tr><td>MSH.11</td><td>Processing ID</td><td></td><td></td></tr><tr><td>MSH.11.1</td><td>Processing ID</td><td>T</td><td>Changeable Data</td></tr><tr><td>MSH.12</td><td>Version ID</td><td></td><td></td></tr><tr><td>MSH.12.1</td><td>Version ID</td><td>2.5.1</td><td>IG Fixed Data</td></tr><tr><td>MSH.15</td><td>Accept Acknowledgment Type</td><td>AL</td><td>IG Fixed Data</td></tr><tr><td>MSH.16</td><td>Application Acknowledgment Type</td><td>AL</td><td>Changeable Data</td></tr><tr><td>MSH.21</td><td>Message Profile Identifier</td><td></td><td></td></tr><tr><td>MSH.21.1</td><td>Entity Identifier</td><td>LOI Common Component</td><td>Test Case Fixed Data</td></tr><tr><td>MSH.21.2</td><td>Namespace ID</td><td>LOI Base Profile</td><td>Changeable Data</td></tr><tr><td>MSH.21.3</td><td>Universal ID</td><td>2.16.840.1.113883.9.66</td><td>Test Case Fixed Data</td></tr><tr><td>MSH.21.4</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr><tr><td>MSH.21[2]</td><td>Message Profile Identifier</td><td></td><td></td></tr><tr><td>MSH.21[2].1</td><td>Entity Identifier</td><td>LOI_GU_Component</td><td>Test Case Fixed Data</td></tr><tr><td>MSH.21[2].2</td><td>Namespace ID</td><td>LOI GU Profile</td><td>Changeable Data</td></tr><tr><td>MSH.21[2].3</td><td>Universal ID</td><td>2.16.840.1.113883.9.78</td><td>Test Case Fixed Data</td></tr><tr><td>MSH.21[2].4</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr><tr><td></td><td></td><td>150</td><td>IO FIXEU Data</td></tr><tr><td>MSH.21[3]</td><td>Message Profile Identifier Entity Identifier</td><td>I AD DDII Commonant</td><td>Test Case Fixed Data</td></tr><tr><td>MSH.21[3].1</td><td>-</td><td>LAB_PRU_Component</td><td></td></tr><tr><td>MSH.21[3].2</td><td>Namespace ID</td><td>LOI PRU Profile</td><td>Changeable Data</td></tr><tr><td>MSH.21[3].3</td><td>Universal ID</td><td>2.16.840.1.113883.9.82</td><td>Test Case Fixed Data</td></tr><tr><td>MSH.21[3].4</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr></tbody></table>	

_		_	_
-1	ח	п	
П		п	

Location	Data Element	Data	Categorization

PID.1Location	Set ID - PID Data Element	1 Data	IG Fixedate orization
PID.3	Patient Identifier List		
PID.3.1	ID Number	PATID1234	Configurable Data
PID.3.4	Assigning Authority		
PID.3.4.1	Namespace ID	NIST MPI	Configurable Data
PID.3.4.2	Universal ID	2.16.840.1.113883.3.72.5.30.2	Configurable Data
PID.3.4.3	Universal ID Type	ISO	IG Fixed Data
PID.3.5	Identifier Type Code	MR	Changeable Data
PID.5	Patient Name		
PID.5.1	Family Name		
PID.5.1.1	Surname	Jones	Changeable Data
PID.5.2	Given Name	William	Changeable Data
PID.5.3	Second and Further Given Names or Initials Thereof	A	Changeable Data
PID.5.4	Suffix (e.g., JR or III)	JR	Changeable Data
PID.5.7	Name Type Code	L	Changeable Data
PID.7	Date/Time of Birth		
PID.7.1	Date/Time	19610615	Changeable Data
PID.8	Administrative Sex	M	Changeable Data
PID.10	Race		
PID.10.1	Identifier	2106-3	Changeable Data
PID.10.2	Text	White	Changeable Data
PID.10.3	Name of Coding System	HL70005	IG Fixed Data
PID.11	Patient Address		
PID.11.1	Street Address		
PID.11.1.1	Street or Mailing Address	2100 Kennwood Ave	Changeable Data
PID.11.2	Other Designation	Apt 41	Changeable Data
PID.11.3	City	Los Angeles	Changeable Data
PID.11.4	State or Province	CA	Changeable Data
PID.11.5	Zip or Postal Code	90067	Changeable Data
PID.11.7	Address Type	Н	Test Case Fixed Data

-ORC-

Location	Data Element	Data	Categorization
ORC.1	Order Control	CA	Test Case Fixed Data
ORC.2	Placer Order Number		
ORC.2.1	Entity Identifier	ORD10	Changeable Data
ORC.2.2	Namespace ID	NIST EHR	Configurable Data
ORC.2.3	Universal ID	2.16.840.1.113883.3.72.5.24	Configurable Data
ORC.2.4	Universal ID Type	ISO	IG Fixed Data
ORC.4	Placer Group Number		
ORC.4.1	Entity Identifier	GORD874211	Changeable Data
ORC.4.2	Namespace ID	NIST EHR	Configurable Data
ORC.4.3	Universal ID	2.16.840.1.113883.3.72.5.24	Configurable Data
ORC.4.4	Universal ID Type	ISO	IG Fixed Data
ORC.9	Date/Time of Transaction		
ORC.9.1	Time	201103311630	System Generated
ORC.12	Ordering Provider		
ORC.12.1	ID Number	5742200012	Changeable Data
ORC.12.2	Family Name		
ORC.12.2.1	Surname	Radon	Changeable Data
ORC.12.3	Given Name	Nicholas	Changeable Data
ORC.12.9	Assigning Authority		
ORC.12.9.1	Namespace ID	NPI	Configurable Data
ORC.12.9.2	Universal ID	2.16.840.1.113883.4.6	Configurable Data
ORC.12.9.3	Universal ID Type	ISO	IG Fixed Data
ORC.12.10	Name Type Code	L	Changeable Data

-OBR				
Location	Data Element	Data	Categorization	
OBR.1	Set ID - OBR	1	IG Fixed Data	
OBR.2	Placer Order Number			
OBR.2.1	Entity Identifier	ORD10	Changeable Data	
OBR.2.2	Namespace ID	NIST EHR	Configurable Data	
OBR.2.3	Universal ID	2.16.840.1.113883.3.72.5.24	Configurable Data	
OBR.2.4	Universal ID Type	ISO	IG Fixed Data	
OBR.4	Universal Service Identifier			
OBR.4.1	Identifier	500	Test Case Fixed Data	
OBR.4.2	Text	Erythrocyte sedimentation rate	Changeable Data	
OBR.4.3	Name of Coding System	99USL	Test Case Fixed Data	
OBR.4.4	Alternate Identifier	30341-2	Test Case Fixed Data	
OBR.4.5	Alternate Text	Erythrocyte sedimentation rate	Changeable Data	
OBR.4.6	Name of Alternate Coding System	LN	Changeable Data	
OBR.4.7	Coding System Version ID	20130421	Changeable Data	
OBR.16	Ordering Provider			
OBR.16.1	ID Number	5742200012	Changeable Data	
OBR.16.2	Family Name			
OBR.16.2.1	Surname	Radon	Changeable Data	
OBR.16.3	Given Name	Nicholas	Changeable Data	
OBR.16.9	Assigning Authority			
OBR.16.9.1	Namespace ID	NPI	Configurable Data	
OBR.16.9.2	Universal ID	2.16.840.1.113883.4.6	Configurable Data	
OBR.16.9.3	Universal ID Type	ISO	IG Fixed Data	
OBR.16.10	Name Type Code	L	Changeable Data	

NPI

Data

Changeahlegorization

Changeable Data

Identifier Typa far Element

Identifier Type Code

OR Cocation

OBR.16.13

N	_	г.
	ш	Н.

NPI

Location	Data Element	Data	Categorization
NTE.1	Set ID - NTE	1	IG Fixed Data
NTE.3	Comment	Cancel Order Request	Changeable Data