-		. •
1 1000	crin	h∩n

Ms. Janet Anderson is a 26 year old white female who presented with fever of unknown origin and generalized joint pain. In order to rule out infection and autoimmune disease, Dr. Nicholas Radon ordered an erythrocyte sedimentation rate (ESR or Sed Rate) blood test to be performed. A blood specimen for the ESR test was collected from the patient and sent to the clinical lab for processing. The analysis could not be performed because the specimen arrived clotted. No result value could be obtained. A Specimen Rejection message for this lab test was generated by the LIS and transmitted to the patient's record in the ambulatory EHR used in Dr. Radon's office practice.

-Comments

This test case tests handling of specimen rejection.

PreCondition-

No Pre-condition.

-PostCondition

No Post Condition.

TestObjectives

- Determine if a valid Rejected Specimen message for an erythrocyte sedimentation rate (ESR) blood test can be created.
- Special case test aspects include ensuring that the system can populate the data elements required for a specimen rejection test message.

Notes to Testers

Note that data elements SPM-21 and SPM-24 is provided in the test data. These elements are optional for the sender and required-but maybe empty for the receiver. These elements are not subject to conformance testing for this test case (sender).

-Test Case Information-

LRI_1.2_1.1-GU - Rejected SED Rate Message (No OBX segment; OBR.25 = X)

Test Case ID LRI_1.2_1.1-GU

	ra:	

Location	Data Element	Data	Categorization
MSH.1	Field Separator		IG Fixed Data
MSH.2	Encoding Characters	\~\&	IG Fixed Data
MSH.3	Sending Application		
MSH.3.2	Universal ID	2.16.840.1.113883.3.72.5.20	Configurable Data
MSH.3.3	Universal ID Type	ISO	IG Fixed Data
MSH.4	Sending Facility		
MSH.4.2	Universal ID	2.16.840.1.113883.3.72.5.21	Configurable Data
MSH.4.3	Universal ID Type	ISO	IG Fixed Data
MSH.6	Receiving Facility		
MSH.6.2	Universal ID	2.16.840.1.113883.3.72.5.23	Configurable Data
MSH.6.3	Universal ID Type	ISO	IG Fixed Data
MSH.7	Date/Time Of Message		
MSH.7.1	Time	20150926140551	System Generated
MSH.9	Message Type		
MSH.9.1	Message Code	ORU	IG Fixed Data
MSH.9.2	Event Type	R01	IG Fixed Data
MSH.9.3	Message Structure	ORU_R01	IG Fixed Data
MSH.10	Message Control ID	LRI_1.2_1.1-GU	System Generated
MSH.11	Processing ID		
MSH.11.1	Processing ID	D	Changeable Data
MSH.12	VersionID		
MSH.12.1	Version ID	2.5.1	IG Fixed Data
MSH.15	Accept Acknowledgment Type	AL	IG Fixed Data
MSH.16	Application Acknowledgment Type	AL	IG Fixed Data
MSH.21	Message Profile Identifier		
MSH.21.1	Entity Identifier	LRI_Common_Component	IG Fixed Data
MSH.21.3	Universal ID	2.16.840.1.113883.9.16	IG Fixed Data
MSH.21.4	Universal ID Type	ISO	IG Fixed Data
MSH.21[2]	Message Profile Identifier		
MSH.21[2].1	Entity Identifier	LRI_GU_Component	IG Fixed Data
MSH.21[2].3	Universal ID	2.16.840.1.113883.9.12	IG Fixed Data
MSH.21[2].4	Universal ID Type	ISO	IG Fixed Data
MSH.21[3]	Message Profile Identifier		
MSH.21[3].1	Entity Identifier	LRI_FRU_Component	IG Fixed Data
MSH.21[3].3	Universal ID	2.16.840.1.113883.9.83	IG Fixed Data
MSH.21[3].4	Universal ID Type	ISO	IG Fixed Data

-PID-

Location	Data Element	Data	Categorization
PID.1	Set ID - PID	1	IG Fixed Data
PID.3	Patient Identifier List		
PID.3.1	ID Number	PATID1236	Configurable Data
PID.3.4	Assigning Authority		
PID.3.4.2	Universal ID	2.16.840.1.113883.3.72.5.30.2	Configurable Data
PID.3.4.3	Universal ID Type	ISO	Changeable Data
PID.3.5	Identifier Type Code	MR	Configurable Data
PID.5	Patient Name		
PID.5.1	Family Name		
PID.5.1.1	Surname	Anderson	Changeable Data

PID & Cation	Given NameData Element	Janet Data	Changeallegorization
PID.5.7	Name Type Code	L	Changeable Data
PID.7	Date/Time of Birth		
PID.7.1	Time	19860930	Changeable Data
PID.8	Administrative Sex	F	Changeable Data
PID.10	Race		
PID.10.1	Identifier	2106-3	Changeable Data
PID.10.2	Text	White	Changeable Data
PID.10.3	Name of the Coding System	HL70005	Changeable Data
PID.18	Patient Account Number		
PID.18.1	ID Number	PATID1236	Configurable Data
PID.18.4	Assigning Authority		
PID.18.4.2	Universal ID	2.16.840.1.113883.3.72.5.30.2	Configurable Data
PID.18.4.3	Universal ID Type	ISO	Changeable Data
PID.18.5	Identifier Type Code	AN	Configurable Data

-ORC-

Logadion	Data Flamont	Dete	Catagorigation
Location	Data Element	Data	Categorization
ORC.1	Order Control	RE	Test Case Fixed Data
ORC.2	Placer Order Number		
ORC.2.1	Entity Identifier	ORD723222-1	Changeable Data
ORC.2.3	Universal ID	2.16.840.1.113883.3.72.5.24	Changeable Data
ORC.2.4	Universal ID Type	ISO	IG Fixed Data
ORC.3	Filler Order Number		
ORC.3.1	Entity Identifier	R-783274-1	Changeable Data
ORC.3.3	Universal ID	2.16.840.1.113883.3.72.5.25	Changeable Data
ORC.3.4	Universal ID Type	ISO	IG Fixed Data
ORC.4	Placer Group Number		
ORC.4.1	Entity Identifier	GORD874222	Changeable Data
ORC.4.3	Universal ID	2.16.840.1.113883.3.72.5.24	Changeable Data
ORC.4.4	Universal ID Type	ISO	IG Fixed Data
ORC.12	Ordering Provider		
ORC.12.1	ID Number	5742200012	Changeable Data
ORC.12.2	Family Name		
ORC.12.2.1	Surname	Radon	Changeable Data
ORC.12.3	Given Name	Nicholas	Changeable Data
ORC.12.9	Assigning Authority		
ORC.12.9.2	Universal ID	2.16.840.1.113883.4.6	Changeable Data
ORC.12.9.3	Universal ID Type	ISO	Changeable Data
ORC.12.10	Name Type Code	L	Changeable Data
ORC.12.13	Identifier Type Code	NPI	Changeable Data

-OBR-

· Date			
Location	Data Element	Data	Categorization
OBR.1	Set ID - OBR	1	IG Fixed Data
OBR.2	Placer Order Number		
OBR.2.1	Entity Identifier	ORD723222-1	Changeable Data
OBR.2.3	Universal ID	2.16.840.1.113883.3.72.5.24	Changeable Data
OBR.2.4	Universal ID Type	ISO	IG Fixed Data
OBR.3	Filler Order Number		
OBR.3.1	Entity Identifier	R-783274-1	Changeable Data
OBR.3.3	Universal ID	2.16.840.1.113883.3.72.5.25	Changeable Data
OBR.3.4	Universal ID Type	ISO	IG Fixed Data
OBR.4	Universal Service Identifier		
OBR.4.1	Identifier	30341-2	Test Case Fixed Data
OBR.4.2	Text	Erythrocyte sedimentation rate	Test Case Fixed Data
OBR.4.3	Name of Coding System	LN	Test Case Fixed Data
	I		I

OBRocation OBR	Alternate IdpatificElement	500 Data	Changeallegorization
OBR.4.5	Alternate Text	Erythrocyte sedimentation rate	Changeable Data
OBR.4.6	Name of Alternate Coding System	99USL	Changeable Data
OBR.4.7	Coding System Version	2.52	Changeable Data
OBR.4.9	Original Text	Erythrocyte sedimentation rate	Changeable Data
OBR.7	Observation Date/Time		
OBR.7.1	Time	201509251400	Changeable Data
OBR.16	Ordering Provider		
OBR.16.1	ID Number	5742200012	Changeable Data
OBR.16.2	Family Name		
OBR.16.2.1	Surname	Radon	Changeable Data
OBR.16.3	Given Name	Nicholas	Changeable Data
OBR.16.9	Assigning Authority		
OBR.16.9.2	Universal ID	2.16.840.1.113883.4.6	Changeable Data
OBR.16.9.3	Universal ID Type	ISO	IG Fixed Data
OBR.16.10	Name Type Code	L	Changeable Data
OBR.16.13	Identifier Type Code	NPI	Changeable Data
OBR.22	Results Rpt/Status Chng - Date/Time		
OBR.22.1	Time	20150926140551	Changeable Data
OBR.25	Result Status	X	Test Case Fixed Data

-OBX-

Location	Data Element	Data	Categorization
OBX.1	Set ID - OBX	1	IG Fixed Data
OBX.2	Value Type	ST	Test Case Fixed Data
OBX.3	Observation Identifier		
OBX.3.1	Identifier	30341-2	Test Case Fixed Data
OBX.3.2	Text	Erythrocyte sedimentation rate	Test Case Fixed Data
OBX.3.3	Name of the Coding System	LN	Test Case Fixed Data
OBX.3.4	Alternate Identifier	815117	Changeable Data
OBX.3.5	Alternate Text	ESR	Changeable Data
OBX.3.6	Name of Alternate Coding System	99USL	Changeable Data
OBX.3.7	Coding System Version	2.52	Changeable Data
OBX.3.9	Original Text	Erythrocyte sedimentation rate	Changeable Data
OBX.4	Observation Sub-ID		
OBX.4.2	Group	1	Test Case Fixed Data
OBX.4.3	Sequence	1	Test Case Fixed Data
OBX.4.4	Identifier	1	Test Case Fixed Data
OBX.5	Observation Value	Test could not be performed, see Note for details	Changeable data
OBX.11	Observation Result Status	X	Test Case Fixed Data
OBX.14	Date/Time of the Observation		
OBX.14.1	Time	201509251400	Changeable Data
OBX.19	Date/Time of the Analysis		
OBX.19.1	Time	20150926130550	Changeable Data
OBX.23	Performing Organization Name		
OBX.23.1	Organization Name	Century Hospital	Changeable Data
OBX.23.6	Assigning Authority		
OBX.23.6.1	Namespace ID	CLIA	Changeable Data
OBX.23.6.2	Universal ID	2.16.840.1.113883.4.7	Changeable Data
OBX.23.6.3	Universal ID Type	ISO	IG Fixed Data
OBX.23.7	Identifier Type Code	XX	Changeable Data
OBX.23.10	Organization Identifier	24D9871327	Changeable Data
OBX.24	Performing Organization Address		
OBX.24.1	Street Address		
OBX.24.1.1	Street or Mailing Address	2070 Test Park	Changeable Data
OBX.24.3	City	Los Angeles	Changeable Data

OB Kocation	State or Propined Element	CA Data	Changeallegorization
OBX.24.5	Zip or Postal Code	90067	Changeable Data
OBX.24.6	Country	USA	Changeable Data
OBX.24.7	Address Type	В	Changeable Data
OBX.24.9	County/Parish Code	06037	Changeable Data
OBX.25	Performing Organization Medical Director		
OBX.25.1	ID Number	5432178916	Changeable Data
OBX.25.2	Family Name		
OBX.25.2.1	Surname	Knowsalot	Changeable Data
OBX.25.3	Given Name	Phil	Changeable Data
OBX.25.4	Second and Further Given Names or Initials Thereof	J.	Changeable Data
OBX.25.5	Suffix (e.g., JR or III)	III	Changeable Data
OBX.25.6	Prefix (e.g., DR)	Dr.	Changeable Data
OBX.25.9	Assigning Authority		
OBX.25.9.1	Namespace ID	NPI	Changeable Data
OBX.25.9.2	Universal ID	2.16.840.1.113883.4.6	Changeable Data
OBX.25.9.3	Universal ID Type	ISO	IG Fixed Data
OBX.25.10	Name Type Code	L	Changeable Data
OBX.25.13	Identifier Type Code	NPI	Changeable Data
OBX.29	Observation Type	RSLT	Test Case Fixed Data
OBX.30	Observation SubType	UNSP	Test Case Fixed Data

-NTE-

Location	Data Element	Data	Categorization
NTE.1	Set ID - NTE	1	IG Fixed Data
NTE.3	Comment	Blood in tube was clotted, resulting in a rejection of the specimen and leaving the lab unable to perform this test. Please resubmit a new specimen, if test is still desired.	Changeable Data

-SPM-

Location	Data Element	Data	Categorization
SPM.1	Set ID - SPM	1	IG Fixed Data
SPM.2	Specimen ID		
SPM.2.1			
SPM.2.1.1		S-2015-66	Configurable Data
SPM.2.1.3		2.16.840.1.113883.3.72.5.24	Configurable Data
SPM.2.1.4		ISO	Changeable Data
SPM.2.2			
SPM.2.2.1		S-9911-33	Changeable Data
SPM.2.2.3		2.16.840.1.113883.3.72.5.25	Changeable Data
SPM.2.2.4		ISO	IG Fixed Data
SPM.4	Specimen Type		
SPM.4.1	Identifier	119297000	Changeable Data
SPM.4.2	Text	BLD	Changeable Data
SPM.4.3	Name of the Coding System	SCT	Changeable Data
SPM.4.4	Alternate Identifier	BldSpc	Changeable Data
SPM.4.5	Alternate Text	Blood	Changeable Data
SPM.4.6	Name of Alternate Coding System	99USL	Changeable Data
SPM.4.7	Coding System Version	201509USEd	Changeable Data
SPM.4.9	Original Text	Blood Specimen	Changeable Data
SPM.17	Specimen Collection Date/Time		
SPM.17.1	Range Start Date/Time		
SPM.17.1.1	Time	201509251400	Changeable Data
SPM.21	Specimen Reject Reason		

SPMocation	Identifier Data Element	RC Data	Changeallegorization	
SPM.21.2	Text	Clotting	Changeable Data	
SPM.21.3	Name of the Coding System	HL70490	Changeable Data	
SPM.21.4	Alternate Identifier	C	Changeable Data	
SPM.21.5	Alternate Text	Clotting	Changeable Data	
SPM.21.6	Name of Alternate Coding System	99USL	Changeable Data	
SPM.21.9	Original Text	Blood specimen clotted	Changeable Data	
SPM.24	Specimen Condition			
SPM.24.1	Identifier	CLOT	Changeable Data	
SPM.24.2	Text	Clotted		
SPM.24.3	Name of the Coding System	HL70493	Changeable Data	
SPM.24.4	Alternate Identifier	C	Changeable Data	
SPM.24.5	Alternate Text	Clotted	Changeable Data	
SPM.24.6	Name of Alternate Coding System	99USL	Changeable Data	
SPM.24.9	Original Text	blood specimen clotted	Changeable Data	

-Patient Information-

Ele me nt	Data
Name	Janet Anderson
Date/Time of Birth	09/30/1986
Administrative Sex	Female
Race	White
Alt Race	

Order Observation

Ordering Provider-

Element	Data
Name	Nicholas Radon
Identifier number	5742200012

Observation Details

Observation Details————————————————————————————————————						
Element	Data					
Observation General Information						
Placer Order Number	ORD723222-1					
Filler Order Number	R-783274-1					
Placer Group Number	GORD874222					
Parent Universal Service Identifier						
Identifier						
Text						
Alt Identifier						
Alt Text						
Original Text						
Observation Details						
Universal Service Identifier	Erythrocyte sedimentation rate					
Observation Date/Time	09/25/2015 2:00 PM					
Observation end Date/Time						
Specimen Action Code						
Relevant Clinical Information						
Relevant Clinical Information Original Text						
Observation Result Information						
Result Status	X					
Results Report/Status Change - Date/Time	09/26/2015 2:05 PM					
Danika Handlina						
Results Handling						
Standard						
Observation Notes						

-Timing/Quantity Information-

Element	Data
Priority	
Start Date/time	
End Date/time	

Results Performing Laboratory			
Element	Data		
Laboratory Name	Century Hospital		
Organization identifier	24D9871327		
Address	2070 Test Park Los Angeles CA 90067 USA		
Director Name	Dr. Phil J. Knowsalot III		
Director identifier	5432178916		

Element	Data	Data		
Specimen Type	BLD			
Alt Specimen Type	Blood			
Specimen Original Text	Blood Specimen			
Start date/time	201509251400			
Specimen Reject Reason	Clotting			
Alt Specimen Reject Reason	Clotting			
Reject Reason Original Text	Blood specimen clotted			
Specimen Condition	Clotted			
Alt Specimen Condition	Clotted			
Condition Original Text	blood specimen clotted			

				—Lab resul	ts —			
Element			Data					
Test performed			Erythrocyte	Erythrocyte sedimentation rate				
Test Report d	late			09/26/2015	14:05			
Result Observation Name	Result	UOM	Range	Abnormal Flag	Status		Date/Time of Analysis	Notes
Erythrocyte sedimentation rate	Test could not be performed, see Note for details				X	09/25/2015 14:00	09/26/2015 13:05	Blood in tub was clotted, resulting in a rejection of the specimer and leaving the lab unable to perform this test. Please resubmit a new specimen, if test is still desired.

 $MSH|^{\sim} \&|^{2.16.840.1.113883.3.72.5.20} |^{2.16.840.1.113883.3.72.5.21} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.113883.3.72.5.23} |^{2.16.840.1.1138} |^{2.16.840.1.1138} |^{2.16.84$

 $GU|D|2.5.1|||AL|AL|||||LRI_Common_Component^{2}.16.840.1.113883.9.16^{I}SO\sim LRI_GU_Component^{2}.16.840.1.113883.9.16^{I}SO\sim LRI_GU_Component^{2}.16.840.1.113883.9.83^{I}SO$

 $PID|1||PATID1236^{^{\circ}}\&2.16.840.1.113883.3.72.5.30.2\&ISO^{^{\circ}}MR||Anderson^{^{\circ}}Janet^{^{\circ}}\&2.16.840.1.113883.3.72.5.30.2\&ISO^{^{\circ}}MR||Anderson^{^{\circ}}Janet^{^{\circ}}\&2.16.840.1.113883.3.72.5.30.2\&ISO^{^{\circ}}AN$

ORC|RE|ORD723222-1^2.16.840.1.113883.3.72.5.24^ISO|R-783274-1^2.16.840.1.113883.3.72.5.25^ISO|GORD874222^2.16.840.1.113883.3.72.5.24^ISO||||||||||5742200012^Radon^Nicholas^2.16.840.1.113883.4.6&ISO^L^^NPI

 $OBR|1|ORD723222-1^{2}.16.840.1.113883.3.72.5.24^{I}SO|R-783274-1^{2}.16.840.1.113883.3.72.5.25^{I}SO|30341-2^{E}rythrocyte sedimentation rate^{LN^500^{E}rythrocyte} sedimentation rate^{99USL^2.52^{E}rythrocyte} sedimentation rate|||201509251400|||||||||5742200012^{Radon^{Nicholas^{^{^{*}}}}}&2.16.840.1.113883.4.6&ISO^{L^{^{*}}}NPI|||||20150926140551||X$

 $OBX[1|ST]30341-2^{Erythrocyte}\ sedimentation\ rate^{LN^815117^{ESR^99USL^2.52^{Erythrocyte}\ sedimentation\ rate}]^{1^11}|Test\ co\ uld\ not\ be\ performed,\ see\ Note\ for\ details|||||X|||201509251400|||||20150926130550||||Century\ Hospital^^^^CLIA&2.16.84\\ 0.1.113883.4.7&ISO^XX^^24D9871327|2070\ Test\ Park^^Los\ Angeles^{CA^90067^USA^B^06037|5432178916^Knowsalot^Phil^J.^III^Dr\\ .^^NPI&2.16.840.1.113883.4.6&ISO^L^^NPI||||RSLT|UNSP$

NTE|1||Blood in tube was clotted, resulting in a rejection of the specimen and leaving the lab unable to perform this t est. Please resubmit a new specimen, if test is still desired.

 $SPM[1|S-2015-66\&\&2.16.840.1.113883.3.72.5.24\&ISO^S-9911-33\&\&2.16.840.1.113883.3.72.5.25\&ISO||119297000^BLD^SCT^BklSpc^B|\\ bod^99USL^201509USEd^Blood Specimen[||||||||||||201509251400||||RC^Clotting^HL70490^C^Clotting^99USL^^Blood specimen clotted|||CLOT^Clotted^HL70493^C^Clotted^99USL^^blood specimen clotted|||CLOT^Clotted^HL70493^C^Clotted^99USL^^Blood specimen clotted|||CLOT^Clotted^HL70493^C^Clotted^99USL^^Blood specimen clotted|||CLOT^Clotted^HL70493^C^Clotted^99USL^^Blood specimen clotted|||CLOT^Clotted^Blood^$