$\mathbf{r}$			. •
· I )	PCCT	าก	tion

Mr. William A. Jones is a 51 year old white male who presented for a cardiology exam. Dr. Nicholas Radon ordered a lipid panel test to be performed to check for hyperlipidemia. A blood specimen for the lipid panel test was collected from the patient and sent to the clinical lab at Century Hospital, 2070 Test Park, Los Angeles, CA, 90067, for processing. The laboratory director at this hospital is Dr. Phil J. Knowsalot, III. The blood specimen was analyzed, and the final result of the lipid panel test showed that all components were within the normal range. The final result report for this lab test was generated by the LIS and transmitted to the patient's record in the ambulatory EHR used in Dr. Radon's office practice. A copy of the test result report was sent to Dr. Pafford M. Hamlin, Sr.

Comments	_
Contraction	
This test case is evaluating the handling of communication of fasting status information as part of a typical Lipid Panel.	
PreCondition—	
FICCOLUMN	
No Pre-condition.	
De de Carallèra	_
PostCondition—	_
No Post Condition.	
1-0-1-001-001-001-001-001-001-001-001-0	
T. (01)	_
TestObjectives—	_
Demonstrate the capability to message a valid LRI lipid panel message that contains fasting status information.	
5 Demonstrate the capacity to message a valid Ere ipid paner message that contains fasting status information.	
	_
Notes to Testers	_

No Notes.

#### -Test Case Information-

## LRI\_3.0\_1.1-NG - Typically Populated Lipid Panel message - Final Results

Test Case ID LRI\_3.0\_1.1-NG

-MSH-

MSH————————————————————————————————————				
Location	Data Element	Data	Categorization	
MSH.1	Field Separator		IG Fixed Data	
MSH.2	<b>Encoding Characters</b>	^~\&#</td><td>IG Fixed Data</td></tr><tr><td>MSH.3</td><td>Sending Application</td><td></td><td></td></tr><tr><td>MSH.3.2</td><td>Universal ID</td><td>372520</td><td>Configurable Data</td></tr><tr><td>MSH.3.3</td><td>Universal ID Type</td><td>L</td><td>Configurable Data</td></tr><tr><td>MSH.4</td><td>Sending Facility</td><td></td><td></td></tr><tr><td>MSH.4.2</td><td>Universal ID</td><td>372521</td><td>Configurable Data</td></tr><tr><td>MSH.4.3</td><td>Universal ID Type</td><td>L</td><td>Configurable Data</td></tr><tr><td>MSH.6</td><td>Receiving Facility</td><td></td><td></td></tr><tr><td>MSH.6.2</td><td>Universal ID</td><td>372523</td><td>Configurable Data</td></tr><tr><td>MSH.6.3</td><td>Universal ID Type</td><td>L</td><td>Configurable Data</td></tr><tr><td>MSH.7</td><td>Date/Time Of Message</td><td></td><td></td></tr><tr><td>MSH.7.1</td><td>Time</td><td>20150926140551</td><td>System Generated</td></tr><tr><td>MSH.9</td><td>Message Type</td><td></td><td></td></tr><tr><td>MSH.9.1</td><td>Message Code</td><td>ORU</td><td>IG Fixed Data</td></tr><tr><td>MSH.9.2</td><td>Event Type</td><td>R01</td><td>IG Fixed Data</td></tr><tr><td>MSH.9.3</td><td>Message Structure</td><td>ORU_R01</td><td>IG Fixed Data</td></tr><tr><td>MSH.10</td><td>Message Control ID</td><td>LRI_3.0_1.1-NG</td><td>System Generated</td></tr><tr><td>MSH.11</td><td>Processing ID</td><td></td><td></td></tr><tr><td>MSH.11.1</td><td>Processing ID</td><td>D</td><td>Changeable Data</td></tr><tr><td>MSH.12</td><td>VersionID</td><td></td><td></td></tr><tr><td>MSH.12.1</td><td>Version ID</td><td>2.5.1</td><td>IG Fixed Data</td></tr><tr><td>MSH.15</td><td>Accept Acknowledgment Type</td><td>AL</td><td>IG Fixed Data</td></tr><tr><td>MSH.16</td><td>Application Acknowledgment Type</td><td>AL</td><td>IG Fixed Data</td></tr><tr><td>MSH.21</td><td>Message Profile Identifier</td><td></td><td></td></tr><tr><td>MSH.21.1</td><td>Entity Identifier</td><td>LRI_Common_Component</td><td>IG Fixed Data</td></tr><tr><td>MSH.21.3</td><td>Universal ID</td><td>2.16.840.1.113883.9.16</td><td>IG Fixed Data</td></tr><tr><td>MSH.21.4</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr><tr><td>MSH.21[2]</td><td>Message Profile Identifier</td><td></td><td></td></tr><tr><td>MSH.21[2].1</td><td>Entity Identifier</td><td>LRI_NG_Component</td><td>IG Fixed Data</td></tr><tr><td>MSH.21[2].3</td><td>Universal ID</td><td>2.16.840.1.113883.9.13</td><td>IG Fixed Data</td></tr><tr><td>MSH.21[2].4</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr><tr><td>MSH.21[3]</td><td>Message Profile Identifier</td><td></td><td></td></tr><tr><td>MSH.21[3].1</td><td>Entity Identifier</td><td>LRI_FRU_Component</td><td>IG Fixed Data</td></tr><tr><td>MSH.21[3].3</td><td>Universal ID</td><td>2.16.840.1.113883.9.83</td><td>IG Fixed Data</td></tr><tr><td>MSH.21[3].4</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr></tbody></table>		

-PID-

Data Element	Data	Categorization
Set ID - PID	1	IG Fixed Data
Patient Identifier List		
ID Number	PATID1234	Configurable Data
Assigning Authority		
Namespace ID	NIST MPI	Changeable Data
Identifier Type Code	MR	Configurable Data
Patient Name		
Family Name		
Surname	Jones	Changeable Data
Given Name	William	Changeable Data
	Set ID - PID  Patient Identifier List  ID Number  Assigning Authority  Namespace ID  Identifier Type Code  Patient Name  Family Name  Surname	Set ID - PID 1  Patient Identifier List  ID Number PATID1234  Assigning Authority  Namespace ID NIST MPI  Identifier Type Code MR  Patient Name  Family Name  Surname Jones

PIDLogation	Second and Further Given Names or Initials Thereof	A Data	Chang Cattegorization
PID.5.7	Name Type Code	L	Changeable Data
PID.7	Date/Time of Birth		
PID.7.1	Time	19610627	Changeable Data
PID.8	Administrative Sex	M	Changeable Data
PID.10	Race		
PID.10.1	Identifier	2106-3	Changeable Data
PID.10.2	Text	White	Changeable Data
PID.10.3	Name of the Coding System	HL70005	Changeable Data
PID.18	Patient Account Number		
PID.18.1	ID Number	PATID1234	Configurable Data
PID.18.4	Assigning Authority		
PID.18.4.1	Namespace ID	NIST MPI	Changeable Data
PID.18.5	Identifier Type Code	AN	Configurable Data

-ORC-

Location	Data Element	Data	Categorization
ORC.1	Order Control	RE	Test Case Fixed Data
ORC.2	Placer Order Number		
ORC.2.1	Entity Identifier	ORD777888	Changeable Data
ORC.2.2	Namespace ID	NIST EHR	Changeable Data
ORC.3	Filler Order Number		
ORC.3.1	Entity Identifier	R-220713	Changeable Data
ORC.3.2	Namespace ID	NIST Lab Filler	Changeable Data
ORC.4	Placer Group Number		
ORC.4.1	Entity Identifier	GORD874244	Changeable Data
ORC.4.2	Namespace ID	NIST EHR	Changeable Data
ORC.12	Ordering Provider		
ORC.12.1	ID Number	5742200012	Changeable Data
ORC.12.2	Family Name		
ORC.12.2.1	Surname	Radon	Changeable Data
ORC.12.3	Given Name	Nicholas	Changeable Data
ORC.12.9	Assigning Authority		
ORC.12.9.1	Namespace ID	NPI	Changeable Data
ORC.12.10	Name Type Code	L	Changeable Data
ORC.12.13	Identifier Type Code	NPI	Changeable Data

-OBR-

Location	Data Element	Data	Categorization
OBR.1	Set ID - OBR	1	IG Fixed Data
OBR.2	Placer Order Number		
OBR.2.1	Entity Identifier	ORD777888	Changeable Data
OBR.2.2	Namespace ID	NIST EHR	Changeable Data
OBR.3	Filler Order Number		
OBR.3.1	Entity Identifier	R-220713	Changeable Data
OBR.3.2	Namespace ID	NIST Lab Filler	Changeable Data
OBR.4	Universal Service Identifier		
OBR.4.1	Identifier	24331-1	Test Case Fixed Data
OBR.4.2	Text	Lipid 1996 panel in Serum or Plasma	Test Case Fixed Data
OBR.4.3	Name of Coding System	LN	Test Case Fixed Data
OBR.4.4	Alternate Identifier	345789	Changeable Data
OBR.4.5	Alternate Text	Lipid Panel	Changeable Data
OBR.4.6	Name of Alternate Coding System	99USL	Changeable Data
OBR.4.7	Coding System Version	2.52	Changeable Data
OBR.4.9	Original Text	Lipid 1996 panel in Serum or Plasma	Changeable Data
OBR.7	Observation Date/Time		

OBRocation	Time Data Element	20150925 <b>Data</b>	Changeable Dataation
OBR.13	Relevant Clinical Information		
OBR.13.1	Identifier	F	Changeable Data
OBR.13.2	Text	Patient was fasting prior to the procedure.	Changeable Data
OBR.13.3	Name of the Coding System	HL70916	Changeable Data
OBR.13.7	Coding System Version	2.7.1	Changeable Data
OBR.13.9	Original Text	fasting 12 hours	Changeable Data
OBR.16	Ordering Provider		
OBR.16.1	ID Number	5742200012	Changeable Data
OBR.16.2	Family Name		
OBR.16.2.1	Surname	Radon	Changeable Data
OBR.16.3	Given Name	Nicholas	Changeable Data
OBR.16.9	Assigning Authority		
OBR.16.9.1	Namespace ID	NPI	Changeable Data
OBR.16.10	Name Type Code	L	Changeable Data
OBR.16.13	Identifier Type Code	NPI	Changeable Data
OBR.22	Results Rpt/Status Chng - Date/Time		
OBR.22.1	Time	20150926140551	Changeable Data
OBR.25	Result Status	F	Test Case Fixed Data
OBR.28	Result Copies To		
OBR.28.1	ID Number	10092000194	Changeable Data
OBR.28.2	Family Name		
OBR.28.2.1	Surname	Hamlin	Changeable Data
OBR.28.3	Given Name	Pafford	Changeable Data
OBR.28.9	Assigning Authority		
OBR.28.9.1	Namespace ID	NPI	Changeable Data
OBR.28.10	Name Type Code	L	Changeable Data
OBR.28.13	Identifier Type Code	NPI	Changeable Data
OBR.49	Results Handling		
OBR.49.1	Identifier	CC	Test Case Fixed Data
OBR.49.2	Text	Copies Requested	Test Case Fixed Data
OBR.49.3	Name of the Coding System	HL70507	IG Fixed Data

Location	Data Element	Data	Categorization
OBX.1	Set ID - OBX	1	IG Fixed Data
OBX.2	Value Type	NM	Test Case Fixed Data
OBX.3	Observation Identifier		
OBX.3.1	Identifier	2093-3	Test Case Fixed Data
OBX.3.2	Text	Cholesterol [Mass/volume] in Serum or Plasma	Test Case Fixed Data
OBX.3.3	Name of the Coding System	LN	Test Case Fixed Data
OBX.3.7	Coding System Version	2.52	Changeable Data
OBX.3.9	Original Text	Cholesterol [Mass/volume] in Serum or Plasma	Test Case Fixed Data
OBX.5	Observation Value	196	Test Case Fixed Data
OBX.6	Units		
OBX.6.1	Identifier	mg/dL	Changeable Data
OBX.6.2	Text	milligrams per deciliter	Changeable Data
OBX.6.3	Name of the Coding System	UCUM	Changeable Data
OBX.6.7	Coding System Version	1.9	Changeable Data
OBX.7	Reference Range	Recommended: <200; Moderate Risk: 200-239; High Risk: >240	Changeable Data
OBX.8	Abnormal Flags	N	Test Case Fixed Data
OBX.11	Observation Result Status	F	Test Case Fixed Data
OBX.14	Date/Time of the Observation		
OBX.14.1	Time	20150925	Changeable Data

OBX <sub>1</sub> 18 <sub>cation</sub>	Date/Time of the Argelyseint	Data	Categorization
OBX.19.1	Time	201509261400	Changeable Data
OBX.23	Performing Organization Name		
OBX.23.1	Organization Name	Century Hospital	Changeable Data
OBX.23.6	Assigning Authority		
OBX.23.6.1	Namespace ID	CLIA	Changeable Data
OBX.23.7	Identifier Type Code	XX	Changeable Data
OBX.23.10	Organization Identifier	24D9871327	Changeable Data
OBX.24	Performing Organization Address		
OBX.24.1	Street Address		
OBX.24.1.1	Street or Mailing Address	2070 Test Park	Changeable Data
OBX.24.3	City	Los Angeles	Changeable Data
OBX.24.4	State or Province	CA	Changeable Data
OBX.24.5	Zip or Postal Code	90067	Changeable Data
OBX.24.7	Address Type	В	Changeable Data
OBX.25	Performing Organization Medical Director		
OBX.25.1	ID Number	5432178916	Changeable Data
OBX.25.2	Family Name		
OBX.25.2.1	Surname	Knowsalot	Changeable Data
OBX.25.3	Given Name	Phil	Changeable Data
OBX.25.6	Prefix (e.g., DR)	Dr.	Changeable Data
OBX.25.9	Assigning Authority		
OBX.25.9.1	Namespace ID	NPI	Changeable Data
OBX.25.10	Name Type Code	L	Changeable Data
OBX.25.13	Identifier Type Code	NPI	Changeable Data
OBX.29	Observation Type	RSLT	Test Case Fixed Data

Location	Data Element	Data	Categorization
OBX.1	Set ID - OBX	2	IG Fixed Data
OBX.2	Value Type	NM	Test Case Fixed Data
OBX.3	Observation Identifier		
OBX.3.1	Identifier	2571-8	Test Case Fixed Data
OBX.3.2	Text	Triglyceride [Mass/volume] in Serum or Plasma	Test Case Fixed Data
OBX.3.3	Name of the Coding System	LN	Test Case Fixed Data
OBX.3.7	Coding System Version	2.52	Changeable Data
OBX.3.9	Original Text	Triglyceride [Mass/volume] in Serum or Plasma	Test Case Fixed Data
OBX.5	Observation Value	100	Test Case Fixed Data
OBX.6	Units		
OBX.6.1	Identifier	mg/dL	Changeable Data
OBX.6.2	Text	milligrams per deciliter	Changeable Data
OBX.6.3	Name of the Coding System	UCUM	Changeable Data
OBX.6.7	Coding System Version	1.9	Changeable Data
OBX.7	Reference Range	40 to 160	Changeable Data
OBX.8	Abnormal Flags	N	Test Case Fixed Data
OBX.11	Observation Result Status	F	Test Case Fixed Data
OBX.14	Date/Time of the Observation		
OBX.14.1	Time	20150925	Changeable Data
OBX.19	Date/Time of the Analysis		
OBX.19.1	Time	201509261400	Changeable Data
OBX.23	Performing Organization Name		
OBX.23.1	Organization Name	Century Hospital	Changeable Data
OBX.23.6	Assigning Authority		
OBX.23.6.1	Namespace ID	CLIA	Changeable Data

OB Location	Identifier Typa factiement	XX Data	Changealdego Tization
OBX.23.10	Organization Identifier	24D9871327	Changeable Data
OBX.24	Performing Organization Address		
OBX.24.1	Street Address		
OBX.24.1.1	Street or Mailing Address	2070 Test Park	Changeable Data
OBX.24.3	City	Los Angeles	Changeable Data
OBX.24.4	State or Province	CA	Changeable Data
OBX.24.5	Zip or Postal Code	90067	Changeable Data
OBX.24.7	Address Type	В	Changeable Data
OBX.25	Performing Organization Medical Director		
OBX.25.1	ID Number	5432178916	Changeable Data
OBX.25.2	Family Name		
OBX.25.2.1	Surname	Knowsalot	Changeable Data
OBX.25.3	Given Name	Phil	Changeable Data
OBX.25.6	Prefix (e.g., DR)	Dr.	Changeable Data
OBX.25.9	Assigning Authority		
OBX.25.9.1	Namespace ID	NPI	Changeable Data
OBX.25.10	Name Type Code	L	Changeable Data
OBX.25.13	Identifier Type Code	NPI	Changeable Data
OBX.29	Observation Type	RSLT	Test Case Fixed Data

Location	Data Element	Data	Categorization
OBX.1	Set ID - OBX	3	IG Fixed Data
OBX.2	Value Type	NM	Test Case Fixed Data
OBX.3	Observation Identifier		
OBX.3.1	Identifier	2085-9	Test Case Fixed Data
OBX.3.2	Text	Cholesterol in HDL [Mass/volume] in Serum or Plasma	Test Case Fixed Data
OBX.3.3	Name of the Coding System	LN	Test Case Fixed Data
OBX.3.7	Coding System Version	2.52	Changeable Data
OBX.3.9	Original Text	Cholesterol in HDL [Mass/volume] in Serum or Plasma	Test Case Fixed Data
OBX.5	Observation Value	60	Test Case Fixed Data
OBX.6	Units		
OBX.6.1	Identifier	mg/dL	Changeable Data
OBX.6.2	Text	milligrams per deciliter	Changeable Data
OBX.6.3	Name of the Coding System	UCUM	Changeable Data
OBX.6.7	Coding System Version	1.9	Changeable Data
OBX.7	Reference Range	29 to 72	Changeable Data
OBX.8	Abnormal Flags	N	Test Case Fixed Data
OBX.11	Observation Result Status	F	Test Case Fixed Data
OBX.14	Date/Time of the Observation		
OBX.14.1	Time	20150925	Changeable Data
OBX.19	Date/Time of the Analysis		
OBX.19.1	Time	201509261400	Changeable Data
OBX.23	Performing Organization Name		
OBX.23.1	Organization Name	Century Hospital	Changeable Data
OBX.23.6	Assigning Authority		
OBX.23.6.1	Namespace ID	CLIA	Changeable Data
OBX.23.7	Identifier Type Code	XX	Changeable Data
OBX.23.10	Organization Identifier	24D9871327	Changeable Data
OBX.24	Performing Organization Address		
OBX.24.1	Street Address		
OBX.24.1.1	Street or Mailing Address	2070 Test Park	Changeable Data
OBX.24.3	City	Los Angeles	Changeable Data

OB Location	State or Propinca Element	CA <b>Data</b>	Changealle gorization
OBX.24.5	Zip or Postal Code	90067	Changeable Data
OBX.25	Performing Organization Medical Director		
OBX.25.1	ID Number	5432178916	Changeable Data
OBX.25.2	Family Name		
OBX.25.2.1	Surname	Knowsalot	Changeable Data
OBX.25.3	Given Name	Phil	Changeable Data
OBX.25.6	Prefix (e.g., DR)	Dr.	Changeable Data
OBX.25.9	Assigning Authority		
OBX.25.9.1	Namespace ID	NPI	Changeable Data
OBX.25.10	Name Type Code	L	Changeable Data
OBX.25.13	Identifier Type Code	NPI	Changeable Data
OBX.29	Observation Type	RSLT	Test Case Fixed Data

Location	Data Element	Data	Categorization
OBX.1	Set ID - OBX	4	IG Fixed Data
OBX.2	Value Type	NM	Test Case Fixed Data
OBX.3	Observation Identifier		
OBX.3.1	Identifier	2089-1	Test Case Fixed Data
OBX.3.2	Text	Cholesterol in LDL [Mass/volume] in Serum or Plasma	Test Case Fixed Data
OBX.3.3	Name of the Coding System	LN	Test Case Fixed Data
OBX.3.7	Coding System Version	2.52	Changeable Data
OBX.3.9	Original Text	Cholesterol in LDL [Mass/volume] in Serum or Plasma	Test Case Fixed Data
OBX.5	Observation Value	116	Test Case Fixed Data
OBX.6	Units		
OBX.6.1	Identifier	mg/dL	Changeable Data
OBX.6.2	Text	milligrams per deciliter	Changeable Data
OBX.6.3	Name of the Coding System	UCUM	Changeable Data
OBX.6.7	Coding System Version	1.9	Changeable Data
OBX.7	Reference Range	Recommended: <130; Moderate Risk: 130-159; High Risk: >160	Changeable Data
OBX.8	Abnormal Flags	N	Test Case Fixed Data
OBX.11	Observation Result Status	F	Test Case Fixed Data
OBX.14	Date/Time of the Observation		
OBX.14.1	Time	20150925	Changeable Data
OBX.19	Date/Time of the Analysis		
OBX.19.1	Time	201509261400	Changeable Data
OBX.23	Performing Organization Name		
OBX.23.1	Organization Name	Century Hospital	Changeable Data
OBX.23.6	Assigning Authority		
OBX.23.6.1	Namespace ID	CLIA	Changeable Data
OBX.23.7	Identifier Type Code	XX	Changeable Data
OBX.23.10	Organization Identifier	24D9871327	Changeable Data
OBX.24	Performing Organization Address		
OBX.24.1	Street Address		
OBX.24.1.1	Street or Mailing Address	2070 Test Park	Changeable Data
OBX.24.3	City	Los Angeles	Changeable Data
OBX.24.4	State or Province	CA	Changeable Data
OBX.24.5	Zip or Postal Code	90067	Changeable Data
OBX.24.7	Address Type	В	Changeable Data
OBX.25	Performing Organization Medical Director		
OBX.25.1	ID Number	5432178916	Changeable Data
OBX.25.2	Family Name		_

OB <del>ocation</del>	Surname Data Element	Knowsalot <b>Data</b>	Changeanegorization
OBX.25.3	Given Name	Phil	Changeable Data
OBX.25.6	Prefix (e.g., DR)	Dr.	Changeable Data
OBX.25.9	Assigning Authority		
OBX.25.9.1	Namespace ID	NPI	Changeable Data
OBX.25.10	Name Type Code	L	Changeable Data
OBX.25.13	Identifier Type Code	NPI	Changeable Data
OBX.29	Observation Type	RSLT	Test Case Fixed Data

-SPM-

Location	Data Element	Data	Categorization
SPM.1	Set ID - SPM	1	IG Fixed Data
SPM.2	Specimen ID		
SPM.2.1			
SPM.2.1.1		S2015-777888	Changeable Data
SPM.2.1.2		NIST EHR	Changeable Data
SPM.2.2			
SPM.2.2.1		S-220713-1	Changeable Data
SPM.2.2.2		NIST Lab Filler	Changeable Data
SPM.4	Specimen Type		
SPM.4.1	Identifier	119297000	Changeable Data
SPM.4.2	Text	BLD	Changeable Data
SPM.4.3	Name of the Coding System	SCT	Changeable Data
SPM.4.7	Coding System Version	201509USEd	Changeable Data
SPM.4.9	Original Text	Blood	Changeable Data
SPM.17	Specimen Collection Date/Time		
SPM.17.1	Range Start Date/Time		
SPM.17.1.1	Time	20150925	Changeable Data

### -Patient Information-

Ele me nt	Data
Name	William A Jones
Date/Time of Birth	06/27/1961
Administrative Sex	Male
Race	White
Alt Race	

#### Order Observation

### Ordering Provider-

Element	Data
Name	Nicholas Radon
Identifier number	5742200012

### Observation Details-

Element	Data
Observation General Information	Data
Placer Order Number	ORD777888
Filler Order Number	R-220713
	GORD874244
Placer Group Number	GORD8/4244
Parent Universal Service Identifier	I
Identifier	
Text	
Alt Identifier	
Alt Text	
Original Text	
Observation Details	
Universal Service Identifier	Lipid 1996 panel in Serum or Plasma
Observation Date/Time	09/25/2015
Observation end Date/Time	07/25/2015
Specimen Action Code	
Relevant Clinical Information	Patient was fasting prior to the procedure.
Relevant Clinical Information Original Text	fasting 12 hours
Observation Result Information	
Result Status	F
Results Report/Status Change - Date/Time	09/26/2015 2:05 PM
Results Copy To	- m 4
Name	Pafford
Identifier	10092000194
Results Handling	
Standard	
Observation Notes	

# Timing/Quantity Information-

Element	Data
---------	------

Priority	
Start Date/time	
End Date/time	

Results Performing Laboratory-

Element	Data
Laboratory Name	Century Hospital
Organization identifier	24D9871327
Address	2070 Test Park Los Angeles CA 90067
Director Name	Dr. Phil Knowsalot
Director identifier	5432178916

-Specimen Information-

Element	Data
Specimen Type	BLD
Alt Specimen Type	
Specimen Original Text	Blood
Start date/time	20150925

Lab results

Element				Data				
Test performed				Lipid 1996 panel in Serum or Plasma				
Test Report date				09/26/2015 14:05				
Result Observation Name	Result	UOM	Range	Abnormal Flag	Status	Date/Time of Observation	Date/Time of Analysis	Notes
Cholesterol [Mass/volume in Serum or Plasma	<sup>1</sup> 196	milligrams per deciliter	Recommender <200; Moderate Risk: 200- 239; High Risk: >240	d: N	F	09/25/2015	09/26/2015 14:00	
Triglyceride [Mass/volume in Serum or Plasma	100	milligrams per deciliter	40 to 160	N	F	09/25/2015	09/26/2015 14:00	
Cholesterol in HDL [Mass/volume in Serum or Plasma		milligrams per deciliter	29 to 72	N	F	09/25/2015	09/26/2015 14:00	
Cholesterol in LDL [Mass/volume in Serum or Plasma		milligrams per deciliter	Recommender <130; Moderate Risk: 130- 159; High Risk: >160	d: N	F	09/25/2015	09/26/2015 14:00	

MSH|\^\&#|\^372520^L|\^372521^L|\^372523^L|\20150926140551|\ORU\^R01\^ORU\_R01|\LRI\_3.0\_1.1-NG|D|2.5.1|\|AL|AL|\|||LRI\_Common\_Component\^2.16.840.1.113883.9.16\^ISO\^LRI\_NG\_Component\^2.16.840.1.113883.9.13\^ISO\^LRI\_FRU\_Component\^2.16.840.1.113883.9.13\^ISO\^LRI\_FRU\_Component\^2.16.840.1.113883.9.13\^ISO\^LRI\_FRU\_Component\^2.16.840.1.113883.9.13\^ISO\^LRI\_FRU\_Component\^2.16.840.1.113883.9.13\^ISO\^LRI\_FRU\_Component\^2.16.840.1.113883.9.13\^ISO\^LRI\_FRU\_COMPONENT\^2.16.840.1.113883.9.13\^ISO\^LRI\_FRU\_CO

PID|1||PATID1234^^^NIST MPI^MR||Jones^William^A^^^^L||19610627|M||2106-3^White^HL70005||||||||PATID1234^^^NIST MPI^AN

ORC|RE|ORD777888^NIST EHR|R-220713^NIST Lab Filler|GORD874244^NIST EHR||||||||5742200012^Radon^Nicholas^^^^NPI^L^^^NPI

 $OBR|1|ORD777888^{N}IST\ EHR|R-220713^{N}IST\ Lab\ Filler|24331-1^{L}ipid\ 1996\ panel\ in\ Serum\ or\ Plasma^{L}N^345789^{L}ipid\ Panel^99U\ SL^2.52^{L}ipid\ 1996\ panel\ in\ Serum\ or\ Plasma|||20150925||||||F^Patient\ was\ fasting\ prior\ to\ the\ procedure.^HL70916^^^2. \\ 7.1^{fasting\ 12\ hours|||5742200012^{Radon^Nicholas^^^^NPI^L^^^NPI||||||20150926140551|||F|||10092000194^{Hamlin^Pafford^^^NPI^L^^^NPI||||||||||||||||||||CC^{Copies\ Requested^HL70507}$ 

 $OBX|1|NM|2093-3^Cholesterol\ [Mass/volume]\ in\ Serum\ or\ Plasma^LN^^^2.52^^Cholesterol\ [Mass/volume]\ in\ Serum\ or\ Plasma||\\ 196|mg/dL^milligrams\ per\ deciliter^UCUM^^^1.9|Recommended: <200;\ Moderate\ Risk: 200-239\ ;\ High\ Risk: >240|N|||F|||20150925|||||201509261400||||Century\ Hospital^^^^CLIA^XX^^24D9871327|2070\ Test\ Park^^Los\ Angeles^CA^90067^^B|5432178916^Kno\ wsalot^Phil^^^Dr.^^NPI^L^^^NPI||||RSLT$ 

 $OBX[3]NM[2085-9^C Cholesterol\ in\ HDL\ [Mass/volume]\ in\ Serum\ or\ Plasma^LN^{^^2}2.52^{^C}Cholesterol\ in\ HDL\ [Mass/volume]\ in\ Serum\ or\ Plasma[60]mg/dL^milligrams\ per\ deciliter^UCUM^{^^1}1.9[29\ to\ 72]N[||F|||20150925][|||201509261400][||Century\ Hospit\ al^{^^2}2178916^Knowsalot^Phil^{^2}Dr.^{^2}NPI^{||RSLT|}2178916^Knowsalot^Phil^{^2}Dr.^{^2}NPI^{||RSLT|}2178916^Knowsalot^{^2}Dr.^{^2}NPI^{^2}Dr.^{^2}NPI^{^2}Dr.^{^2}NPI^{^2}Dr.^{^2}NPI^{^2}Dr.^{^2}NPI^{^2}Dr.^{^2}Dr$ 

 $OBX|4|NM|2089-1^Cholesterol\ in\ LDL\ [Mass/volume]\ in\ Serum\ or\ Plasma^LN^^^2.52^^Cholesterol\ in\ LDL\ [Mass/volume]\ in\ Serum\ or\ Plasma^LN^^2.52^^Cholesterol\ in\ LDL\ [Mass/volume]\ in\ Serum\ or\ Plasma^LN^^2.52^^Cholesterol\ in\ Serum\ or\ Plasma^LN^2.52^^Cholesterol\ in\ Serum\ or\ Plasma^LN^2.52^^Cholesterol\ in\ Serum\ or\ Plasma^LN^2.52^^Choleste$