

Test Case Information

EDOS 2.4 3.1-M04 GU - M04 Update combo	
Test Case ID	EDOS_2.4_3.1-M04_GU

MSH :

Location	Data Element	Data	Categorization
MSH.1[1]	Field Separator		IG Fixed Data
MSH.2[1]	Encoding Characters	^~\&#	IG Fixed Data
MSH.3[1]	Sending Application		
MSH.3[1].1	Namespace ID	NIST Test Lab APP	Configurable Data
MSH.3[1].2	Universal ID	2.16.840.1.113883.3.72.5.20	Configurable Data
MSH.3[1].3	Universal ID Type	ISO	IG Fixed Data
MSH.4[1]	Sending Facility		
MSH.4[1].1	Namespace ID	NIST Lab Facility	Configurable Data
MSH.4[1].2	Universal ID	2.16.840.1.113883.3.72.5.21	Configurable Data
MSH.4[1].3	Universal ID Type	ISO	IG Fixed Data
MSH.6[1]	Receiving Facility		
MSH.6[1].1	Namespace ID	NIST EHR Facility	Configurable Data
MSH.6[1].2	Universal ID	2.16.840.1.113883.3.72.5.23	Configurable Data
MSH.6[1].3	Universal ID Type	ISO	IG Fixed Data
MSH.7[1]	Date/Time Of Message		
MSH.7[1].1	Time	20130421113601-0700	System Generated
MSH.9[1]	Message Type		
MSH.9[1].1	Message Code	MFN	IG Fixed Data
MSH.9[1].2	Trigger Event	M04	IG Fixed Data
MSH.9[1].3	Message Structure	MFN_M04	IG Fixed Data
MSH.10[1]	Message Control ID	EDOS_2.4_3.1-M04_GU	System Generated
MSH.11[1]	Processing ID		
MSH.11[1].1	Processing ID	T	Changeable Data
MSH.12[1]	VersionID		
MSH.12[1].1	Version ID	2.5.1	IG Fixed Data
MSH.15[1]	Accept Acknowledgment Type	AL	Changeable Data
MSH.16[1]	Application Acknowledgment Type	NE	Changeable Data
MSH.21[1]	Message Profile Identifier		
MSH.21[1].1	Entity Identifier	EDOS_Common_Component	Test Case Fixed Data
MSH.21[1].2	Namespace ID	EDOS Base Profile	Changeable Data
MSH.21[1].3	Universal ID	2.16.840.1.113883.9.67	Test Case Fixed Data
MSH.21[1].4	Universal ID Type	ISO	IG Fixed Data
MSH.21[2]	Message Profile Identifier		
MSH.21[2].1	Entity Identifier	EDOS_GU_Component	Test Case Fixed Data
MSH.21[2].2	Namespace ID	EDOS GU Profile	Changeable Data
MSH.21[2].3	Universal ID	2.16.840.1.113883.9.68	Test Case Fixed Data
MSH.21[2].4	Universal ID Type	ISO	IG Fixed Data

## MFI :

Location	Data Element	Data	Categorization
<b>MFI.1[1]</b>	<b>Master File Identifier</b>		
MFI.1[1].1	Identifier	CDM	IG Fixed Data
MFI.1[1].2	Text	Charge description master file	Changeable Data
MFI.1[1].3	Name of Coding System	HL70175	IG Fixed Data
MFI.1[1].7	Coding System Version ID	2.5.1	Test Case Fixed Data
<b>MFI.3[1]</b>	<b>File-Level Event Code</b>	UPD	Test Case Fixed Data
<b>MFI.6[1]</b>	<b>Response Level Code</b>	NE	IG Fixed Data

## MFE :

Location	Data Element	Data	Categorization
<b>MFE.1[1]</b>	<b>Record-Level Event Code</b>	MDC	Test Case Fixed Data
<b>MFE.3[1]</b>	<b>Effective Date/Time</b>		
MFE.3[1].1	Time	20131219145310	System Generated
<b>MFE.4[1]</b>	<b>Primary Key Value - MFE</b>		
MFE.4[1].1	Identifier	1300	Changeable Data
MFE.4[1].2	Text	Arbovirus IgG and IgM Panel (DNG, WNV) in Serum	Changeable Data
MFE.4[1].3	Name of Coding System	99USI	Changeable Data
MFE.4[1].7	Coding System Version ID	20130421	Changeable Data
<b>MFE.5[1]</b>	<b>Primary Key Value Type</b>	CWE	IG Fixed Data

## CDM :

Location	Data Element	Data	Categorization
<b>CDM.1[1]</b>	<b>Primary Key Value - CDM</b>		
CDM.1[1].1	Identifier	1300	Test Case Fixed Data
CDM.1[1].2	Text	Arbovirus IgG and IgM Panel (DNG, WNV) in Serum	Changeable Data
CDM.1[1].3	Name of Coding System	99USI	Changeable Data
CDM.1[1].7	Coding System Version ID	20130421	Changeable Data
<b>CDM.3[1]</b>	<b>Charge Description Short</b>	Arbovirus IgG and IgM Panel	Changeable Data
<b>CDM.7[1]</b>	<b>Procedure Code</b>		
CDM.7[1].1	Identifier	86788	Changeable Data
CDM.7[1].2	Text	West Nile virus antibody, IgM	Changeable Data
CDM.7[1].3	Name of Coding System	C4	IG Fixed Data
CDM.7[1].7	Coding System Version ID	19.2	Changeable Data
<b>CDM.7[2]</b>	<b>Procedure Code</b>		
CDM.7[2].1	Identifier	86789	Changeable Data
CDM.7[2].2	Text	West Nile virus antibody, IgG	Changeable Data
CDM.7[2].3	Name of Coding System	C4	IG Fixed Data
CDM.7[2].7	Coding System Version ID	19.2	Changeable Data
<b>CDM.7[3]</b>	<b>Procedure Code</b>		
CDM.7[3].1	Identifier	86790	Changeable Data
CDM.7[3].2	Text	DengueFever antibody, IgG and IgM	Changeable Data
CDM.7[3].3	Name of Coding System	C4	IG Fixed Data
CDM.7[3].7	Coding System Version ID	19.2	Changeable Data
<b>CDM.7[4]</b>	<b>Procedure Code</b>		
CDM.7[4].1	Identifier	86790	Changeable Data
CDM.7[4].2	Text	DengueFever antibody, IgG and IgM	Changeable Data
CDM.7[4].3	Name of Coding System	C4	IG Fixed Data
CDM.7[4].7	Coding System Version ID	19.2	Changeable Data

## MFE :

Location	Data Element	Data	Categorization
<b>MFE.1[1]</b>	<b>Record-Level Event Code</b>	MAD	Test Case Fixed Data
<b>MFE.3[1]</b>	<b>Effective Date/Time</b>		
MFE.3[1].1	Time	20131219145310	System Generated
<b>MFE.4[1]</b>	<b>Primary Key Value - MFE</b>		
MFE.4[1].1	Identifier	400	Changeable Data
MFE.4[1].2	Text	Lipid Panel	Changeable Data
MFE.4[1].3	Name of Coding System	99USI	Changeable Data
MFE.4[1].7	Coding System Version ID	20130421	Changeable Data
<b>MFE.5[1]</b>	<b>Primary Key Value Type</b>	CWE	IG Fixed Data

## CDM :

Location	Data Element	Data	Categorization
<b>CDM.1[1]</b>	<b>Primary Key Value - CDM</b>		
CDM.1[1].1	Identifier	400	Test Case Fixed Data
CDM.1[1].2	Text	Lipid Panel	Changeable Data
CDM.1[1].3	Name of Coding System	99USI	Changeable Data
CDM.1[1].7	Coding System Version ID	20130421	Changeable Data
<b>CDM.3[1]</b>	<b>Charge Description Short</b>	Lipid Panel	Changeable Data
<b>CDM.7[1]</b>	<b>Procedure Code</b>		
CDM.7[1].1	Identifier	80061	Changeable Data
CDM.7[1].2	Text	Lipid panel This panel must include the following: Cholesterol, serum, total (82465) Lipoprotein, direct measurement, high density cholesterol (HDL cholesterol) (83718) Triglycerides (84478)	Changeable Data
CDM.7[1].3	Name of Coding System	C4	IG Fixed Data
CDM.7[1].7	Coding System Version ID	19.2	Changeable Data

## MFE :

Location	Data Element	Data	Categorization
<b>MFE.1[1]</b>	<b>Record-Level Event Code</b>	MAD	Test Case Fixed Data
<b>MFE.3[1]</b>	<b>Effective Date/Time</b>		
MFE.3[1].1	Time	20131219145310	System Generated
<b>MFE.4[1]</b>	<b>Primary Key Value - MFE</b>		
MFE.4[1].1	Identifier	410	Changeable Data
MFE.4[1].2	Text	Low density lipoprotein cholesterol, serum (LDL) - measured	Changeable Data
MFE.4[1].3	Name of Coding System	99USI	Changeable Data
MFE.4[1].7	Coding System Version ID	20130421	Changeable Data
<b>MFE.5[1]</b>	<b>Primary Key Value Type</b>	CWE	IG Fixed Data

## CDM :

Location	Data Element	Data	Categorization
<b>CDM.1[1]</b>	<b>Primary Key Value - CDM</b>		
CDM.1[1].1	Identifier	410	Test Case Fixed Data
CDM.1[1].2	Text	Low density lipoprotein cholesterol, serum (LDL) - measured	Changeable Data
CDM.1[1].3	Name of Coding System	99USI	Changeable Data
CDM.1[1].7	Coding System Version ID	20130421	Changeable Data
<b>CDM.3[1]</b>	<b>Charge Description Short</b>	Low density lipoprotein cholesterol, serum (LDL)	Changeable Data
<b>CDM.7[1]</b>	<b>Procedure Code</b>		
CDM.7[1].1	Identifier	87721	Changeable Data
CDM.7[1].2	Text	lipoprotein, direct measurement; low density cholesterol (ldl cholesterol)	Changeable Data
CDM.7[1].3	Name of Coding System	C4	IG Fixed Data
CDM.7[1].7	Coding System Version ID	19.2	Changeable Data

## MFE :

Location	Data Element	Data	Categorization
MFE.1[1]	Record-Level Event Code	MUP	Test Case Fixed Data
MFE.3[1]	Effective Date/Time		
MFE.3[1].1	Time	20131219145310	System Generated
MFE.4[1]	Primary Key Value - MFE		
MFE.4[1].1	Identifier	600	Changeable Data
MFE.4[1].2	Text	Prostate Biopsy Pathology Report	Changeable Data
MFE.4[1].3	Name of Coding System	99USI	Changeable Data
MFE.4[1].7	Coding System Version ID	20130421	Changeable Data
MFE.5[1]	Primary Key Value Type	CWE	IG Fixed Data

## CDM :

Location	Data Element	Data	Categorization
CDM.1[1]	Primary Key Value - CDM		
CDM.1[1].1	Identifier	600	Test Case Fixed Data
CDM.1[1].2	Text	Prostate Biopsy Pathology Report	Changeable Data
CDM.1[1].3	Name of Coding System	99USI	Changeable Data
CDM.1[1].7	Coding System Version ID	20130421	Changeable Data
CDM.3[1]	Charge Description Short	Prostate Biopsy Pathology Report	Changeable Data
CDM.7[1]	Procedure Code		
CDM.7[1].1	Identifier	G0416	Changeable Data
CDM.7[1].2	Text	Surgical pathology, gross and micro exam for prostate needle saturation biopsy sampling 1-20 specimens	Changeable Data
CDM.7[1].3	Name of Coding System	C4	IG Fixed Data
CDM.7[1].7	Coding System Version ID	19.2	Changeable Data

## MFE :

Location	Data Element	Data	Categorization
MFE.1[1]	Record-Level Event Code	MDC	Test Case Fixed Data
MFE.3[1]	Effective Date/Time		
MFE.3[1].1	Time	20131219145310	System Generated
MFE.4[1]	Primary Key Value - MFE		
MFE.4[1].1	Identifier	1305	Changeable Data
MFE.4[1].2	Text	SLE IgG Titer Serum	Changeable Data
MFE.4[1].3	Name of Coding System	99USI	Changeable Data
MFE.4[1].7	Coding System Version ID	20130421	Changeable Data
MFE.5[1]	Primary Key Value Type	CWE	IG Fixed Data

CDM :

Location	Data Element	Data	Categorization
<b>CDM.1[1]</b>	<b>Primary Key Value - CDM</b>		
CDM.1[1].1	Identifier	1305	Changeable Data
CDM.1[1].2	Text	SLE IgG Titer Serum	Changeable Data
CDM.1[1].3	Name of Coding System	99USI	Changeable Data
CDM.1[1].7	Coding System Version ID	20130421	Changeable Data
<b>CDM.3[1]</b>	<b>Charge Description Short</b>	SLE IgG	Changeable Data
<b>CDM.7[1]</b>	<b>Procedure Code</b>		
CDM.7[1].1	Identifier	86653	Changeable Data
CDM.7[1].2	Text	St. Louis encephalitis antibody, IgG and IgM	Changeable Data
CDM.7[1].3	Name of Coding System	C4	IG Fixed Data
CDM.7[1].7	Coding System Version ID	19.2	Changeable Data