

Master File Identification

Element name	Data
--------------	------

[1]

Master File Entry

Element name	Data
--------------	------

General Segment

Element name	Data
Observations Required to Interpret the Observation	Fasting status [Presence] - reported
Observations Required to Interpret the Observation	Race
Interpretation of Observations	Test used to measure blood sugar, electrolytes and fluid balance, kidney and liver function.
Patient Preparation	Patient fasting required for 12 hours.
Service/Test/Observation Performance Schedule	Daily
Other Names	CMP
Expected Turn-Around Time	1 day

Observation batteries

Element name	Data
Test/Observations Included Within an Ordered Test Battery [1]	Serum Glucose
Test/Observations Included Within an Ordered Test Battery [2]	Blood Urea Nitrogen (BUN)
Test/Observations Included Within an Ordered Test Battery [3]	Creatinine
Test/Observations Included Within an Ordered Test Battery [4]	BUN/Creatinine Ratio
Test/Observations Included Within an Ordered Test Battery [5]	GFR, calculated
Test/Observations Included Within an Ordered Test Battery [6]	Calcium
Test/Observations Included Within an Ordered Test Battery [7]	Total protein, serum
Test/Observations Included Within an Ordered Test Battery [8]	Albumin
Test/Observations Included Within an Ordered Test Battery [9]	Globulin
Test/Observations Included Within an Ordered Test Battery [10]	Albumin/globulin ratio
Test/Observations Included Within an Ordered Test Battery [11]	Total bilirubin, serum
Test/Observations Included Within an Ordered Test Battery [12]	Alkaline phosphatase (ALP)
Test/Observations Included Within an Ordered Test Battery [13]	Alanine aminotransferase (ALT)
Test/Observations Included Within an Ordered Test Battery [14]	Aspartate aminotransferase (ASP)
Test/Observations Included Within an Ordered Test Battery [15]	Sodium, serum

Test/Observations Included Within an Ordered Test Battery [16]	Potassium, serum
Test/Observations Included Within an Ordered Test Battery [17]	Chloride, serum
Test/Observations Included Within an Ordered Test Battery [18]	Carbon dioxide, serum
Test/Observations Included Within an Ordered Test Battery [19]	Anion gap

Observations that Require Specimens

Element name	Data
Container Description	Gold Serum Separator tube
Container Description	Red, No Additive tube
Container Volume	5.0
Container Volume	5.0
Container Units	mililiter
Container Units	mililiter
Specimen Handling Code	Refrigerated temperature

Observations that Require Specimens

Element name	Data
Container Description	Green Lithium Heparin tube
Container Volume	3.0
Container Units	mililiter
Specimen Handling Code	Refrigerated temperature