NIST HL7 V2 Validation Tool Test Data Electronic Directory Of Service (eDOS)

Version 1.0

13th March, 2015

eDOS Test Case MFK_1_1.2-MFK_M10_GU

 $EDOS_1.0_2.1\text{-}M10_GU$

Message Content Data Sheet

Test Case Information

MFK_1_1.2-GU - EI	OOS_1.0_2.1-M10_GU
Test Case ID	MFK_1_1.2-MFK_M10_GU

MSH:

Location	Data Element	Data	Categorization
MSH.1	Field Separator		IG Fixed Data
MSH.2	Encoding Characters	^~\&#</td><td>IG Fixed Data</td></tr><tr><td>MSH.3</td><td>Sending Application</td><td></td><td></td></tr><tr><td>MSH.3.1</td><td>Namespace ID</td><td>NIST EHR</td><td>Configurable Data</td></tr><tr><td>MSH.3.2</td><td>Universal ID</td><td>2.16.840.1.113883.3.72.5.22</td><td>Configurable Data</td></tr><tr><td>MSH.3.3</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr><tr><td>MSH.4</td><td>Sending Facility</td><td></td><td></td></tr><tr><td>MSH.4.1</td><td>Namespace ID</td><td>NIST EHR Facility</td><td>Configurable Data</td></tr><tr><td>MSH.4.2</td><td>Universal ID</td><td>2.16.840.1.113883.3.72.5.23</td><td>Configurable Data</td></tr><tr><td>MSH.4.3</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr><tr><td>MSH.6</td><td>Receiving Facility</td><td></td><td></td></tr><tr><td>MSH.6.1</td><td>Namespace ID</td><td>NIST Test Lab APP</td><td>Configurable Data</td></tr><tr><td>MSH.6.2</td><td>Universal ID</td><td>2.16.840.1.113883.3.72.5.20</td><td>Configurable Data</td></tr><tr><td>MSH.6.3</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr><tr><td>MSH.7</td><td>Date/Time Of Message</td><td></td><td></td></tr><tr><td>MSH.7.1</td><td>Time</td><td>201302201110015-0800</td><td>System Generated</td></tr><tr><td>MSH.9</td><td>Message Type</td><td></td><td></td></tr><tr><td>MSH.9.1</td><td>Message Code</td><td>MFK</td><td>IG Fixed Data</td></tr><tr><td>MSH.9.2</td><td>Trigger Event</td><td>M10</td><td>IG Fixed Data</td></tr><tr><td>MSH.9.3</td><td>Message Structure</td><td>MFK_M01</td><td>IG Fixed Data</td></tr><tr><td>MSH.10</td><td>Message Control ID</td><td>MFK_1_1.2-GU</td><td>System Generated</td></tr><tr><td>MSH.11</td><td>Processing ID</td><td></td><td></td></tr><tr><td>MSH.11.1</td><td>Processing ID</td><td>Т</td><td>Changeable Data</td></tr><tr><td>MSH.12</td><td>VersionID</td><td></td><td></td></tr><tr><td>MSH.12.1</td><td>Version ID</td><td>2.5.1</td><td>IG Fixed Data</td></tr><tr><td>MSH.15</td><td>Accept Acknowledgment Typ</td><td>e NE</td><td>Changeable Data</td></tr></tbody></table>	

NIST HL7 V2 eDOS Validation Tool Test Data

MSH.16	Application Acknowledgment Type	NE	Changeable Data
MSH.21[1]	Message Profile Identifier		
MSH.21[1].1	Entity Identifier	EDOS_GU_RESPONSE_PROFII	Eest Case Fixed Data
MSH.21[1].2	Namespace ID	EDOS Base Profile	Changeable Data
MSH.21[1].3	Universal ID	2.16.840.1.113883.9.75	Test Case Fixed Data
MSH.21[1].4	Universal ID Type	ISO	IG Fixed Data

MSA:

Location	Data Element	Data	Categorization
MSA.1	Acknowledgement Code	CA	Test Case Fixed Data
MSA.2	Message Control ID	EDOS_1.0_2.1-M10_GU	Test Case Fixed Data

MFI:

Location	Data Element	Data	Categorization
MFI.1	Master File Identifier		
MFI.1.1	Identifier	OMC	IG Fixed Data
MFI.1.2	Text	Observation batteries master file	Changeable Data
MFI.1.3	Name of Coding System	HL70175	IG Fixed Data
MFI.1.7	Coding System Version ID	2.5.1	Test Case Fixed Data
MFI.3	File-Level Event Code	REP	Test Case Fixed Data
MFI.6	Response Level Code	NE	IG Fixed Data

Test Data Specification

Master File Identification

Element name	Data
Hile-Level Event Code	Replace current version of this master file with the version contained in this message

UI 7 Magaga
HL7 Message
MSH ^~\&# NIST EHR^2.16.840.1.113883.3.72.5.22^ISO NIST EHR Facility^2.16.840.1.113883.3.72.5.23^ISO NIST Test Lab APP ^2.16.840.1.113883.3.72.5.20^ISO 201302201110015-0800 MFK^M10^MFK_M01 MFK_1_1.2-GU T 2.5.1 NE NE EDOS_GU_RESPONSE _PROFILE^EDOS Base Profile^2.16.840.1.113883.9.75^ISO</td></tr><tr><td>MSA CA EDOS_1.0_2.1-M10_GU</td></tr><tr><td>MFI OMC^Observation batteries master file^HL70175^^^2.5.1 REP NE</td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr><tr><td></td></tr></tbody></table>