Element	Data
Identifier	OMC
Event code	UPD

Panel/Test battery Information - GHP Record information

Element	Data
Event code	MAC
Effective date and time	20131219145310
Test/Panel Name	GHP
Test/Panel identifier (code system)	800 (99USL)

General information-

Element	Data	
Specimen required ?	Y	
Preferred Report Name	General Health Profile	
Is the test/panel orderable ?	Y	
Nature	S	
Interpretation	This blood test is used to determine general health status and to screen for and monitor a variety of disorders. This profile includes a complete metabolic profile, comprehensive CBC, Urinalysis and total Thyrotropin (T4).	
Patient preparation	Patient fasting required for 12 hours.	
Factors that may affect the test	Insufficient specimen, Gross hemolysis, Improper labeling.	
Schedule	Daily	
Is the test exclusive?	N	
Diagnostic Service Sector	LAB	
Expected turn around time	1 day	

- Panel/Test battery Details -

Element	Data
Test included in panel/test battery	CMP
Test included in panel/test battery (code & code system)	100 (99USL)
Test included in panel/test battery	CBC_diff
Test included in panel/test battery (code & code system)	200 (99USL)
Test included in panel/test battery	TSH
Test included in panel/test battery (code & code system)	700 (99USL)
Test included in panel/test battery	Comprehensive Urinalysis
Test included in panel/test battery (code & code system)	300 (99USL)

Specimen information details-

~F		
Element	Data	
Sequence ID	1	
Container description	Gold Serum Separator tube	
Container description	Red, No Additive tube	
Container volume	5.0	
Container volume	5.0	
Container units	milliliter	
Container units (code & code system)	mL (UCUM)	
Container units	milliliter	
Container units (code & code system)	mL (UCUM)	
Specimen (code & code system)	119364003 (SCT)	
Normal collection volume	1	
Normal collection volume unit	milliliter	
Normal collection volume unit (code & code system)	mL (UCUM)	

Normal collection volume	0.5
Normal collection volume unit	milliliter
Normal collection volume unit (code & code system)	mL (UCUM)
Specimen Requirements	Protect from light. Allow serum tube to clot completely at room temperature. Separate serum or plasma from cells ASAP or within 30 minutes of collection.
Specimen Handling Code	Refrigerated temperature
Specimen Handling Code (code & code system)	REF (HL70376)
Specimen Preference	P

-Specimen information details-

Element	Data	
Sequence ID	1.2	
Container description	Lavender Top (EDTA) tube	
Container description	Pink Top (K2EDTA) tube	
Container volume	3.0	
Container volume	3.0	
Container units	milliliters	
Container units (code & code system)	mL (UCUM)	
Container units	milliliters	
Container units (code & code system)	mL (UCUM)	
Specimen (code & code system)	119297000 (SCT)	
Specimen (code & code system)	EDTK (HL70371)	
Normal collection volume	3	
Normal collection volume unit	milliliters	
Normal collection volume unit (code & code system)	mL (UCUM)	
Normal collection volume	0.5	
Normal collection volume unit	milliliters	
Normal collection volume unit (code & code system)	mL (UCUM)	
Specimen Requirements	Refrigeration is required if specimen is not brought immediat to laboratory. Two blood smears should be prepared if sam is not delivered to the laboratory within 4 hrs. Sample should be analyzed within 6 hours at room temperature and 24 hrs when stored at 4 degrees C.	
Specimen Handling Code	Critical refrigerated	
Specimen Handling Code (code & code system)	CREF (HL70376)	
Specimen Preference	P	

Specimen information details-

Element	Data
Sequence ID	1.3
Container description	Sterile, plastic, leak proof container
Container volume	4
Container units	fluid ounce (US)
Container units (code & code system)	[foz_us] (UCUM)
Specimen (code & code system)	122575003 (SCT)
Normal collection volume	20
Normal collection volume unit	milliliter
Normal collection volume unit (code & code system)	mL (UCUM)
Normal collection volume	4
Normal collection volume unit	milliliter
Normal collection volume unit (code & code system)	mL (UCUM)
Specimen Requirements	Keep refrigerated
Specimen Handling Code	Refrigerated temperature

Specimen Handling Code (code & code system)	REF (HL70376)	
Specimen Preference	P	