

Description

The initial Laboratory Test Compendium built by Century Hospital Clinical Laboratory consists of a single Prothrombin Time panel ("PT+INR"). This compendium is delivered electronically to Dr. Nicholas Radon's practice EHR, as agreed to between Laboratory and Dr. Radon's Office and EHR-S Vendor. The EHR-S will integrate the eDOS into its test directory and use it to allow Dr. Radon to place orders electronically to Century Hospital's Clinical Laboratory. The initial laboratory test compendium is composed of four messages. Century Hospital's LIS transmit the first message which provides information about the PT and INR observations.

Comments

This is the initial "smoke testing" to reveal simple failures severe enough to preclude further testing. Initial upload of 2 individual observations for a new compendium.
Since this is a replace message, all data in this message will replace information from a previous eDOS M08 message.

PreCondition

No Pre-Condition.

PostCondition

eDOS M08 data elements are incorporated appropriately into the EHR-S.

TestObjectives

- Demonstrate capability to support a minimal M08 initial load message.
- Demonstrate capability to support a single occurrence of all required ("R") elements in the MSH segment.
- Demonstrate capability to support a single occurrence of all required ("R") elements in the MFI segment.
- Demonstrate capability to support a single occurrence of all required ("R") elements in the MFE segment.
- Demonstrate capability to support a single occurrence of all required ("R") elements in the OM1 segment.

Test Case Information

EDOS_0.0_1.1-M08_GU - M08_New_Smoke_Test

Test Case ID	EDOS_0.0_1.1-M08_GU
--------------	---------------------

MSH

Location	Data Element	Data	Categorization
MSH.1	Field Separator		IG Fixed Data
MSH.2	Encoding Characters	~\&	IG Fixed Data
MSH.4	Sending Facility		
MSH.4.2	Universal ID	2.16.840.1.113883.3.72.5.21	Configurable Data
MSH.4.3	Universal ID Type	ISO	IG Fixed Data
MSH.7	Date/Time Of Message		
MSH.7.1	Time	20130421113601-0700	System Generated
MSH.9	Message Type		
MSH.9.1	Message Code	MFN	IG Fixed Data
MSH.9.2	Trigger Event	M08	IG Fixed Data
MSH.9.3	Message Structure	MFN_M08	IG Fixed Data
MSH.10	Message Control ID	EDOS_0.0_1.1-M08_GU	System Generated
MSH.11	Processing ID		
MSH.11.1	Processing ID	D	Changeable Data
MSH.12	VersionID		
MSH.12.1	Version ID	2.5.1	IG Fixed Data
MSH.21	Message Profile Identifier		
MSH.21.1	Entity Identifier	EDOS_GU_Profile	Test Case Fixed Data
MSH.21.3	Universal ID	2.16.840.1.113883.9.70	Test Case Fixed Data
MSH.21.4	Universal ID Type	ISO	IG Fixed Data

MFI

Location	Data Element	Data	Categorization
MFI.1	Master File Identifier		
MFI.1.1	Identifier	OMM	IG Fixed Data
MFI.1.3	Name of Coding System	HL70175	IG Fixed Data
MFI.3	File-Level Event Code	REP	Test Case Fixed Data
MFI.6	Response Level Code	NE	IG Fixed Data

MFE

Location	Data Element	Data	Categorization
MFE.1	Record-Level Event Code	MAD	Test Case Fixed Data
MFE.3	Effective Date/Time		
MFE.3.1	Time	20131219145310	System Generated
MFE.4	Primary Key Value - MFE		
MFE.4.1	Identifier	11	Changeable Data
MFE.4.2	Text	Prothrombin Time, PT	Changeable Data
MFE.4.3	Name of Coding System	99USL	Changeable Data
MFE.5	Primary Key Value Type	CWE	IG Fixed Data

OM1

Location	Data Element	Data	Categorization
OM1.1	Sequence Number - Test/Observation Master File	1	IG Fixed Data
OM1.2	Producer's Service/Test/Observation ID		
OM1.2.1	Identifier	11	Changeable Data
OM1.2.2	Text	Prothrombin Time, PT	Changeable Data
OM1.2.3	Name of Coding System	99USL	Changeable Data
OM1.4	Specimen Required	Y	Changeable Data

Location	Producer ID Data Element	Data	Categorization
OM1.5.2	Text	Century Hospital Clinical Laboratory	Changeable Data
OM1.11	Preferred Long Name for the Observation	Prothrombin Time	Changeable Data
OM1.12	Orderability	N	Test Case Fixed Data
OM1.18	Nature of Service/Test/Observation	A	Test Case Fixed Data

MFE

Location	Data Element	Data	Categorization
MFE.1	Record-Level Event Code	MAD	Test Case Fixed Data
MFE.3	Effective Date/Time		
MFE.3.1	Time	20131219145310	System Generated
MFE.4	Primary Key Value - MFE		
MFE.4.1	Identifier	12	Changeable Data
MFE.4.2	Text	INR	Changeable Data
MFE.4.3	Name of Coding System	99USL	Changeable Data
MFE.5	Primary Key Value Type	CWE	IG Fixed Data

OM1

Location	Data Element	Data	Categorization
OM1.1	Sequence Number - Test/Observation Master File	2	IG Fixed Data
OM1.2	Producer's Service/Test/Observation ID		
OM1.2.1	Identifier	12	Changeable Data
OM1.2.2	Text	INR	Changeable Data
OM1.2.3	Name of Coding System	99USL	Changeable Data
OM1.4	Specimen Required	N	Changeable Data
OM1.5	Producer ID		
OM1.5.2	Text	Century Hospital Clinical Laboratory	Changeable Data
OM1.11	Preferred Long Name for the Observation	International Normalized Ratio	Changeable Data
OM1.12	Orderability	N	Test Case Fixed Data
OM1.18	Nature of Service/Test/Observation	C	Test Case Fixed Data

Element	Data
Identifier	OMM
Event code	REP

Test Details - Prothrombin Time, PT
Record information

Element	Data
Event code	MAD
Effective date and time	20131219145310
Test/Panel Name	Prothrombin Time, PT
Test/Panel identifier (code system)	11 (99USL)

General information

Element	Data
Specimen required ?	Y
Preferred Long Name	Prothrombin Time
Is the test/panel orderable ?	N
Nature	A

Test Details - INR
Record information

Element	Data
Event code	MAD
Effective date and time	20131219145310
Test/Panel Name	INR
Test/Panel identifier (code system)	12 (99USL)

General information

Element	Data
Specimen required ?	N
Preferred Long Name	International Normalized Ratio
Is the test/panel orderable ?	N
Nature	C

MSH|^~\&|^2.16.840.1.113883.3.72.5.21^ISO|||20130421113601-0700||MFN^M08^MFN_M08|EDOS_0.0_1.1-M08_GU|D|2.5.1|||||||E
DOS_GU_Profile^^2.16.840.1.113883.9.70^ISO

MFI|OMM^^HL70175||REP|||NE

MFE|MAD||20131219145310|11^Prothrombin Time, PT^99USL|CWE

OM1|1|11^Prothrombin Time, PT^99USL||Y|^Century Hospital Clinical Laboratory|||||Prothrombin Time|N|||||A

MFE|MAD||20131219145310|12^INR^99USL|CWE

OM1|2|12^INR^99USL||N|^Century Hospital Clinical Laboratory|||||International Normalized Ratio|N|||||C