

### Description

Mr. William A. Jones is a 51 year old white male who presented with a complaint of diarrhea times 3 days. Dr. Nicholas Radon ordered a stool culture test to be performed. A stool specimen for the microbiology test was collected from the patient and sent to the clinical lab at Century Hospital, 2070 Test Park, Los Angeles, CA, 90067, for processing. The stool specimen was analyzed and preliminary result are sent 24 hours later, identifying E. coli O157:H7, Salmonella and Shigella flexneri.

Approximately 48 hours later the Final results, including susceptibility testing of Shigella flexneri are available. The final result report was generated by the LIS and transmitted to the patient's record in the ambulatory EHR used in Dr. Radon's office practice.

Five hours later the laboratorian realizes that the result in the susceptibility suite for Shigella Flexneri were missing results for additional antibiotic.

### Comments

This is a variation of a follow up report to LRI\_4.1\_2.1-GU\_FRU

This test case is evaluating correct use of the status codes in OBR-25 and OBX-11.

Special case test aspects include ensuring that the system can correctly populate the data elements for linking parent/child results in accordance with the requirements specified in the implementation guide for the FRN and FRU profiles.

### PreCondition

LRI\_4.0\_1.1-GU AND LRI\_4.1\_2.1-GU\_FRU have been sent previously.

### PostCondition

No Post Condition.

### TestObjectives

- Determine if the system can correctly create a valid Parent-Child message for an appended Stool Culture/Susceptibility microbiology test report.

### Notes to Testers

This is a variation of a follow up report to LRI\_4.1\_2.1-GU\_FRU

Check that the parent-child linkage is correctly implemented:

The second OBR has OBR-26 populated with the result that spawned the respective susceptibility testing, e.g. the second group, first sequence of the Stool culture ( Salmonella I, group O:4)

The third OBR has OBR-26 populated with the result that spawned the respective susceptibility testing, e.g. the third group, first sequence Stool culture (Shigella flexneri).

In both cases the OBR-29 references the Placer and Filler Order Number of the Stool culture OBR.

In the third order group additional results are now available.

Test Case Information

**LRI\_4.1\_4.1-GU\_FRU - Culture-Escherichia coli, Salmonella, Shigella - Parent/Child Susceptibility - Correction**

Test Case ID LRI\_4.1\_4.1-GU\_FRU

MSH

Location	Data Element	Data	Categorization
<b>MSH.1</b>	<b>Field Separator</b>		IG Fixed Data
<b>MSH.2</b>	<b>Encoding Characters</b>	~\&	IG Fixed Data
<b>MSH.3</b>	<b>Sending Application</b>		
MSH.3.2	Universal ID	2.16.840.1.113883.3.72.5.20	Configurable Data
MSH.3.3	Universal ID Type	ISO	IG Fixed Data
<b>MSH.4</b>	<b>Sending Facility</b>		
MSH.4.2	Universal ID	2.16.840.1.113883.3.72.5.21	Configurable Data
MSH.4.3	Universal ID Type	ISO	IG Fixed Data
<b>MSH.6</b>	<b>Receiving Facility</b>		
MSH.6.2	Universal ID	2.16.840.1.113883.3.72.5.23	Configurable Data
MSH.6.3	Universal ID Type	ISO	IG Fixed Data
<b>MSH.7</b>	<b>Date/Time Of Message</b>		
MSH.7.1	Time	20150926140551	System Generated
<b>MSH.9</b>	<b>Message Type</b>		
MSH.9.1	Message Code	ORU	IG Fixed Data
MSH.9.2	Event Type	R01	IG Fixed Data
MSH.9.3	Message Structure	ORU_R01	IG Fixed Data
<b>MSH.10</b>	<b>Message Control ID</b>	LRI_4.1_4.1-GU_FRU	System Generated
<b>MSH.11</b>	<b>Processing ID</b>		
MSH.11.1	Processing ID	D	Changeable Data
<b>MSH.12</b>	<b>VersionID</b>		
MSH.12.1	Version ID	2.5.1	IG Fixed Data
<b>MSH.15</b>	<b>Accept Acknowledgment Type</b>	AL	IG Fixed Data
<b>MSH.16</b>	<b>Application Acknowledgment Type</b>	AL	IG Fixed Data
<b>MSH.21</b>	<b>Message Profile Identifier</b>		
MSH.21.1	Entity Identifier	LRI_Common_Component	IG Fixed Data
MSH.21.3	Universal ID	2.16.840.1.113883.9.16	IG Fixed Data
MSH.21.4	Universal ID Type	ISO	IG Fixed Data
<b>MSH.21[2]</b>	<b>Message Profile Identifier</b>		
MSH.21[2].1	Entity Identifier	LRI_GU_Component	IG Fixed Data
MSH.21[2].3	Universal ID	2.16.840.1.113883.9.12	IG Fixed Data
MSH.21[2].4	Universal ID Type	ISO	IG Fixed Data
<b>MSH.21[3]</b>	<b>Message Profile Identifier</b>		
MSH.21[3].1	Entity Identifier	LRI_FRU_Component	IG Fixed Data
MSH.21[3].3	Universal ID	2.16.840.1.113883.9.83	IG Fixed Data
MSH.21[3].4	Universal ID Type	ISO	IG Fixed Data

PID

Location	Data Element	Data	Categorization
<b>PID.1</b>	<b>Set ID - PID</b>	1	IG Fixed Data
<b>PID.3</b>	<b>Patient Identifier List</b>		
PID.3.1	ID Number	PATID1234	Configurable Data
PID.3.4	Assigning Authority		
PID.3.4.2	Universal ID	2.16.840.1.113883.3.72.5.30.2	Configurable Data
PID.3.4.3	Universal ID Type	ISO	IG Fixed Data
PID.3.5	Identifier Type Code	MR	Configurable Data
<b>PID.5</b>	<b>Patient Name</b>		
PID.5.1	Family Name		
PID.5.1.1	Surname	Jones	Changeable Data

<b>PID.5.2</b>	<b>Location</b>	<b>Data Element</b>	<b>William</b>	<b>Data</b>	<b>Changeable Data</b>	<b>Categorization</b>
PID.5.3		Second and Further Given Names or Initials Thereof	A			Changeable Data
PID.5.7		Name Type Code	L			Changeable Data
<b>PID.7</b>		<b>Date/Time of Birth</b>				
PID.7.1		Time	19610615			Changeable Data
<b>PID.8</b>		<b>Administrative Sex</b>	M			Changeable Data
<b>PID.10</b>		<b>Race</b>				
PID.10.1		Identifier	2106-3			Changeable Data
PID.10.2		Text	White			Changeable Data
PID.10.3		Name of the Coding System	HL70005			Changeable Data
<b>PID.18</b>		<b>Patient Account Number</b>				
PID.18.1		ID Number	PATID1234			Configurable Data
PID.18.4		Assigning Authority				
PID.18.4.2		Universal ID	2.16.840.1.113883.3.72.5.30.2			Configurable Data
PID.18.4.3		Universal ID Type	ISO			IG Fixed Data
PID.18.5		Identifier Type Code	AN			Configurable Data

#### ORC

<b>Location</b>	<b>Data Element</b>	<b>Data</b>	<b>Categorization</b>
<b>ORC.1</b>	<b>Order Control</b>	RE	Test Case Fixed Data
<b>ORC.2</b>	<b>Placer Order Number</b>		
ORC.2.1	Entity Identifier	ORD723222-4	Changeable Data
ORC.2.3	Universal ID	2.16.840.1.113883.3.72.5.24	Changeable Data
ORC.2.4	Universal ID Type	ISO	IG Fixed Data
<b>ORC.3</b>	<b>Filler Order Number</b>		
ORC.3.1	Entity Identifier	R-783274-4	Changeable Data
ORC.3.3	Universal ID	2.16.840.1.113883.3.72.5.25	Changeable Data
ORC.3.4	Universal ID Type	ISO	IG Fixed Data
<b>ORC.4</b>	<b>Placer Group Number</b>		
ORC.4.1	Entity Identifier	GORD874211	Changeable Data
ORC.4.3	Universal ID	2.16.840.1.113883.3.72.5.24	Changeable Data
ORC.4.4	Universal ID Type	ISO	IG Fixed Data
<b>ORC.12</b>	<b>Ordering Provider</b>		
ORC.12.1	ID Number	5742200012	Changeable Data
ORC.12.2	Family Name		
ORC.12.2.1	Surname	Radon	Changeable Data
ORC.12.3	Given Name	Nicholas	Changeable Data
ORC.12.9	Assigning Authority		
ORC.12.9.2	Universal ID	2.16.840.1.113883.4.6	Changeable Data
ORC.12.9.3	Universal ID Type	ISO	Changeable Data
ORC.12.10	Name Type Code	L	Changeable Data
ORC.12.13	Identifier Type Code	NPI	Changeable Data

#### OBR

<b>Location</b>	<b>Data Element</b>	<b>Data</b>	<b>Categorization</b>
<b>OBR.1</b>	<b>Set ID - OBR</b>	1	IG Fixed Data
<b>OBR.2</b>	<b>Placer Order Number</b>		
OBR.2.1	Entity Identifier	ORD723222-4	Changeable Data
OBR.2.3	Universal ID	2.16.840.1.113883.3.72.5.24	Changeable Data
OBR.2.4	Universal ID Type	ISO	Changeable Data
<b>OBR.3</b>	<b>Filler Order Number</b>		
OBR.3.1	Entity Identifier	R-783274-4	Changeable Data
OBR.3.3	Universal ID	2.16.840.1.113883.3.72.5.25	Changeable Data
OBR.3.4	Universal ID Type	ISO	Changeable Data
<b>OBR.4</b>	<b>Universal Service Identifier</b>		
OBR.4.1	Identifier	625-4	Test Case Fixed Data

<b>OBR.4.2</b>	<b>Location</b>	<b>Data Element</b>	<b>Bacteria identified in Stool by Culture</b>	<b>Test Case Fixed Data</b>
OBR.4.3		Name of Coding System	LN	Test Case Fixed Data
OBR.4.4		Alternate Identifier	3456543	Changeable Data
OBR.4.5		Alternate Text	CULTURE STOOL	Changeable Data
OBR.4.6		Name of Alternate Coding System	99USL	Changeable Data
OBR.4.7		Coding System Version	2.52	Changeable Data
OBR.4.9		Original Text	Stool Culture	Changeable Data
<b>OBR.7</b>		<b>Observation Date/Time</b>		
OBR.7.1		Time	201509231400	Changeable Data
<b>OBR.16</b>		<b>Ordering Provider</b>		
OBR.16.1		ID Number	5742200012	Changeable Data
OBR.16.2		Family Name		
OBR.16.2.1		Surname	Radon	Changeable Data
OBR.16.3		Given Name	Nicholas	Changeable Data
OBR.16.9		Assigning Authority		
OBR.16.9.2		Universal ID	2.16.840.1.113883.4.6	Changeable Data
OBR.16.9.3		Universal ID Type	ISO	Changeable Data
OBR.16.10		Name Type Code	L	Changeable Data
OBR.16.13		Identifier Type Code	NPI	Changeable Data
<b>OBR.22</b>		<b>Results Rpt/Status Chng - Date/Time</b>		
OBR.22.1		Time	20150926140551	Changeable Data
<b>OBR.25</b>		<b>Result Status</b>	F	Test Case Fixed Data
<b>OBR.28</b>		<b>Result Copies To</b>		
OBR.28.1		ID Number	10092000194	Changeable Data
OBR.28.2		Family Name		
OBR.28.2.1		Surname	Hamlin	Changeable Data
OBR.28.3		Given Name	Pafford	Changeable Data
OBR.28.9		Assigning Authority		
OBR.28.9.2		Universal ID	2.16.840.1.113883.4.6	Changeable Data
OBR.28.9.3		Universal ID Type	ISO	Changeable Data
OBR.28.10		Name Type Code	L	Changeable Data
OBR.28.13		Identifier Type Code	NPI	Changeable Data
<b>OBR.47</b>		<b>Filler Supplemental Service Information</b>		
OBR.47.1		Identifier	MIC	Test Case Fixed Data
OBR.47.2		Text	Observation of type microbiology	Test Case Fixed Data
OBR.47.3		Name of the Coding System	HL70411	IG Fixed Data
OBR.47.7		Coding System Version	2.5.1	IG Fixed Data
<b>OBR.49</b>		<b>Results Handling</b>		
OBR.49.1		Identifier	CC	Test Case Fixed Data
OBR.49.2		Text	Copies Requested	Test Case Fixed Data
OBR.49.3		Name of the Coding System	HL70507	IG Fixed Data

#### OBX

<b>Location</b>	<b>Data Element</b>	<b>Data</b>	<b>Categorization</b>
<b>OBX.1</b>	<b>Set ID - OBX</b>	1	IG Fixed Data
<b>OBX.2</b>	<b>Value Type</b>	CWE	Test Case Fixed Data
<b>OBX.3</b>	<b>Observation Identifier</b>		
OBX.3.1	Identifier	625-4	Test Case Fixed Data
OBX.3.2	Text	Bacteria identified in Stool by Culture	Test Case Fixed Data
OBX.3.3	Name of the Coding System	LN	Test Case Fixed Data
OBX.3.7	Coding System Version	2.52	Changeable Data
OBX.3.9	Original Text	Stool Culture	Test Case Fixed Data
<b>OBX.4</b>	<b>Observation Sub-ID</b>		
OBX.4.2	Group	1	Changeable Data
OBX.4.3	Sequence	1	Changeable Data
OBX.4.4	Identifier	Islt-1	Changeable Data

Location	Data Element	Data	Categorization
OBX.5.1	Identifier	103429008	Changeable Data
OBX.5.2	Text	Enterohemorrhagic Escherichia coli, serotype O157:H7	Changeable Data
OBX.5.3	Name of the Coding System	SCT	Changeable Data
OBX.5.7	Coding System Version	201509USED	Changeable Data
OBX.5.9	Original Text	Shiga toxin producing E. coli O157:H7 isolated	Changeable Data
OBX.8	Abnormal Flags	A	Test Case Fixed Data
OBX.11	Observation Result Status	F	Test Case Fixed Data
OBX.14	Date/Time of the Observation		
OBX.14.1	Time	201509231400	Changeable Data
OBX.19	Date/Time of the Analysis		
OBX.19.1	Time	201509251930	Changeable Data
OBX.23	Performing Organization Name		
OBX.23.1	Organization Name	Century Hospital	Changeable Data
OBX.23.6	Assigning Authority		
OBX.23.6.2	Universal ID	2.16.840.1.113883.4.7	Changeable Data
OBX.23.6.3	Universal ID Type	ISO	IG Fixed Data
OBX.23.7	Identifier Type Code	XX	Changeable Data
OBX.23.10	Organization Identifier	24D9871327	Changeable Data
OBX.24	Performing Organization Address		
OBX.24.1	Street Address		
OBX.24.1.1	Street or Mailing Address	2070 Test Park	Changeable Data
OBX.24.3	City	Los Angeles	Changeable Data
OBX.24.4	State or Province	CA	Changeable Data
OBX.24.5	Zip or Postal Code	90067	Changeable Data
OBX.24.7	Address Type	B	Changeable Data
OBX.25	Performing Organization Medical Director		
OBX.25.1	ID Number	5432178916	Changeable Data
OBX.25.2	Family Name		
OBX.25.2.1	Surname	Knowsalot	Changeable Data
OBX.25.3	Given Name	Phil	Changeable Data
OBX.25.4	Second and Further Given Names or Initials Thereof	J.	Changeable Data
OBX.25.9	Assigning Authority		
OBX.25.9.2	Universal ID	2.16.840.1.113883.4.6	Changeable Data
OBX.25.9.3	Universal ID Type	ISO	IG Fixed Data
OBX.25.10	Name Type Code	L	Changeable Data
OBX.25.13	Identifier Type Code	NPI	Changeable Data
OBX.29	Observation Type	RSLT	Test Case Fixed Data

#### NTE

Location	Data Element	Data	Categorization
NTE.1	Set ID - NTE	1	IG Fixed Data
NTE.3	Comment	Susceptibility testing for E.coli is not performed, because antibiotics should not be used to treat this infection. There is no evidence that treatment with antibiotics is helpful, and taking antibiotics may increase the risk of hemolytic-uremic syndrome (HUS). Antidiarrheal agents like Imodium may also increase that risk. Non-specific supportive therapy, including hydration, is important.	Changeable Data

#### OBX

Location	Data Element	Data	Categorization
OBX.1	Set ID - OBX	2	IG Fixed Data

OBX.2	Location	Value Type	Data Element	CWE	Data	Test Case Fixed Data	Categorization
OBX.3	Observation Identifier						
OBX.3.1	Identifier			625-4		Test Case Fixed Data	
OBX.3.2	Text			Bacteria identified in Stool by Culture		Test Case Fixed Data	
OBX.3.3	Name of the Coding System			LN		Test Case Fixed Data	
OBX.3.7	Coding System Version			2.52		Changeable Data	
OBX.3.9	Original Text			Stool Culture		Test Case Fixed Data	
OBX.4	Observation Sub-ID						
OBX.4.2	Group			2		Changeable Data	
OBX.4.3	Sequence			1		Changeable Data	
OBX.4.4	Identifier			Islt-2		Changeable Data	
OBX.5	Observation Value						
OBX.5.1	Identifier			398567006		Changeable Data	
OBX.5.2	Text			Salmonella I, group O:4		Changeable Data	
OBX.5.3	Name of the Coding System			SCT		Changeable Data	
OBX.5.7	Coding System Version			201509USED		Changeable Data	
OBX.5.9	Original Text			Salmonella I, group O:4 isolated		Changeable Data	
OBX.8	Abnormal Flags			A		Test Case Fixed Data	
OBX.11	Observation Result Status			F		Test Case Fixed Data	
OBX.14	Date/Time of the Observation						
OBX.14.1	Time			201509231400		Changeable Data	
OBX.19	Date/Time of the Analysis						
OBX.19.1	Time			201509251930		Changeable Data	
OBX.23	Performing Organization Name						
OBX.23.1	Organization Name			Century Hospital		Changeable Data	
OBX.23.6	Assigning Authority						
OBX.23.6.2	Universal ID			2.16.840.1.113883.4.7		Changeable Data	
OBX.23.6.3	Universal ID Type			ISO		IG Fixed Data	
OBX.23.7	Identifier Type Code			XX		Changeable Data	
OBX.23.10	Organization Identifier			24D9871327		Changeable Data	
OBX.24	Performing Organization Address						
OBX.24.1	Street Address						
OBX.24.1.1	Street or Mailing Address			2070 Test Park		Changeable Data	
OBX.24.3	City			Los Angeles		Changeable Data	
OBX.24.4	State or Province			CA		Changeable Data	
OBX.24.5	Zip or Postal Code			90067		Changeable Data	
OBX.24.7	Address Type			B		Changeable Data	
OBX.25	Performing Organization Medical Director						
OBX.25.1	ID Number			5432178916		Changeable Data	
OBX.25.2	Family Name						
OBX.25.2.1	Surname			Knowsalot		Changeable Data	
OBX.25.3	Given Name			Phil		Changeable Data	
OBX.25.4	Second and Further Given Names or Initials Thereof			J.		Changeable Data	
OBX.25.9	Assigning Authority						
OBX.25.9.2	Universal ID			2.16.840.1.113883.4.6		Changeable Data	
OBX.25.9.3	Universal ID Type			ISO		IG Fixed Data	
OBX.25.10	Name Type Code			L		Changeable Data	
OBX.25.13	Identifier Type Code			NPI		Changeable Data	
OBX.29	Observation Type			RSLT		Test Case Fixed Data	

NTE

Location	Data Element	Data	Categorization
NTE.1	Set ID - NTE	1	IG Fixed Data
		Salmonella gastrointestinal infections usually resolve in 5-7 days and most do not require treatment other than oral fluids.	

Location	Data Element	Data	Categorization
NTE.3	Comment	Persons with severe diarrhea may require rehydration with intravenous fluids. Antibiotic therapy can prolong the duration of excretion of non-typhoidal Salmonella and is recommended only for patients with severe illness (e.g., those with severe diarrhea, high fever, bloodstream infection, or who need hospitalization) or those at risk of severe disease or complications, including young infants, older adults (over 65 years old) and immunocompromised persons. Antibiotic resistance is increasing among some Salmonella bacteria; therefore, susceptibility testing can help guide appropriate therapy. Choices for antibiotic therapy for severe infections include fluoroquinolones, third-generation cephalosporins, and ampicillin (for susceptible infections).	Changeable Data

#### OBX

Location	Data Element	Data	Categorization
OBX.1	Set ID - OBX	3	IG Fixed Data
OBX.2	Value Type	CWE	Test Case Fixed Data
OBX.3	Observation Identifier		
OBX.3.1	Identifier	625-4	Test Case Fixed Data
OBX.3.2	Text	Bacteria identified in Stool by Culture	Test Case Fixed Data
OBX.3.3	Name of the Coding System	LN	Test Case Fixed Data
OBX.3.7	Coding System Version	2.52	Changeable Data
OBX.3.9	Original Text	Stool Culture	Test Case Fixed Data
OBX.4	Observation Sub-ID		
OBX.4.2	Group	3	Changeable Data
OBX.4.3	Sequence	1	Changeable Data
OBX.4.4	Identifier	IsIt-3	Changeable Data
OBX.5	Observation Value		
OBX.5.1	Identifier	85729005	Changeable Data
OBX.5.2	Text	Shigella flexneri	Changeable Data
OBX.5.3	Name of the Coding System	SCT	Changeable Data
OBX.5.7	Coding System Version	201509USED	Changeable Data
OBX.5.9	Original Text	Shigella flexneri isolated	Changeable Data
OBX.8	Abnormal Flags	A	Test Case Fixed Data
OBX.11	Observation Result Status	F	Test Case Fixed Data
OBX.14	Date/Time of the Observation		
OBX.14.1	Time	201509231400	Changeable Data
OBX.19	Date/Time of the Analysis		
OBX.19.1	Time	201509251930	Changeable Data
OBX.23	Performing Organization Name		
OBX.23.1	Organization Name	Century Hospital	Changeable Data
OBX.23.6	Assigning Authority		
OBX.23.6.2	Universal ID	2.16.840.1.113883.4.7	Changeable Data
OBX.23.6.3	Universal ID Type	ISO	IG Fixed Data
OBX.23.7	Identifier Type Code	XX	Changeable Data
OBX.23.10	Organization Identifier	24D9871327	Changeable Data
OBX.24	Performing Organization Address		
OBX.24.1	Street Address		
OBX.24.1.1	Street or Mailing Address	2070 Test Park	Changeable Data
OBX.24.3	City	Los Angeles	Changeable Data
OBX.24.4	State or Province	CA	Changeable Data
OBX.24.5	Zip or Postal Code	90067	Changeable Data
OBX.24.7	Address Type	B	Changeable Data
OBX.25	Performing Organization Medical		



Location	Director	Data Element	Data	Categorization
OBX.25.1	ID Number		5432178916	Changeable Data
OBX.25.2	Family Name			
OBX.25.2.1	Surname		Knowsalot	Changeable Data
OBX.25.3	Given Name		Phil	Changeable Data
OBX.25.4	Second and Further Given Names or Initials Thereof		J.	Changeable Data
OBX.25.9	Assigning Authority			
OBX.25.9.2	Universal ID		2.16.840.1.113883.4.6	Changeable Data
OBX.25.9.3	Universal ID Type		ISO	IG Fixed Data
OBX.25.10	Name Type Code		L	Changeable Data
OBX.25.13	Identifier Type Code		NPI	Changeable Data
<b>OBX.29</b>	<b>Observation Type</b>		RSLT	Test Case Fixed Data

#### NTE

Location	Data Element	Data	Categorization
<b>NTE.1</b>	<b>Set ID - NTE</b>	1	IG Fixed Data
<b>NTE.3</b>	<b>Comment</b>	Antibiotic treatment is recommended for patients with severe disease, bloody diarrhea, or compromised immune systems. Resistance to traditional first-line drugs like ampicillin and trimethoprim-sulfamethoxazole is common, and resistance to some other antibiotics is increasing. With this in mind, antibiotic susceptibility testing can help guide appropriate therapy. When an ampicillin- or trimethoprim-sulfamethoxazole-resistant strain is isolated, choices for therapy include fluoroquinolones, ceftriaxone, and azithromycin.	Changeable Data

#### SPM

Location	Data Element	Data	Categorization
<b>SPM.1</b>	<b>Set ID - SPM</b>	1	IG Fixed Data
<b>SPM.2</b>	<b>Specimen ID</b>		
SPM.2.1			
SPM.2.1.1		S-2015-66	Configurable Data
SPM.2.1.3		2.16.840.1.113883.3.72.5.24	Configurable Data
SPM.2.1.4		ISO	Changeable Data
SPM.2.2			
SPM.2.2.1		S-9911-33	Changeable Data
SPM.2.2.3		2.16.840.1.113883.3.72.5.25	Changeable Data
SPM.2.2.4		ISO	IG Fixed Data
<b>SPM.4</b>	<b>Specimen Type</b>		
SPM.4.1	Identifier	119339001	Changeable Data
SPM.4.2	Text	Stool specimen	Changeable Data
SPM.4.3	Name of the Coding System	SCT	Changeable Data
SPM.4.7	Coding System Version	201509USEd	Changeable Data
SPM.4.9	Original Text	Stool	Changeable Data
<b>SPM.17</b>	<b>Specimen Collection Date/Time</b>		
SPM.17.1	Range Start Date/Time		
SPM.17.1.1	Time	201509231400	Changeable Data

#### ORC

Location	Data Element	Data	Categorization
<b>ORC.1</b>	<b>Order Control</b>	RE	Test Case Fixed Data
<b>ORC.3</b>	<b>Filler Order Number</b>		
ORC.3.1	Entity Identifier	R-783274-6	Changeable Data
ORC.3.3	Universal ID	2.16.840.1.113883.3.72.5.25	Changeable Data



Location	Data Element	Data	Categorization
<b>ORC.4</b>	<b>Placer Group Number</b>		
ORC.4.1	Entity Identifier	GORD874211	Changeable Data
ORC.4.3	Universal ID	2.16.840.1.113883.3.72.5.24	Changeable Data
ORC.4.4	Universal ID Type	ISO	IG Fixed Data
<b>ORC.12</b>	<b>Ordering Provider</b>		
ORC.12.1	ID Number	5742200012	Changeable Data
ORC.12.2	Family Name		
ORC.12.2.1	Surname	Radon	Changeable Data
ORC.12.3	Given Name	Nicholas	Changeable Data
ORC.12.9	Assigning Authority		
ORC.12.9.2	Universal ID	2.16.840.1.113883.4.6	Changeable Data
ORC.12.9.3	Universal ID Type	ISO	Changeable Data
ORC.12.10	Name Type Code	L	Changeable Data
ORC.12.13	Identifier Type Code	NPI	Changeable Data

#### OBR

Location	Data Element	Data	Categorization
<b>OBR.1</b>	<b>Set ID - OBR</b>	2	IG Fixed Data
<b>OBR.3</b>	<b>Filler Order Number</b>		
OBR.3.1	Entity Identifier	R-783274-6	Changeable Data
OBR.3.3	Universal ID	2.16.840.1.113883.3.72.5.25	Changeable Data
OBR.3.4	Universal ID Type	ISO	Changeable Data
<b>OBR.4</b>	<b>Universal Service Identifier</b>		
OBR.4.1	Identifier	50545-3	Test Case Fixed Data
OBR.4.2	Text	Bacterial susceptibility panel in Isolate by Minimum inhibitory concentration (MIC)	Test Case Fixed Data
OBR.4.3	Name of Coding System	LN	Test Case Fixed Data
OBR.4.7	Coding System Version	2.52	Changeable Data
OBR.4.9	Original Text	Bacteria susceptibility	Changeable Data
<b>OBR.7</b>	<b>Observation Date/Time</b>		
OBR.7.1	Time	201509231400	Changeable Data
<b>OBR.11</b>	<b>Specimen Action Code</b>	G	Test Case Fixed Data
<b>OBR.16</b>	<b>Ordering Provider</b>		
OBR.16.1	ID Number	5742200012	Changeable Data
OBR.16.2	Family Name		
OBR.16.2.1	Surname	Radon	Changeable Data
OBR.16.3	Given Name	Nicholas	Changeable Data
OBR.16.9	Assigning Authority		
OBR.16.9.2	Universal ID	2.16.840.1.113883.4.6	Changeable Data
OBR.16.9.3	Universal ID Type	ISO	Changeable Data
OBR.16.10	Name Type Code	L	Changeable Data
OBR.16.13	Identifier Type Code	NPI	Changeable Data
<b>OBR.22</b>	<b>Results Rpt/Status Chng - Date/Time</b>		
OBR.22.1	Time	20150927112054	Changeable Data
<b>OBR.25</b>	<b>Result Status</b>	F	Test Case Fixed Data
<b>OBR.26</b>	<b>Parent Result</b>		
OBR.26.1	Parent Observation Identifier		
OBR.26.1.1	Identifier	625-4	Test Case Fixed Data
OBR.26.1.2	Text	Bacteria identified in Stool by Culture	Test Case Fixed Data
OBR.26.1.3	Name of the Coding System	LN	Test Case Fixed Data
OBR.26.1.7	Coding System Version	2.52	Changeable Data
OBR.26.1.9	Original Text	Stool Culture	Changeable Data
OBR.26.2	Parent Observation Sub-Identifier		
OBR.26.2.2	Group	2	Test Case Fixed Data
OBR.26.2.3	Sequence	1	Test Case Fixed Data
OBR.26.2.4	Identifier	Islt-2	Test Case Fixed Data

<b>OBR.28</b>	<b>Location</b>	<b>Result Copies To Data Element</b>	<b>Data</b>	<b>Categorization</b>
OBR.28.1	ID Number		10092000194	Changeable Data
OBR.28.2	Family Name			
OBR.28.2.1	Surname		Hamlin	Changeable Data
OBR.28.3	Given Name		Pafford	Changeable Data
OBR.28.9	Assigning Authority			
OBR.28.9.2	Universal ID		2.16.840.1.113883.4.6	Changeable Data
OBR.28.9.3	Universal ID Type		ISO	Changeable Data
OBR.28.10	Name Type Code		L	Changeable Data
OBR.28.13	Identifier Type Code		NPI	Changeable Data
<b>OBR.29</b>	<b>Parent</b>			
OBR.29.1	Placer Assigned Identifier			
OBR.29.1.1	Entity Identifier		ORD723222-4	Changeable Data
OBR.29.1.3	Universal ID		2.16.840.1.113883.3.72.5.24	Changeable Data
OBR.29.1.4	Universal ID Type		ISO	Changeable Data
OBR.29.2	Filler Assigned Identifier			
OBR.29.2.1	Entity Identifier		R-783274-4	Changeable Data
OBR.29.2.3	Universal ID		2.16.840.1.113883.3.72.5.25	Changeable Data
OBR.29.2.4	Universal ID Type		ISO	Changeable Data
<b>OBR.47</b>	<b>Filler Supplemental Service Information</b>			
OBR.47.1	Identifier		MIC	Test Case Fixed Data
OBR.47.2	Text		Observation of type microbiology	Test Case Fixed Data
OBR.47.3	Name of the Coding System		HL70411	IG Fixed Data
OBR.47.7	Coding System Version		2.5.1	IG Fixed Data
<b>OBR.49</b>	<b>Results Handling</b>			
OBR.49.1	Identifier		CC	Test Case Fixed Data
OBR.49.2	Text		Copies Requested	Test Case Fixed Data
OBR.49.3	Name of the Coding System		HL70507	IG Fixed Data

#### OBX

<b>Location</b>	<b>Data Element</b>	<b>Data</b>	<b>Categorization</b>
<b>OBX.1</b>	<b>Set ID - OBX</b>	1	IG Fixed Data
<b>OBX.2</b>	<b>Value Type</b>	SN	Test Case Fixed Data
<b>OBX.3</b>	<b>Observation Identifier</b>		
OBX.3.1	Identifier	28-1	Test Case Fixed Data
OBX.3.2	Text	Ampicillin [Susceptibility] by Minimum inhibitory concentration (MIC)	Test Case Fixed Data
OBX.3.3	Name of the Coding System	LN	Test Case Fixed Data
OBX.3.7	Coding System Version	2.52	Changeable Data
<b>OBX.4</b>	<b>Observation Sub-ID</b>		
OBX.4.2	Group	1	Test Case Fixed Data
OBX.4.3	Sequence	1	Test Case Fixed Data
OBX.4.4	Identifier	IsIt-2	Test Case Fixed Data
<b>OBX.5</b>	<b>Observation Value</b>		
OBX.5.1	Comparator	<	Test Case Fixed Data
OBX.5.2	Num1	0.06	Test Case Fixed Data
<b>OBX.6</b>	<b>Units</b>		
OBX.6.1	Identifier	ug/mL	Changeable Data
OBX.6.3	Name of the Coding System	UCUM	Changeable Data
OBX.6.7	Coding System Version	1.9	Changeable Data
<b>OBX.8</b>	<b>Abnormal Flags</b>	S	Test Case Fixed Data
<b>OBX.11</b>	<b>Observation Result Status</b>	F	Test Case Fixed Data
<b>OBX.14</b>	<b>Date/Time of the Observation</b>		
OBX.14.1	Time	201509231400	Changeable Data
<b>OBX.19</b>	<b>Date/Time of the Analysis</b>		
OBX.19.1	Time	201509261100	Changeable Data

OBX.23	Location	Performing Organization Name	Data	Categorization
	OBX.23.1	Organization Name	Century Hospital	Changeable Data
	OBX.23.6	Assigning Authority		
	OBX.23.6.2	Universal ID	2.16.840.1.113883.4.7	Changeable Data
	OBX.23.6.3	Universal ID Type	ISO	IG Fixed Data
	OBX.23.7	Identifier Type Code	XX	Changeable Data
	OBX.23.10	Organization Identifier	24D9871327	Changeable Data
<b>OBX.24</b>	<b>Performing Organization Address</b>			
	OBX.24.1	Street Address		
	OBX.24.1.1	Street or Mailing Address	2070 Test Park	Changeable Data
	OBX.24.3	City	Los Angeles	Changeable Data
	OBX.24.4	State or Province	CA	Changeable Data
	OBX.24.5	Zip or Postal Code	90067	Changeable Data
	OBX.24.7	Address Type	B	Changeable Data
<b>OBX.25</b>	<b>Performing Organization Medical Director</b>			
	OBX.25.1	ID Number	5432178916	Changeable Data
	OBX.25.2	Family Name		
	OBX.25.2.1	Surname	Knowsalot	Changeable Data
	OBX.25.3	Given Name	Phil	Changeable Data
	OBX.25.4	Second and Further Given Names or Initials Thereof	J.	Changeable Data
	OBX.25.9	Assigning Authority		
	OBX.25.9.2	Universal ID	2.16.840.1.113883.4.6	Changeable Data
	OBX.25.9.3	Universal ID Type	ISO	IG Fixed Data
	OBX.25.10	Name Type Code	L	Changeable Data
	OBX.25.13	Identifier Type Code	NPI	Changeable Data
<b>OBX.29</b>	<b>Observation Type</b>		RSLT	Test Case Fixed Data
<b>OBX.30</b>	<b>Observation SubType</b>		SUR	Test Case Fixed Data

OBX

	Location	Data Element	Data	Categorization
<b>OBX.1</b>	<b>Set ID - OBX</b>		2	IG Fixed Data
<b>OBX.2</b>	<b>Value Type</b>		SN	Test Case Fixed Data
<b>OBX.3</b>	<b>Observation Identifier</b>			
	OBX.3.1	Identifier	267-5	Test Case Fixed Data
	OBX.3.2	Text	Gentamicin [Susceptibility] by Minimum inhibitory concentration (MIC)	Test Case Fixed Data
	OBX.3.3	Name of the Coding System	LN	Test Case Fixed Data
	OBX.3.7	Coding System Version	2.52	Changeable Data
<b>OBX.4</b>	<b>Observation Sub-ID</b>			
	OBX.4.2	Group	1	Test Case Fixed Data
	OBX.4.3	Sequence	1	Test Case Fixed Data
	OBX.4.4	Identifier	IsIt-2	Test Case Fixed Data
<b>OBX.5</b>	<b>Observation Value</b>			
	OBX.5.2	Num1	0.05	Test Case Fixed Data
<b>OBX.6</b>	<b>Units</b>			
	OBX.6.1	Identifier	ug/mL	Changeable Data
	OBX.6.3	Name of the Coding System	UCUM	Changeable Data
	OBX.6.7	Coding System Version	1.9	Changeable Data
<b>OBX.8</b>	<b>Abnormal Flags</b>		S	Test Case Fixed Data
<b>OBX.11</b>	<b>Observation Result Status</b>		F	Test Case Fixed Data
<b>OBX.14</b>	<b>Date/Time of the Observation</b>			
	OBX.14.1	Time	201509231400	Changeable Data
<b>OBX.19</b>	<b>Date/Time of the Analysis</b>			
	OBX.19.1	Time	201509261100	Changeable Data
<b>OBX.23</b>	<b>Performing Organization Name</b>			

OBX.23.1	Organization Name	Century Hospital	Changeable Data
Location	Data Element	Data	Categorization
OBX.23.6	Assigning Authority		
OBX.23.6.2	Universal ID	2.16.840.1.113883.4.7	Changeable Data
OBX.23.6.3	Universal ID Type	ISO	IG Fixed Data
OBX.23.7	Identifier Type Code	XX	Changeable Data
OBX.23.10	Organization Identifier	24D9871327	Changeable Data
<b>OBX.24</b>	<b>Performing Organization Address</b>		
OBX.24.1	Street Address		
OBX.24.1.1	Street or Mailing Address	2070 Test Park	Changeable Data
OBX.24.3	City	Los Angeles	Changeable Data
OBX.24.4	State or Province	CA	Changeable Data
OBX.24.5	Zip or Postal Code	90067	Changeable Data
OBX.24.7	Address Type	B	Changeable Data
<b>OBX.25</b>	<b>Performing Organization Medical Director</b>		
OBX.25.1	ID Number	5432178916	Changeable Data
OBX.25.2	Family Name		
OBX.25.2.1	Surname	Knowsalot	Changeable Data
OBX.25.3	Given Name	Phil	Changeable Data
OBX.25.4	Second and Further Given Names or Initials Thereof	J.	Changeable Data
OBX.25.9	Assigning Authority		
OBX.25.9.2	Universal ID	2.16.840.1.113883.4.6	Changeable Data
OBX.25.9.3	Universal ID Type	ISO	IG Fixed Data
OBX.25.10	Name Type Code	L	Changeable Data
OBX.25.13	Identifier Type Code	NPI	Changeable Data
<b>OBX.29</b>	<b>Observation Type</b>	RSLT	Test Case Fixed Data
<b>OBX.30</b>	<b>Observation SubType</b>	SUR	Test Case Fixed Data

#### OBX

Location	Data Element	Data	Categorization
<b>OBX.1</b>	<b>Set ID - OBX</b>	3	IG Fixed Data
<b>OBX.2</b>	<b>Value Type</b>	SN	Test Case Fixed Data
<b>OBX.3</b>	<b>Observation Identifier</b>		
OBX.3.1	Identifier	185-9	Test Case Fixed Data
OBX.3.2	Text	Ciprofloxacin [Susceptibility] by Minimum inhibitory concentration (MIC)	Test Case Fixed Data
OBX.3.3	Name of the Coding System	LN	Test Case Fixed Data
OBX.3.7	Coding System Version	2.52	Changeable Data
<b>OBX.4</b>	<b>Observation Sub-ID</b>		
OBX.4.2	Group	1	Test Case Fixed Data
OBX.4.3	Sequence	1	Test Case Fixed Data
OBX.4.4	Identifier	Islt-2	Test Case Fixed Data
<b>OBX.5</b>	<b>Observation Value</b>		
OBX.5.2	Num1	0.05	Test Case Fixed Data
<b>OBX.6</b>	<b>Units</b>		
OBX.6.1	Identifier	ug/mL	Changeable Data
OBX.6.3	Name of the Coding System	UCUM	Changeable Data
OBX.6.7	Coding System Version	1.9	Changeable Data
<b>OBX.8</b>	<b>Abnormal Flags</b>	S	Test Case Fixed Data
<b>OBX.11</b>	<b>Observation Result Status</b>	F	Test Case Fixed Data
<b>OBX.14</b>	<b>Date/Time of the Observation</b>		
OBX.14.1	Time	201509231400	Changeable Data
<b>OBX.19</b>	<b>Date/Time of the Analysis</b>		
OBX.19.1	Time	201509261100	Changeable Data
<b>OBX.23</b>	<b>Performing Organization Name</b>		
OBX.23.1	Organization Name	Century Hospital	Changeable Data

Location	Data Element	Data	Categorization
OBX.23.6.2	Universal ID	2.16.840.1.113883.4.7	Changeable Data
OBX.23.6.3	Universal ID Type	ISO	IG Fixed Data
OBX.23.7	Identifier Type Code	XX	Changeable Data
OBX.23.10	Organization Identifier	24D9871327	Changeable Data
<b>OBX.24</b>	<b>Performing Organization Address</b>		
OBX.24.1	Street Address		
OBX.24.1.1	Street or Mailing Address	2070 Test Park	Changeable Data
OBX.24.3	City	Los Angeles	Changeable Data
OBX.24.4	State or Province	CA	Changeable Data
OBX.24.5	Zip or Postal Code	90067	Changeable Data
OBX.24.7	Address Type	B	Changeable Data
<b>OBX.25</b>	<b>Performing Organization Medical Director</b>		
OBX.25.1	ID Number	5432178916	Changeable Data
OBX.25.2	Family Name		
OBX.25.2.1	Surname	Knowsalot	Changeable Data
OBX.25.3	Given Name	Phil	Changeable Data
OBX.25.4	Second and Further Given Names or Initials Thereof	J.	Changeable Data
OBX.25.9	Assigning Authority		
OBX.25.9.2	Universal ID	2.16.840.1.113883.4.6	Changeable Data
OBX.25.9.3	Universal ID Type	ISO	IG Fixed Data
OBX.25.10	Name Type Code	L	Changeable Data
OBX.25.13	Identifier Type Code	NPI	Changeable Data
<b>OBX.29</b>	<b>Observation Type</b>	RSLT	Test Case Fixed Data
<b>OBX.30</b>	<b>Observation SubType</b>	SUR	Test Case Fixed Data

#### ORC

Location	Data Element	Data	Categorization
<b>ORC.1</b>	<b>Order Control</b>	RE	Test Case Fixed Data
<b>ORC.3</b>	<b>Filler Order Number</b>		
ORC.3.1	Entity Identifier	R-783274-7	Changeable Data
ORC.3.3	Universal ID	2.16.840.1.113883.3.72.5.25	Changeable Data
ORC.3.4	Universal ID Type	ISO	IG Fixed Data
<b>ORC.4</b>	<b>Placer Group Number</b>		
ORC.4.1	Entity Identifier	GORD874211	Changeable Data
ORC.4.3	Universal ID	2.16.840.1.113883.3.72.5.24	Changeable Data
ORC.4.4	Universal ID Type	ISO	IG Fixed Data
<b>ORC.12</b>	<b>Ordering Provider</b>		
ORC.12.1	ID Number	5742200012	Changeable Data
ORC.12.2	Family Name		
ORC.12.2.1	Surname	Radon	Changeable Data
ORC.12.3	Given Name	Nicholas	Changeable Data
ORC.12.9	Assigning Authority		
ORC.12.9.2	Universal ID	2.16.840.1.113883.4.6	Changeable Data
ORC.12.9.3	Universal ID Type	ISO	Changeable Data
ORC.12.10	Name Type Code	L	Changeable Data
ORC.12.13	Identifier Type Code	NPI	Changeable Data

#### OBR

Location	Data Element	Data	Categorization
<b>OBR.1</b>	<b>Set ID - OBR</b>	3	IG Fixed Data
<b>OBR.3</b>	<b>Filler Order Number</b>		
OBR.3.1	Entity Identifier	R-783274-7	Changeable Data
OBR.3.3	Universal ID	2.16.840.1.113883.3.72.5.25	Changeable Data
OBR.3.4	Universal ID Type	ISO	Changeable Data

<b>OBR.4</b>	<b>Location</b>	<b>Universal Service Data Identifier</b>	<b>Data</b>	<b>Categorization</b>
OBR.4.1	Identifier		50545-3	Test Case Fixed Data
OBR.4.2	Text		Bacterial susceptibility panel in Isolate by Minimum inhibitory concentration (MIC)	Test Case Fixed Data
OBR.4.3	Name of Coding System		LN	Test Case Fixed Data
OBR.4.7	Coding System Version		2.52	Changeable Data
OBR.4.9	Original Text		Bacteria susceptibility	Changeable Data
<b>OBR.7</b>	<b>Observation Date/Time</b>			
OBR.7.1	Time		201509231400	Changeable Data
<b>OBR.11</b>	<b>Specimen Action Code</b>		G	Test Case Fixed Data
<b>OBR.16</b>	<b>Ordering Provider</b>			
OBR.16.1	ID Number		5742200012	Changeable Data
OBR.16.2	Family Name			
OBR.16.2.1	Surname		Radon	Changeable Data
OBR.16.3	Given Name		Nicholas	Changeable Data
OBR.16.9	Assigning Authority			
OBR.16.9.2	Universal ID		2.16.840.1.113883.4.6	Changeable Data
OBR.16.9.3	Universal ID Type		ISO	Changeable Data
OBR.16.10	Name Type Code		L	Changeable Data
OBR.16.13	Identifier Type Code		NPI	Changeable Data
<b>OBR.22</b>	<b>Results Rpt/Status Chng - Date/Time</b>			
OBR.22.1	Time		20150927164251	Changeable Data
<b>OBR.25</b>	<b>Result Status</b>		C	Test Case Fixed Data
<b>OBR.26</b>	<b>Parent Result</b>			
OBR.26.1	Parent Observation Identifier			
OBR.26.1.1	Identifier		625-4	Test Case Fixed Data
OBR.26.1.2	Text		Bacteria identified in Stool by Culture	Test Case Fixed Data
OBR.26.1.3	Name of the Coding System		LN	Test Case Fixed Data
OBR.26.1.7	Coding System Version		2.52	Changeable Data
OBR.26.1.9	Original Text		Stool Culture	Changeable Data
OBR.26.2	Parent Observation Sub-Identifier			
OBR.26.2.2	Group		3	Test Case Fixed Data
OBR.26.2.3	Sequence		1	Test Case Fixed Data
OBR.26.2.4	Identifier		Islt-3	Test Case Fixed Data
<b>OBR.28</b>	<b>Result Copies To</b>			
OBR.28.1	ID Number		10092000194	Changeable Data
OBR.28.2	Family Name			
OBR.28.2.1	Surname		Hamlin	Changeable Data
OBR.28.3	Given Name		Pafford	Changeable Data
OBR.28.9	Assigning Authority			
OBR.28.9.2	Universal ID		2.16.840.1.113883.4.6	Changeable Data
OBR.28.9.3	Universal ID Type		ISO	Changeable Data
OBR.28.10	Name Type Code		L	Changeable Data
OBR.28.13	Identifier Type Code		NPI	Changeable Data
<b>OBR.29</b>	<b>Parent</b>			
OBR.29.1	Placer Assigned Identifier			
OBR.29.1.1	Entity Identifier		ORD723222-4	Changeable Data
OBR.29.1.3	Universal ID		2.16.840.1.113883.3.72.5.24	Changeable Data
OBR.29.1.4	Universal ID Type		ISO	Changeable Data
OBR.29.2	Filler Assigned Identifier			
OBR.29.2.1	Entity Identifier		R-783274-4	Changeable Data
OBR.29.2.3	Universal ID		2.16.840.1.113883.3.72.5.25	Changeable Data
OBR.29.2.4	Universal ID Type		ISO	Changeable Data
<b>OBR.47</b>	<b>Filler Supplemental Service Information</b>			
OBR.47.1	Identifier		MIC	Test Case Fixed Data
OBR.47.2	Text		Observation of type microbiology	Test Case Fixed Data



OBR.47.2	Name of the Coding System	HL70411	IG Fixed Data
OBR.47.7	Coding System Version	2.5.1	IG Fixed Data
<b>OBR.49</b>	<b>Results Handling</b>		
OBR.49.1	Identifier	CC	Test Case Fixed Data
OBR.49.2	Text	Copies Requested	Test Case Fixed Data
OBR.49.3	Name of the Coding System	HL70507	IG Fixed Data

#### OBX

Location	Data Element	Data	Categorization
<b>OBX.1</b>	<b>Set ID - OBX</b>	1	IG Fixed Data
<b>OBX.2</b>	<b>Value Type</b>	SN	Test Case Fixed Data
<b>OBX.3</b>	<b>Observation Identifier</b>		
OBX.3.1	Identifier	28-1	Test Case Fixed Data
OBX.3.2	Text	Ampicillin [Susceptibility] by Minimum inhibitory concentration (MIC)	Test Case Fixed Data
OBX.3.3	Name of the Coding System	LN	Test Case Fixed Data
OBX.3.7	Coding System Version	2.52	Changeable Data
<b>OBX.4</b>	<b>Observation Sub-ID</b>		
OBX.4.2	Group	1	Test Case Fixed Data
OBX.4.3	Sequence	1	Test Case Fixed Data
OBX.4.4	Identifier	IsIt-3	Test Case Fixed Data
<b>OBX.5</b>	<b>Observation Value</b>		
OBX.5.1	Comparator	<	Test Case Fixed Data
OBX.5.2	Num1	16	Test Case Fixed Data
<b>OBX.6</b>	<b>Units</b>		
OBX.6.1	Identifier	ug/mL	Changeable Data
OBX.6.3	Name of the Coding System	UCUM	Changeable Data
OBX.6.7	Coding System Version	1.9	Changeable Data
<b>OBX.8</b>	<b>Abnormal Flags</b>	I	Test Case Fixed Data
<b>OBX.11</b>	<b>Observation Result Status</b>	F	Test Case Fixed Data
<b>OBX.14</b>	<b>Date/Time of the Observation</b>		
OBX.14.1	Time	201509231400	Changeable Data
<b>OBX.19</b>	<b>Date/Time of the Analysis</b>		
OBX.19.1	Time	201509261100	Changeable Data
<b>OBX.23</b>	<b>Performing Organization Name</b>		
OBX.23.1	Organization Name	Century Hospital	Changeable Data
OBX.23.6	Assigning Authority		
OBX.23.6.2	Universal ID	2.16.840.1.113883.4.7	Changeable Data
OBX.23.6.3	Universal ID Type	ISO	IG Fixed Data
OBX.23.7	Identifier Type Code	XX	Changeable Data
OBX.23.10	Organization Identifier	24D9871327	Changeable Data
<b>OBX.24</b>	<b>Performing Organization Address</b>		
OBX.24.1	Street Address		
OBX.24.1.1	Street or Mailing Address	2070 Test Park	Changeable Data
OBX.24.3	City	Los Angeles	Changeable Data
OBX.24.4	State or Province	CA	Changeable Data
OBX.24.5	Zip or Postal Code	90067	Changeable Data
OBX.24.7	Address Type	B	Changeable Data
<b>OBX.25</b>	<b>Performing Organization Medical Director</b>		
OBX.25.1	ID Number	5432178916	Changeable Data
OBX.25.2	Family Name		
OBX.25.2.1	Surname	Knowsalot	Changeable Data
OBX.25.3	Given Name	Phil	Changeable Data
OBX.25.4	Second and Further Given Names or Initials Thereof	J.	Changeable Data
OBX.25.9	Assigning Authority		
OBX.25.9.2	Universal ID	2.16.840.1.113883.4.6	Changeable Data



OBX.25.9.3 Location	Universal ID Type Data Element	ISO Data	IG Fixed Data Categorization
OBX.25.10	Name Type Code	L	Changeable Data
OBX.25.13	Identifier Type Code	NPI	Changeable Data
<b>OBX.29</b>	<b>Observation Type</b>	RSLT	Test Case Fixed Data
<b>OBX.30</b>	<b>Observation SubType</b>	SUR	Test Case Fixed Data

#### OBX

Location	Data Element	Data	Categorization
<b>OBX.1</b>	<b>Set ID - OBX</b>	2	IG Fixed Data
<b>OBX.2</b>	<b>Value Type</b>	SN	Test Case Fixed Data
<b>OBX.3</b>	<b>Observation Identifier</b>		
OBX.3.1	Identifier	516-5	Test Case Fixed Data
OBX.3.2	Text	Trimethoprim+Sulfamethoxazole [Susceptibility] by Minimum inhibitory concentration (MIC)	Test Case Fixed Data
OBX.3.3	Name of the Coding System	LN	Test Case Fixed Data
OBX.3.7	Coding System Version	2.52	Changeable Data
<b>OBX.4</b>	<b>Observation Sub-ID</b>		
OBX.4.2	Group	1	Test Case Fixed Data
OBX.4.3	Sequence	1	Test Case Fixed Data
OBX.4.4	Identifier	Islt-3	Test Case Fixed Data
<b>OBX.5</b>	<b>Observation Value</b>		
OBX.5.2	Num1	2	Test Case Fixed Data
OBX.5.3	Separator/Suffix	/	
OBX.5.4	Num2	38	
<b>OBX.6</b>	<b>Units</b>		
OBX.6.1	Identifier	ug/mL	Changeable Data
OBX.6.3	Name of the Coding System	UCUM	Changeable Data
OBX.6.7	Coding System Version	1.9	Changeable Data
<b>OBX.8</b>	<b>Abnormal Flags</b>	S	Test Case Fixed Data
<b>OBX.11</b>	<b>Observation Result Status</b>	B	Test Case Fixed Data
<b>OBX.14</b>	<b>Date/Time of the Observation</b>		
OBX.14.1	Time	201509231400	Changeable Data
<b>OBX.19</b>	<b>Date/Time of the Analysis</b>		
OBX.19.1	Time	201509271120	Changeable Data
<b>OBX.23</b>	<b>Performing Organization Name</b>		
OBX.23.1	Organization Name	Century Hospital	Changeable Data
OBX.23.6	Assigning Authority		
OBX.23.6.2	Universal ID	2.16.840.1.113883.4.7	Changeable Data
OBX.23.6.3	Universal ID Type	ISO	IG Fixed Data
OBX.23.7	Identifier Type Code	XX	Changeable Data
OBX.23.10	Organization Identifier	24D9871327	Changeable Data
<b>OBX.24</b>	<b>Performing Organization Address</b>		
OBX.24.1	Street Address		
OBX.24.1.1	Street or Mailing Address	2070 Test Park	Changeable Data
OBX.24.3	City	Los Angeles	Changeable Data
OBX.24.4	State or Province	CA	Changeable Data
OBX.24.5	Zip or Postal Code	90067	Changeable Data
OBX.24.7	Address Type	B	Changeable Data
<b>OBX.25</b>	<b>Performing Organization Medical Director</b>		
OBX.25.1	ID Number	5432178916	Changeable Data
OBX.25.2	Family Name		
OBX.25.2.1	Surname	Knowsalot	Changeable Data
OBX.25.3	Given Name	Phil	Changeable Data
OBX.25.4	Second and Further Given Names or Initials Thereof	J.	Changeable Data
OBX.25.9	Assigning Authority		

Location	Data Element	Data	Categorization
OBX.25.9.3	Universal ID Type	ISO	IG Fixed Data
OBX.25.10	Name Type Code	L	Changeable Data
OBX.25.13	Identifier Type Code	NPI	Changeable Data
OBX.29	Observation Type	RSLT	Test Case Fixed Data
OBX.30	Observation SubType	SUR	Test Case Fixed Data

NTE

Location	Data Element	Data	Categorization
NTE.1	Set ID - NTE	1	IG Fixed Data
NTE.3	Comment	Due to the indeterminate amoxicillin test result, additional antibiotics were tested and appended to the previous report.	Changeable Data

OBX

Location	Data Element	Data	Categorization
OBX.1	Set ID - OBX	3	IG Fixed Data
OBX.2	Value Type	SN	Test Case Fixed Data
OBX.3	Observation Identifier		
OBX.3.1	Identifier	185-9	Test Case Fixed Data
OBX.3.2	Text	Ciprofloxacin [Susceptibility] by Minimum inhibitory concentration (MIC)	Test Case Fixed Data
OBX.3.3	Name of the Coding System	LN	Test Case Fixed Data
OBX.3.7	Coding System Version	2.52	Changeable Data
OBX.4	Observation Sub-ID		
OBX.4.2	Group	1	Test Case Fixed Data
OBX.4.3	Sequence	1	Test Case Fixed Data
OBX.4.4	Identifier	IsIt-3	Test Case Fixed Data
OBX.5	Observation Value		
OBX.5.2	Num1	0.05	Test Case Fixed Data
OBX.6	Units		
OBX.6.1	Identifier	ug/mL	Changeable Data
OBX.6.3	Name of the Coding System	UCUM	Changeable Data
OBX.6.7	Coding System Version	1.9	Changeable Data
OBX.8	Abnormal Flags	S	Test Case Fixed Data
OBX.11	Observation Result Status	B	Test Case Fixed Data
OBX.14	Date/Time of the Observation		
OBX.14.1	Time	201509231400	Changeable Data
OBX.19	Date/Time of the Analysis		
OBX.19.1	Time	201509271120	Changeable Data
OBX.23	Performing Organization Name		
OBX.23.1	Organization Name	Century Hospital	Changeable Data
OBX.23.6	Assigning Authority		
OBX.23.6.2	Universal ID	2.16.840.1.113883.4.7	Changeable Data
OBX.23.6.3	Universal ID Type	ISO	IG Fixed Data
OBX.23.7	Identifier Type Code	XX	Changeable Data
OBX.23.10	Organization Identifier	24D9871327	Changeable Data
OBX.24	Performing Organization Address		
OBX.24.1	Street Address		
OBX.24.1.1	Street or Mailing Address	2070 Test Park	Changeable Data
OBX.24.3	City	Los Angeles	Changeable Data
OBX.24.4	State or Province	CA	Changeable Data
OBX.24.5	Zip or Postal Code	90067	Changeable Data
OBX.24.7	Address Type	B	Changeable Data
OBX.25	Performing Organization Medical Director		
OBX.25.1	ID Number	5432178916	Changeable Data
OBX.25.2	Family Name		

OBX.25.2.1 Location	Surname Data Element	Knowsalot Data	Changeable Data Categorization
OBX.25.3	Given Name	Phil	Changeable Data
OBX.25.4	Second and Further Given Names or Initials Thereof	J.	Changeable Data
OBX.25.9	Assigning Authority		
OBX.25.9.2	Universal ID	2.16.840.1.113883.4.6	Changeable Data
OBX.25.9.3	Universal ID Type	ISO	IG Fixed Data
OBX.25.10	Name Type Code	L	Changeable Data
OBX.25.13	Identifier Type Code	NPI	Changeable Data
OBX.29	Observation Type	RSLT	Test Case Fixed Data
OBX.30	Observation SubType	SUR	Test Case Fixed Data

NTE

Location	Data Element	Data	Categorization
NTE.1	Set ID - NTE	1	IG Fixed Data
NTE.3	Comment	Due to the indeterminate amoxicillin test result, additional antibiotics were tested and appended to the previous report.	Changeable Data

**Patient Information**

Element	Data
Name	William A Jones
Date/Time of Birth	06/15/1961
Administrative Sex	Male
Race	White
Alt Race	

**Order Observation[\*]****Order Observation****Ordering Provider**

Element	Data
Name	Nicholas Radon
Identifier number	5742200012

**Observation Details**

Element	Data
Observation General Information	
Placer Order Number	ORD723222-4
Filler Order Number	R-783274-4
Placer Group Number	GORD874211
Parent Universal Service Identifier	
Identifier	
Text	
Alt Identifier	
Alt Text	
Original Text	
Observation Details	
Universal Service Identifier	Bacteria identified in Stool by Culture
Observation Date/Time	09/23/2015 2:00 PM
Observation end Date/Time	
Specimen Action Code	
Relevant Clinical Information	
Relevant Clinical Information Original Text	
Observation Result Information	
Result Status	F
Results Report/Status Change - Date/Time	09/26/2015 2:05 PM
Results Copy To	
Name	Pafford
Identifier	10092000194
Results Handling	
Standard	
Observation Notes	

**Timing/Quantity Information**

Element	Data
Priority	
Start Date/time	
End Date/time	

#### Results Performing Laboratory

Element	Data
Laboratory Name	Century Hospital
Organization identifier	24D9871327
Address	2070 Test Park Los Angeles CA 90067
Director Name	Phil J. Knowsalot
Director identifier	5432178916

#### Specimen Information

Element	Data
Specimen Type	Stool specimen
Alt Specimen Type	
Specimen Original Text	Stool
Start date/time	201509231400

#### Lab results

Element				Data				
Test performed				Stool Culture				
Test Report date				09/26/2015 14:05				
Result Observation Name	Result	UOM	Range	Abnormal Flag	Status	Date/Time of Observation	Date/Time of Analysis	Notes

Stool Culture	Shiga toxin producing E. coli O157:H7 isolated			A	F	09/23/2015 14:00	09/25/2015 19:30	<p>Susceptibility testing for E.coli is not performed, because antibiotics should not be used to treat this infection. There is no evidence that treatment with antibiotics is helpful, and taking antibiotics may increase the risk of hemolytic-uremic syndrome (HUS). Antidiarrheal agents like Imodium may also increase that risk. Non-specific supportive therapy, including hydration, is important.</p>
								<p>Salmonella gastrointestinal infections usually resolve in 5-7 days and most do not require treatment other than oral fluids. Persons with severe diarrhea may require rehydration with intravenous fluids. Antibiotic therapy can prolong the duration of excretion of non-typhoidal Salmonella and is recommended only for</p>

Stool Culture	Salmonella I, group O:4 isolated		A	F	09/23/2015 14:00	09/25/2015 19:30	<p>patients with severe illness (e.g., those with severe diarrhea, high fever, bloodstream infection, or who need hospitalization) or those at risk of severe disease or complications, including young infants, older adults (over 65 years old) and immunocompromised persons. Antibiotic resistance is increasing among some Salmonella bacteria; therefore, susceptibility testing can help guide appropriate therapy. Choices for antibiotic therapy for severe infections include fluoroquinolones, third-generation cephalosporins, and ampicillin (for susceptible infections).</p>
---------------	----------------------------------	--	---	---	------------------	------------------	---



Stool Culture	Shigella flexneri isolated			A	F	09/23/2015 14:00	09/25/2015 19:30	Antibiotic treatment is recommended for patients with severe disease, bloody diarrhea, or compromised immune systems. Resistance to traditional first-line drugs like ampicillin and trimethoprim-sulfamethoxazole is common, and resistance to some other antibiotics is increasing. With this in mind, antibiotic susceptibility testing can help guide appropriate therapy. When an ampicillin- or trimethoprim-sulfamethoxazole-resistant strain is isolated, choices for therapy include fluoroquinolones, ceftriaxone, and azithromycin.
---------------	----------------------------	--	--	---	---	------------------	------------------	--

#### Order Observation

#### Ordering Provider

Element	Data
Name	Nicholas Radon
Identifier number	5742200012

#### Observation Details

Element	Data
Observation General Information	
Placer Order Number	
Filler Order Number	R-783274-6
Placer Group Number	GORD874211

Parent Universal Service Identifier	
Identifier	
Text	
Alt Identifier	
Alt Text	
Original Text	
Observation Details	
Universal Service Identifier	Bacterial susceptibility panel in Isolate by Minimum inhibitory concentration (MIC)
Observation Date/Time	09/23/2015 2:00 PM
Observation end Date/Time	
Specimen Action Code	G
Relevant Clinical Information	
Relevant Clinical Information Original Text	
Observation Result Information	
Result Status	F
Results Report/Status Change - Date/Time	09/27/2015 11:20 AM
Results Copy To	
Name	Pafford
Identifier	10092000194
Results Handling	
Standard	
Observation Notes	

#### Timing/Quantity Information

Element	Data
Priority	
Start Date/time	
End Date/time	

#### Results Performing Laboratory

Element	Data
Laboratory Name	Century Hospital
Organization identifier	24D9871327
Address	2070 Test Park Los Angeles CA 90067
Director Name	Phil J. Knowsalot
Director identifier	5432178916

#### Lab results

Element		Data						
Test performed		Bacteria susceptibility						
Test Report date		09/27/2015 11:20						
Result Observation Name	Result	UOM	Range	Abnormal Flag	Status	Date/Time of Observation	Date/Time of Analysis	Notes

Ampicillin [Susceptibility] by Minimum inhibitory concentration (MIC)	< 0.06	ug/mL		S	F	09/23/2015 14:00	09/26/2015 11:00	
Gentamicin [Susceptibility] by Minimum inhibitory concentration (MIC)	0.05	ug/mL		S	F	09/23/2015 14:00	09/26/2015 11:00	
Ciprofloxacin [Susceptibility] by Minimum inhibitory concentration (MIC)	0.05	ug/mL		S	F	09/23/2015 14:00	09/26/2015 11:00	

### Order Observation

### Ordering Provider

Element	Data
Name	Nicholas Radon
Identifier number	5742200012

### Observation Details

Element	Data
Observation General Information	
Placer Order Number	
Filler Order Number	R-783274-7
Placer Group Number	GORD874211
Parent Universal Service Identifier	
Identifier	
Text	
Alt Identifier	
Alt Text	
Original Text	
Observation Details	
Universal Service Identifier	Bacterial susceptibility panel in Isolate by Minimum inhibitory concentration (MIC)
Observation Date/Time	09/23/2015 2:00 PM
Observation end Date/Time	
Specimen Action Code	G
Relevant Clinical Information	
Relevant Clinical Information Original Text	
Observation Result Information	
Result Status	C
Results Report/Status Change - Date/Time	09/27/2015 4:42 PM
Results Copy To	
Name	Pafford

Identifier	10092000194
Results Handling	
Standard	
Observation Notes	

#### Timing/Quantity Information

Element	Data
Priority	
Start Date/time	
End Date/time	

#### Results Performing Laboratory

Element	Data
Laboratory Name	Century Hospital
Organization identifier	24D9871327
Address	2070 Test Park Los Angeles CA 90067
Director Name	Phil J. Knowsalot
Director identifier	5432178916

#### Lab results

Element		Data						
Test performed		Bacteria susceptibility						
Test Report date		09/27/2015 16:42						
Result Observation Name	Result	UOM	Range	Abnormal Flag	Status	Date/Time of Observation	Date/Time of Analysis	Notes
Ampicillin [Susceptibility] by Minimum inhibitory concentration (MIC)	< 16	ug/mL		I	F	09/23/2015 14:00	09/26/2015 11:00	
Trimethoprim+Sulfamethoxazole [Susceptibility] by Minimum inhibitory concentration (MIC)	2 / 38	ug/mL		S	B	09/23/2015 14:00	09/27/2015 11:20	Due to the indeterminate amoxicillin test result, additional antibiotics were tested and appended to the previous report.
Ciprofloxacin [Susceptibility] by Minimum inhibitory concentration (MIC)	0.05	ug/mL		S	B	09/23/2015 14:00	09/27/2015 11:20	Due to the indeterminate amoxicillin test result, additional antibiotics were tested and appended to the previous report.



MSH|^~&|^2.16.840.1.113883.3.72.5.20^ISO|^2.16.840.1.113883.3.72.5.21^ISO||^2.16.840.1.113883.3.72.5.23^ISO|20150926140551||ORU^R01^ORU\_R01|LRI\_4.1\_4.1-

GU\_FRU|D|2.5.1||AL|AL||||LRI\_Common\_Component^^2.16.840.1.113883.9.16^ISO~LRI\_GU\_Component^^2.16.840.1.113883.9.12^ISO~LRI\_FRU\_Component^^2.16.840.1.113883.9.83^ISO

PID|1||PATID|234^^&2.16.840.1.113883.3.72.5.30.2&ISO^MR||Jones^William^A^^^^L||19610615|M||2106-3^White^HL70005|||||PATID|234^^&2.16.840.1.113883.3.72.5.30.2&ISO^AN

ORC|RE|ORD723222-4^^2.16.840.1.113883.3.72.5.24^ISO|R-783274-4^^2.16.840.1.113883.3.72.5.25^ISO|GORD874211^^2.16.840.1.113883.3.72.5.24^ISO|||||5742200012^Radon^Nicholas^^^^^^&2.16.840.1.113883.4.6&ISO^L^^NPI

OBR|1|ORD723222-4^^2.16.840.1.113883.3.72.5.24^ISO|R-783274-4^^2.16.840.1.113883.3.72.5.25^ISO|625-4^Bacteria identified in Stool by Culture^LN^3456543^CULTURE STOOL^99USL^2.52^^Stool Culture||201509231400|||||5742200012^Radon^Nicholas^^^^^^&2.16.840.1.113883.4.6&ISO^L^^NPI|||||20150926140551||F||10092000194^Hamlin^Pafford^^^^^^&2.16.840.1.113883.4.6&ISO^L^^NPI|||||MIC^Observation of type microbiology^HL70411^^^^2.5.1||CC^Copies Requested^HL70507

OBX|1|CWE|625-4^Bacteria identified in Stool by Culture^LN^^^^2.52^^Stool Culture|^1^1^Isl-1|103429008^Enterohemorrhagic Escherichia coli, serotype O157:H7^SCT^^^^201509USED^^Shiga toxin producing E. coli O157:H7 isolated||A||F||201509231400||||201509251930|||Century Hospital^^^^&2.16.840.1.113883.4.7&ISO^XX^^24D9871327|2070 Test Park^^Los Angeles^CA^90067^^B|5432178916^Knowsalot^Phil^J.^^^^&2.16.840.1.113883.4.6&ISO^L^^NPI||||RSLT

NTE|1||Susceptibility testing for E.coli is not performed, because antibiotics should not be used to treat this infection. There is no evidence that treatment with antibiotics is helpful, and taking antibiotics may increase the risk of hemolytic-uremic syndrome (HUS). Antidiarrheal agents like Imodium may also increase that risk. Non-specific supportive therapy, including hydration, is important.

OBX|2|CWE|625-4^Bacteria identified in Stool by Culture^LN^^^^2.52^^Stool Culture|^2^1^Isl-2|398567006^Salmonella I, group O:4^SCT^^^^201509USED^^Salmonella I, group O:4 isolated||A||F||201509231400||||201509251930|||Century Hospital^^^^&2.16.840.1.113883.4.7&ISO^XX^^24D9871327|2070 Test Park^^Los Angeles^CA^90067^^B|5432178916^Knowsalot^Phil^J.^^^^&2.16.840.1.113883.4.6&ISO^L^^NPI||||RSLT

NTE|1||Salmonella gastrointestinal infections usually resolve in 5-7 days and most do not require treatment other than oral fluids. Persons with severe diarrhea may require rehydration with intravenous fluids. Antibiotic therapy can prolong the duration of excretion of non-typhoidal Salmonella and is recommended only for patients with severe illness (e.g., those with severe diarrhea, high fever, bloodstream infection, or who need hospitalization) or those at risk of severe disease or complications, including young infants, older adults (over 65 years old) and immunocompromised persons. Antibiotic resistance is increasing among some Salmonella bacteria; therefore, susceptibility testing can help guide appropriate therapy. Choices for antibiotic therapy for severe infections include fluoroquinolones, third-generation cephalosporins, and ampicillin (for susceptible infections).

OBX|3|CWE|625-4^Bacteria identified in Stool by Culture^LN^^^^2.52^^Stool Culture|^3^1^Isl-3|85729005^Shigella flexneri^SCT^^^^201509USED^^Shigella flexneri isolated||A||F||201509231400||||201509251930|||Century Hospital^^^^&2.16.840.1.113883.4.7&ISO^XX^^24D9871327|2070 Test Park^^Los Angeles^CA^90067^^B|5432178916^Knowsalot^Phil^J.^^^^&2.16.840.1.113883.4.6&ISO^L^^NPI||||RSLT

NTE|1||Antibiotic treatment is recommended for patients with severe disease, bloody diarrhea, or compromised immune systems. Resistance to traditional first-line drugs like ampicillin and trimethoprim-sulfamethoxazole is common, and resistance to some other antibiotics is increasing. With this in mind, antibiotic susceptibility testing can help guide appropriate therapy. When an ampicillin- or trimethoprim-sulfamethoxazole-resistant strain is isolated, choices for therapy include fluoroquinolones, ceftriaxone, and azithromycin.

SPM|1|S-2015-66&&2.16.840.1.113883.3.72.5.24&ISO^S-9911-33&&2.16.840.1.113883.3.72.5.25&ISO||119339001^Stool specimen^SCT^^^^201509USED^^Stool|||||201509231400

ORC|RE|R-783274-6^^2.16.840.1.113883.3.72.5.25^ISO|GORD874211^^2.16.840.1.113883.3.72.5.24^ISO|||||5742200012^Radon^Nicholas^^^^^^&2.16.840.1.113883.4.6&ISO^L^^NPI

OBR|2||R-783274-6^^2.16.840.1.113883.3.72.5.25^ISO|50545-3^Bacterial susceptibility panel in Isolate by Minimum inhibitory concentration (MIC)^LN^^^^2.52^^Bacteria susceptibility||201509231400||G|||5742200012^Radon^Nicholas^^^^^^&2.16.840.1.113883.4.6&ISO^L^^NPI|||||20150927112054||F|625-4^Bacteria identified in Stool by Culture&LN&&&2.52&&Stool Culture^&2&1&Isl-2||10092000194^Hamlin^Pafford^^^^^^&2.16.840.1.113883.4.6&ISO^L^^NPI|ORD723222-4&&2.16.840.1.113883.3.72.5.24&ISO^R-783274-4&&2.16.840.1.113883.3.72.5.25&ISO|||||MIC^Observation of type microbiology^HL70411^^^^2.5.1||CC^Copies Requested^HL70507

OBX|1|SN|28-1^Ampicillin [Susceptibility] by Minimum inhibitory concentration (MIC)^LN^^^^2.52|^1^1^IsIt-2|<^0.06|ug/mL  
^^UCUM^^^^1.9||S||F|||201509231400|||201509261100|||Century Hospital^^^^&2.16.840.1.113883.4.7&ISO^XX^^24D9871327|  
2070 Test Park^^Los Angeles^CA^90067^^B|5432178916^Knowsallot^Phil^J.^^^^&2.16.840.1.113883.4.6&ISO^L^^^NPI|||RSLT|SUR

OBX|2|SN|267-5^Gentamicin [Susceptibility] by Minimum inhibitory concentration (MIC)^LN^^^^2.52|^1^1^IsIt-2|^0.05|ug/mL  
^^UCUM^^^^1.9||S||F|||201509231400|||201509261100|||Century Hospital^^^^&2.16.840.1.113883.4.7&ISO^XX^^24D9871327|  
2070 Test Park^^Los Angeles^CA^90067^^B|5432178916^Knowsallot^Phil^J.^^^^&2.16.840.1.113883.4.6&ISO^L^^^NPI|||RSLT|SUR

OBX|3|SN|185-9^Ciprofloxacin [Susceptibility] by Minimum inhibitory concentration (MIC)^LN^^^^2.52|^1^1^IsIt-2|^0.05|ug  
/mL^^UCUM^^^^1.9||S||F|||201509231400|||201509261100|||Century Hospital^^^^&2.16.840.1.113883.4.7&ISO^XX^^24D98713  
27|2070 Test Park^^Los Angeles^CA^90067^^B|5432178916^Knowsallot^Phil^J.^^^^&2.16.840.1.113883.4.6&ISO^L^^^NPI|||RSLT|S  
UR

ORC|RE||R-783274-7^^2.16.840.1.113883.3.72.5.25^ISO|GORD874211^^2.16.840.1.113883.3.72.5.24^ISO|||||5742200012^Radon  
^Nicholas^^^^&2.16.840.1.113883.4.6&ISO^L^^^NPI

OBR|3||R-783274-7^^2.16.840.1.113883.3.72.5.25^ISO|50545-3^Bacterial susceptibility panel in Isolate by Minimum inhibit  
ory concentration (MIC)^LN^^^^2.52^^Bacteria susceptibility|||201509231400|||G||||5742200012^Radon^Nicholas^^^^&2.16  
.840.1.113883.4.6&ISO^L^^^NPI|||||20150927164251||C|625-4&Bacteria identified in Stool by Culture&LN&&&2.52&&Stool Cu  
lture^&3&1&IsIt-3||10092000194^Hamlin^Pafford^^^^&2.16.840.1.113883.4.6&ISO^L^^^NPI|ORD723222-  
4&&2.16.840.1.113883.3.7  
2.5.24&ISO^R-783274-4&&2.16.840.1.113883.3.72.5.25&ISO|||||MIC^Observation of type microbiology^HL70411^^^^  
2.5.1||CC^Copies Requested^HL70507

OBX|1|SN|28-1^Ampicillin [Susceptibility] by Minimum inhibitory concentration (MIC)^LN^^^^2.52|^1^1^IsIt-3|<^16|ug/mL^^  
UCUM^^^^1.9||I||F|||201509231400|||201509261100|||Century Hospital^^^^&2.16.840.1.113883.4.7&ISO^XX^^24D9871327|20  
70 Test Park^^Los Angeles^CA^90067^^B|5432178916^Knowsallot^Phil^J.^^^^&2.16.840.1.113883.4.6&ISO^L^^^NPI|||RSLT|SUR

OBX|2|SN|516-5^Trimethoprim+Sulfamethoxazole [Susceptibility] by Minimum inhibitory concentration (MIC)^LN^^^^2.52|^1^1  
^IsIt-3|^2^/38|ug/mL^^UCUM^^^^1.9||S||B|||201509231400|||201509271120|||Century Hospital^^^^&2.16.840.1.113883.4.7  
&ISO^XX^^24D9871327|2070 Test Park^^Los  
Angeles^CA^90067^^B|5432178916^Knowsallot^Phil^J.^^^^&2.16.840.1.113883.4.6&ISO  
^L^^^NPI|||RSLT|SUR

NTE|1||Due to the indeterminate amoxicillin test result, additional antibiotics were tested and appended to the previous report.

OBX|3|SN|185-9^Ciprofloxacin [Susceptibility] by Minimum inhibitory concentration (MIC)^LN^^^^2.52|^1^1^IsIt-3|^0.05|ug  
/mL^^UCUM^^^^1.9||S||B|||201509231400|||201509271120|||Century Hospital^^^^&2.16.840.1.113883.4.7&ISO^XX^^24D98713  
27|2070 Test Park^^Los Angeles^CA^90067^^B|5432178916^Knowsallot^Phil^J.^^^^&2.16.840.1.113883.4.6&ISO^L^^^NPI|||RSLT|S  
UR

NTE|1||Due to the indeterminate amoxicillin test result, additional antibiotics were tested and appended to the previous report.