Description

Mr. William A. Jones is a 51 year old white male who presented with pain, swelling, and redness in several joints including his elbows, wrists, knees, and hips. In order to screen for rheumatoid diseases, Dr. Nicholas Radon ordered an erythrocyte sedimentation rate (ESR or Sed Rate) blood test using the Good Health Clinic's EHR system, specifying that the laboratory's patient service center was to draw the blood specimen from the patient later that same day. Dr Radon also requests that the results be CC'd to Dr Pafford M. Hamlin Sr, a specialist to which he is referring Mr. Jones as well as Dr. Daniel Davison who is Mr. Jones' regular doctor.

The laboratory performs the requested test and the final result was 10 mm/hr, a normal value with the normal reference range being 0-17 mm/hr. The final result report for this lab test was generated by the LIS and transmitted to the patient's record in the Good health Clinic's EHR system used in Dr. Radon's office practice on the same day.

PreCondition No Pre-condition.

-PostCondition

No Post Condition.

TestObjectives

- Maximally populated message.
- Demonstrate ability to import and incorporate all supported relevant (R,RE,C(a/b)) data elements for the test case, including multiple
 occurrences:
- Demonstrate the capability to support Timing Qty Group and Timing/Quantity Segment (TQ1).
- Demonstrate the capability to support Notes and Comments Segments (NTE) that follow the Observation Request Segment (OBR), including multiple occurrences.
- Maximally populated MSH Segment Demonstrate the capability of the EHR to support all supported data elements including repeating elements in the MSH segment.
- Maximally populated PID Segment Demonstrate the capability of the EHR to support all supported data elements including repeating elements in the PID segment.
- Maximally populated NTE segment Demonstrate the capability of the EHR to support all supported data elements including support of formatted text.
- Maximally populated ORC Segment Demonstrate the capability of the EHR to support all supported data elements including repeating elements in the ORC segment.
- Maximally populated OBR Segment Demonstrate the capability of the EHR to support all supported data elements including repeating elements in the OBR segment.
- Maximally populated SPM Segment Demonstrate the capability of the EHR to support all supported data elements in the SPM segment.
- Maximally populated TQ1 Segment Demonstrate the capability of the EHR to support all supported data elements including repeating elements in the TQ1 segment .
- Specific test points include the use of the NM data type for OBX.2, the value of "10" for results (OBX.5), the value of "N" for abnormal flags (OBX.8), and the observation results status is "F" (OBX.11).

Notes to Testers

Note that data element SPM-24 is provided in the test data. This element is optional for the sender and required-but maybe empty for the receiver. This element is not subject to conformance testing for this test case (sender).

-Test Case Information-

LRI_1.0_1.1-GU - Maximally Populated SED Rate message - Final Results

Test Case ID LRI_1.0_1.1-GU

-MSH-

Location	Data Element	Data	Categorization
MSH.1	Field Separator		IG Fixed Data
MSH.2	Encoding Characters	~\&#</td><td>IG Fixed Data</td></tr><tr><td>MSH.3</td><td>Sending Application</td><td></td><td></td></tr><tr><td>MSH.3.1</td><td>Namespace ID</td><td>NIST Test Lab APP</td><td>Configurable Data</td></tr><tr><td>MSH.3.2</td><td>Universal ID</td><td>2.16.840.1.113883.3.72.5.20</td><td>Configurable Data</td></tr><tr><td>MSH.3.3</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr><tr><td>MSH.4</td><td>Sending Facility</td><td></td><td></td></tr><tr><td>MSH.4.1</td><td>Namespace ID</td><td>NIST Lab Facility</td><td>Configurable Data</td></tr><tr><td>MSH.4.2</td><td>Universal ID</td><td>2.16.840.1.113883.3.72.5.21</td><td>Configurable Data</td></tr><tr><td>MSH.4.3</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr><tr><td>MSH.6</td><td>Receiving Facility</td><td></td><td></td></tr><tr><td>MSH.6.1</td><td>Namespace ID</td><td>NIST EHR Facility</td><td>Configurable Data</td></tr><tr><td>MSH.6.2</td><td>Universal ID</td><td>2.16.840.1.113883.3.72.5.23</td><td>Configurable Data</td></tr><tr><td>MSH.6.3</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr><tr><td>MSH.7</td><td>Date/Time Of Message</td><td></td><td></td></tr><tr><td>MSH.7.1</td><td>Time</td><td>20150926140551</td><td>System Generated</td></tr><tr><td>MSH.9</td><td>Message Type</td><td></td><td></td></tr><tr><td>MSH.9.1</td><td>Message Code</td><td>ORU</td><td>IG Fixed Data</td></tr><tr><td>MSH.9.2</td><td>Event Type</td><td>R01</td><td>IG Fixed Data</td></tr><tr><td>MSH.9.3</td><td>Message Structure</td><td>ORU_R01</td><td>IG Fixed Data</td></tr><tr><td>MSH.10</td><td>Message Control ID</td><td>LRI_1.0_1.1-GU</td><td>System Generated</td></tr><tr><td>MSH.11</td><td>Processing ID</td><td></td><td></td></tr><tr><td>MSH.11.1</td><td>Processing ID</td><td>D</td><td>Changeable Data</td></tr><tr><td>MSH.12</td><td>VersionID</td><td></td><td></td></tr><tr><td>MSH.12.1</td><td>Version ID</td><td>2.5.1</td><td>IG Fixed Data</td></tr><tr><td>MSH.15</td><td>Accept Acknowledgment Type</td><td>AL</td><td>IG Fixed Data</td></tr><tr><td>MSH.16</td><td>Application Acknowledgment Type</td><td>AL</td><td>IG Fixed Data</td></tr><tr><td>MSH.21</td><td>Message Profile Identifier</td><td></td><td></td></tr><tr><td>MSH.21.1</td><td>Entity Identifier</td><td>LRI Common Component</td><td>IG Fixed Data</td></tr><tr><td>MSH.21.2</td><td>Namespace ID</td><td>Profile Component</td><td>Changeable Data</td></tr><tr><td>MSH.21.3</td><td>Universal ID</td><td>2.16.840.1.113883.9.16</td><td>IG Fixed Data</td></tr><tr><td>MSH.21.4</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr><tr><td>MSH.21[2]</td><td>Message Profile Identifier</td><td></td><td></td></tr><tr><td>MSH.21[2].1</td><td>Entity Identifier</td><td>LRI_GU_Component</td><td>IG Fixed Data</td></tr><tr><td>MSH.21[2].2</td><td>Namespace ID</td><td>Profile Component</td><td>Changeable Data</td></tr><tr><td>MSH.21[2].3</td><td>Universal ID</td><td>2.16.840.1.113883.9.12</td><td>IG Fixed Data</td></tr><tr><td>MSH.21[2].4</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr><tr><td>MSH.21[3]</td><td>Message Profile Identifier</td><td></td><td></td></tr><tr><td>MSH.21[3].1</td><td>Entity Identifier</td><td>LRI_FRU_Component</td><td>IG Fixed Data</td></tr><tr><td>MSH.21[3].2</td><td>Namespace ID</td><td>Profile Component</td><td>Changeable Data</td></tr><tr><td>MSH.21[3].3</td><td>Universal ID</td><td>2.16.840.1.113883.9.83</td><td>IG Fixed Data</td></tr><tr><td>MSH.21[3].4</td><td>Universal ID Type</td><td>ISO</td><td>IG Fixed Data</td></tr></tbody></table>	

-PID-

Location	Data Element	Data	Categorization
PID.1	Set ID - PID	1	IG Fixed Data
PID.3	Patient Identifier List		
PID.3.1	ID Number	PATID1234	Configurable Data
PID.3.4	Assigning Authority		

P Doeation	Namespa pata Element	NIST MPI Data	Changealle gorization
PID.3.4.2	Universal ID	2.16.840.1.113883.3.72.5.30.2	Configurable Data
PID.3.4.3	Universal ID Type	ISO	IG Fixed Data
PID.3.5	Identifier Type Code	MR	Configurable Data
PID.5	Patient Name		
PID.5.1	Family Name		
PID.5.1.1	Surname	Jones	Changeable Data
PID.5.2	Given Name	William	Changeable Data
PID.5.3	Second and Further Given Names or Initials Thereof	A	Changeable Data
PID.5.7	Name Type Code	L	Changeable Data
PID.7	Date/Time of Birth		
PID.7.1	Time	19610615	Changeable Data
PID.8	Administrative Sex	M	Changeable Data
PID.10	Race		
PID.10.1	Identifier	2106-3	Changeable Data
PID.10.2	Text	White	Changeable Data
PID.10.3	Name of the Coding System	HL70005	Changeable Data
PID.10.9	Original Text	White	Changeable Data
PID.10[2]	Race		
PID.10[2].1	Identifier	1002-5	Changeable Data
PID.10[2].2	Text	American Indian or Alaska Native	Changeable Data
PID.10[2].3	Name of the Coding System	HL70005	Changeable Data
PID.10[2].9	Original Text	American Indian	Changeable Data
PID.18	Patient Account Number		
PID.18.1	ID Number	PATACC1234	Configurable Data
PID.18.4	Assigning Authority		
PID.18.4.1	Namespace ID	NIST MPI	Changeable Data
PID.18.4.2	Universal ID	2.16.840.1.113883.3.72.5.30.2	Configurable Data
PID.18.4.3	Universal ID Type	ISO	IG Fixed Data
PID.18.5	Identifier Type Code	AN	Configurable Data

-ORC-

Location	Data Element	Data	Categorization
ORC.1	Order Control	RE	Test Case Fixed Data
ORC.2	Placer Order Number		
ORC.2.1	Entity Identifier	ORD723222	Changeable Data
ORC.2.2	Namespace ID	NIST EHR	Changeable Data
ORC.2.3	Universal ID	2.16.840.1.113883.3.72.5.24	Changeable Data
ORC.2.4	Universal ID Type	ISO	IG Fixed Data
ORC.3	Filler Order Number		
ORC.3.1	Entity Identifier	R-783274	Changeable Data
ORC.3.2	Namespace ID	NIST Lab Filler	Changeable Data
ORC.3.3	Universal ID	2.16.840.1.113883.3.72.5.25	Changeable Data
ORC.3.4	Universal ID Type	ISO	IG Fixed Data
ORC.4	Placer Group Number		
ORC.4.1	Entity Identifier	GORD874211	Changeable Data
ORC.4.2	Namespace ID	NIST EHR	Changeable Data
ORC.4.3	Universal ID	2.16.840.1.113883.3.72.5.24	Changeable Data
ORC.4.4	Universal ID Type	ISO	IG Fixed Data
ORC.12	Ordering Provider		
ORC.12.1	ID Number	5742200012	Changeable Data
ORC.12.2	Family Name		
ORC.12.2.1	Surname	Radon	Changeable Data
ORC.12.3	Given Name	Nicholas	Changeable Data
ORC.12.4	Second and Further Given Names or Initials Thereof	M	Changeable Data

OR Location	Suffix (e.g., patar Element	^{JR} Data	Changeallegorization
ORC.12.6	Prefix (e.g., DR)	DR	Changeable Data
ORC.12.9	Assigning Authority		
ORC.12.9.1	Namespace ID	NPI	Changeable Data
ORC.12.9.2	Universal ID	2.16.840.1.113883.4.6	Changeable Data
ORC.12.9.3	Universal ID Type	ISO	IG Fixed Data
ORC.12.10	Name Type Code	L	Changeable Data
ORC.12.13	Identifier Type Code	NPI	Changeable Data

-OBR-

Location	Data Element	Data	Categorization
OBR.1	Set ID - OBR	1	IG Fixed Data
OBR.2	Placer Order Number		
OBR.2.1	Entity Identifier	ORD723222	Changeable Data
OBR.2.2	Namespace ID	NIST EHR	Changeable Data
OBR.2.3	Universal ID	2.16.840.1.113883.3.72.5.24	Changeable Data
OBR.2.4	Universal ID Type	ISO	Changeable Data
OBR.3	Filler Order Number		
OBR.3.1	Entity Identifier	R-783274	Changeable Data
OBR.3.2	Namespace ID	NIST Lab Filler	Changeable Data
OBR.3.3	Universal ID	2.16.840.1.113883.3.72.5.25	Changeable Data
OBR.3.4	Universal ID Type	ISO	IG Fixed Data
OBR.4	Universal Service Identifier		
OBR.4.1	Identifier	30341-2	Test Case Fixed Data
OBR.4.2	Text	Erythrocyte sedimentation rate	Test Case Fixed Data
OBR.4.3	Name of Coding System	LN	Test Case Fixed Data
OBR.4.4	Alternate Identifier	500	Changeable Data
OBR.4.5	Alternate Text	Erythrocyte sedimentation rate	Changeable Data
OBR.4.6	Name of Alternate Coding System	99USL	Changeable Data
OBR.4.7	Coding System Version	2.52	Changeable Data
OBR.4.9	Original Text	Erythrocyte sedimentation rate	Changeable Data
OBR.7	Observation Date/Time		
OBR.7.1	Time	201509251400	Changeable Data
OBR.11	Specimen Action Code	L	Changeable Data
OBR.16	Ordering Provider		
OBR.16.1	ID Number	5742200012	Changeable Data
OBR.16.2	Family Name		
OBR.16.2.1	Surname	Radon	Changeable Data
OBR.16.3	Given Name	Nicholas	Changeable Data
OBR.16.4	Second and Further Given Names or Initials Thereof	M	Changeable Data
OBR.16.5	Suffix (e.g., JR or III)	JR	Changeable Data
OBR.16.6	Prefix (e.g., DR)	DR	Changeable Data
OBR.16.9	Assigning Authority		
OBR.16.9.1	Namespace ID	NPI	Changeable Data
OBR.16.9.2	Universal ID	2.16.840.1.113883.4.6	Changeable Data
OBR.16.9.3	Universal ID Type	ISO	IG Fixed Data
OBR.16.10	Name Type Code	L	Changeable Data
OBR.16.13	Identifier Type Code	NPI	Changeable Data
OBR.22	Results Rpt/Status Chng - Date/Time		
OBR.22.1	Time	20150926140551	Changeable Data
OBR.25	Result Status	F	Test Case Fixed Data
OBR.28	Result Copies To		
OBR.28.1	ID Number	10092000194	Changeable Data
OBR.28.2	Family Name		S
	Surname	Hamlin	Changeable Data
OBR.28.2.1	Surname	1141111111	

OBROSSION	Second and Further Given Names or Initials Thereof	M Data	Chang Cattle goaization
OBR.28.5	Suffix (e.g., JR or III)	Sr.	Changeable Data
OBR.28.6	Prefix (e.g., DR)	Dr.	Changeable Data
OBR.28.9	Assigning Authority		
OBR.28.9.1	Namespace ID	NPI	Changeable Data
OBR.28.9.2	Universal ID	2.16.840.1.113883.4.6	Changeable Data
OBR.28.9.3	Universal ID Type	ISO	Changeable Data
OBR.28.10	Name Type Code	L	Changeable Data
OBR.28.13	Identifier Type Code	NPI	Changeable Data
OBR.28[2]	Result Copies To		
OBR.28[2].1	ID Number	2129416824	Changeable Data
OBR.28[2].2	Family Name		
OBR.28[2].2.1	Surname	Davison	Changeable Data
OBR.28[2].3	Given Name	Daniel	Changeable Data
OBR.28[2].4	Second and Further Given Names or Initials Thereof	D	Changeable Data
OBR.28[2].5	Suffix (e.g., JR or III)	III	Changeable Data
OBR.28[2].6	Prefix (e.g., DR)	Dr.	Changeable Data
OBR.28[2].9	Assigning Authority		
OBR.28[2].9.1	Namespace ID	NPI	Changeable Data
OBR.28[2].9.2	Universal ID	2.16.840.1.113883.4.6	Changeable Data
OBR.28[2].9.3	Universal ID Type	ISO	IG Fixed Data
OBR.28[2].10	Name Type Code	L	Changeable Data
OBR.28[2].13	Identifier Type Code	NPI	Changeable Data
OBR.49	Results Handling		
OBR.49.1	Identifier	CC	Test Case Fixed Data
OBR.49.2	Text	Copies Requested	Test Case Fixed Data
OBR.49.3	Name of the Coding System	HL70507	IG Fixed Data

-NTE-

Location	Data Element	Data	Categorization
NTE.1	Set ID - NTE	1	IG Fixed Data
NTE.3	Comment	Patient is extremely anxious about needles used for drawing blood.\br\If patient is overly frightened, nervous, or anxious please reschedule blood draw.	Changeable Data

-NTE-

Location	Data Element	Data	Categorization
NTE.1	Set ID - NTE	2	IG Fixed Data
NTE.3	Comment	Patient is allergic to latex	Changeable Data

-TQ1-

· · · · · · · · · · · · · · · · · · ·			
Location	Data Element	Data	Categorization
TQ1.1	Set ID - TQ1	1	IG Fixed Data
TQ1.7	Start date/time		
TQ1.7.1	Time	201509251400	Changeable Data
TQ1.8	End date/time		
TQ1.8.1	Time	201509261400	Changeable Data
TQ1.9	Priority		
TQ1.9.1	Identifier	R	Changeable Data
TQ1.9.2	Text	Routine	Changeable Data
TQ1.9.3	Name of the Coding System	HL70485	Changeable Data
TQ1.9.9	Original Text	Routine	Changeable Data

-OBX-

Location	Data Element	Data	Categorization

OBX _{Location}	Set ID - OB Data Element	Data Data	IG Fixe alegorization
OBX.2	Value Type	NM	Test Case Fixed Data
OBX.3	Observation Identifier		
OBX.3.1	Identifier	30341-2	Test Case Fixed Data
OBX.3.2	Text	Erythrocyte sedimentation rate	Test Case Fixed Data
OBX.3.3	Name of the Coding System	LN	Test Case Fixed Data
OBX.3.4	Alternate Identifier	815117	Changeable Data
OBX.3.5	Alternate Text	ESR	Changeable Data
OBX.3.6	Name of Alternate Coding System	99USL	Changeable Data
OBX.3.7	Coding System Version	2.52	Changeable Data
OBX.3.9	Original Text	Erythrocyte sedimentation rate	Changeable Data
OBX.4	Observation Sub-ID		
OBX.4.2	Group	1	Test Case Fixed Data
OBX.4.3	Sequence	1	Test Case Fixed Data
OBX.4.4	Identifier	1	Test Case Fixed Data
OBX.5	Observation Value	10	Test Case Fixed Data
OBX.6	Units		
OBX.6.1	Identifier	mm/h	Changeable Data
OBX.6.2	Text	millimeter per hour	Changeable Data
OBX.6.3	Name of the Coding System	UCUM	Changeable Data
OBX.6.7	Coding System Version	1.9	Changeable Data
OBX.7	Reference Range	0 to 17	Changeable Data
OBX.8	Abnormal Flags	N	Test Case Fixed Data
OBX.11	Observation Result Status	F	Test Case Fixed Data
OBX.11	Date/Time of the Observation		Cast I area Batta
OBX.14.1	Time	201509251400	Changeable Data
OBX.19	Date/Time of the Analysis	201309231100	Changeaole Data
OBX.19.1	Time	20150926130550	Changeable Data
OBX.13.1	Performing Organization Name	20130720130330	Changeaoic Data
OBX.23.1	Organization Name	Century Hospital	Changeable Data
OBX.23.6	Assigning Authority	Century Prospitar	Changeaoic Data
OBX.23.6.1	Namespace ID	CLIA	Changachla Data
OBX.23.6.1 OBX.23.6.2	Universal ID		Changeable Data
		2.16.840.1.113883.4.7 ISO	Changeable Data IG Fixed Data
OBX.23.6.3	Universal ID Type		
OBX.23.7	Identifier Type Code	XX	Changeable Data
OBX.23.10	Organization Identifier	24D9871327	Changeable Data
OBX.24	Performing Organization Address		
OBX.24.1	Street Address	2050 F P	(d) 11 =
OBX.24.1.1	Street or Mailing Address	2070 Test Park	Changeable Data
OBX.24.3	City	Los Angeles	Changeable Data
OBX.24.4	State or Province	CA	Changeable Data
OBX.24.5	Zip or Postal Code	90067	Changeable Data
OBX.24.6	Country	USA	Changeable Data
OBX.24.7	Address Type	В	Changeable Data
OBX.24.9	County/Parish Code	06037	Changeable Data
OBX.25	Performing Organization Medical Director		
OBX.25.1	ID Number	5432178916	Changeable Data
OBX.25.2	Family Name		
OBX.25.2.1	Surname	Knowsalot	Changeable Data
OBX.25.3	Given Name	Phil	Changeable Data
OBX.25.4	Second and Further Given Names or Initials Thereof	J.	Changeable Data
0071077	Suffix (e.g., JR or III)	III	Changeable Data
OBX.25.5			-
OBX.25.5	Prefix (e.g. DR)	Dr	Changeable Data
OBX.25.6 OBX.25.9	Prefix (e.g., DR) Assigning Authority	Dr.	Changeable Data

OBX.25.9.3	Universal Data Element Universal ID Type	2.16.840.1.113883. 15ata ISO	Chang eallegorization IG Fixed Data		
OBX.25.10	Name Type Code	L	Changeable Data		
OBX.25.13	Identifier Type Code	NPI	Changeable Data		
OBX.29	Observation Type	RSLT	Test Case Fixed Data		
OBX.30	Observation SubType	UNSP	Test Case Fixed Data		

-SPM-

Location	Data Element	Data	Categorization	
SPM.1	Set ID - SPM	1	IG Fixed Data	
SPM.2	Specimen ID			
SPM.2.1				
SPM.2.1.1		S-2015-66	Configurable Data	
SPM.2.1.2		GoodHealthC_EHR	Configurable Data	
SPM.2.1.3		2.16.840.1.113883.3.72.5.24	Configurable Data	
SPM.2.1.4		ISO	Changeable Data	
SPM.2.2				
SPM.2.2.1		S-9911-33	Changeable Data	
SPM.2.2.2		NIST Lab Filler	Changeable Data	
SPM.2.2.3		2.16.840.1.113883.3.72.5.25	Changeable Data	
SPM.2.2.4		ISO	IG Fixed Data	
SPM.4	Specimen Type			
SPM.4.1	Identifier	119297000	Changeable Data	
SPM.4.2	Text	BLD	Changeable Data	
SPM.4.3	Name of the Coding System	SCT	Changeable Data	
SPM.4.4	Alternate Identifier	BldSpc	Changeable Data	
SPM.4.5	Alternate Text	Blood	Changeable Data	
SPM.4.6	Name of Alternate Coding System	99USL	Changeable Data	
SPM.4.7	Coding System Version	201509USEd	Changeable Data	
SPM.4.9	Original Text	Blood Specimen	Changeable Data	
SPM.17	Specimen Collection Date/Time			
SPM.17.1	Range Start Date/Time			
SPM.17.1.1	Time	201509251400	Changeable Data	
SPM.24	Specimen Condition			
SPM.24.1	Identifier	COOL	Changeable Data	
SPM.24.2	Text	Cool	Changeable Data	
SPM.24.3	Name of the Coding System	HL70493	Changeable Data	
SPM.24.4	Alternate Identifier	CL	Changeable Data	
SPM.24.5	Alternate Text	Cool	Changeable Data	
SPM.24.6	Name of Alternate Coding System	99USL	Changeable Data	
SPM.24.9	Original Text	Cool	Changeable Data	

	Fatient information	
Element	Data	
Name	William A Jones	
Date/Time of Birth	06/15/1961	
Administrative Sex	Male	
Race	White	
Alt Race		
Race[2]	American Indian or Alaska Native	

Order Observation

Alt Race[2]

Ordering Provider-

Element	Data
Name	DR Nicholas M Radon JR
Identifier number	5742200012

Observation Details

Element	Data			
Observation General Information				
Placer Order Number	ORD723222			
Filler Order Number	R-783274			
Placer Group Number	GORD874211			
Parent Universal Service Identifier				
Identifier				
Text				
Alt Identifier				
Alt Text				
Original Text				
Observation Details				
Universal Service Identifier	Erythrocyte sedimentation rate			
Observation Date/Time	09/25/2015 2:00 PM			
Observation end Date/Time				
Specimen Action Code	L			
Relevant Clinical Information				
Relevant Clinical Information Original Text				
Observation Result Information				
Result Status	F			
Results Report/Status Change - Date/Time	09/26/2015 2:05 PM			
Results Report/Status Change - Date/Time	09/20/2013 2.03 FIVI			
Results Copy To	'			
Name	Dr. Pafford M Sr.			
Identifier	10092000194			
Results Copy To[2]				
Name	Dr. Daniel D III			
Identifier	2129416824			
Results Handling				

Standard	
Observation Notes	
Notes and comments	Patient is extremely anxious about needles used for drawing blood.\br\lf patient is overly frightened, nervous, or anxious please reschedule blood draw.
Notes and comments	Patient is allergic to latex

Timing/Quantity Information-

Element	Data
Priority	Routine
Start Date/time	09/25/2015 2:00 PM
End Date/time	09/26/2015 2:00 PM

Results Performing Laboratory

Element	Data		
Laboratory Name	Century Hospital		
Organization identifier	24D9871327		
Address	2070 Test Park Los Angeles CA 90067 USA		
Director Name	Dr. Phil J. Knowsalot III		
Director identifier	5432178916		

Specimen Information-

Element	Data			
Specimen Type	BLD			
Alt Specimen Type	Blood			
Specimen Original Text	Blood Specimen	Blood Specimen		
Start date/time	201509251400			
Specimen Condition	Cool			
Alt Specimen Condition	Cool			
Condition Original Text	Cool			

Lab results

Element			Data					
Test performed			Erythrocyte sedimentation rate					
Test Report date			09/26/2015 14:05					
Result Observation Name	Result	UOM	Range	Abnormal Flag	Status	Date/Time of Observation	Date/Time of Analysis	Notes
Erythrocyte sedimentation rate	10	millimeter per hour	0 to 17	N	F	09/25/2015 14:00	09/26/2015	

 $MSH/^{\width=0.113883.3.72.5.20^{\sc ility}^2.16.840.1.113883.3.72.5.20^{\sc ility}^2.16.840.1.113883.3.72.5.21^{\sc ility}^2.16.840.1.113883.3.72.5.23^{\sc ility}^2.16.840.1.113883.9.12^{\sc ility}^2.16.840.1$

 $PID|1||PATID1234^{^{N}}IST\ MPI\&2.16.840.1.113883.3.72.5.30.2\&ISO^{MR}||Jones^{William^{A^{^{N}}}L}||19610615|M||2106-3^{White^{HL7000}}5^{^{^{N}}}White^{1002-5^{American Indian}}OS^{^{N}}American Indian|||||||||PATACC1234^{^{N}}IST\ MPI\&2.16.840.\\ 1.113883.3.72.5.30.2\&ISO^{AN}$

 $ORC|RE|ORD723222^{NIST}EHR^{2}.16.840.1.113883.3.72.5.24^{ISO}|R-783274^{NIST}Lab\ Filler^{2}.16.840.1.113883.3.72.5.25^{ISO}|GORD-874211^{NIST}EHR^{2}.16.840.1.113883.3.72.5.24^{ISO}|||||||5742200012^{Radon^{Nicholas^{M^{P}}}}PR^{^{N^{N}}}PI\&2.16.840.1.113883.4.6\&ISO^{L^{N^{N}}}PI\&2.16.840.1.113883.4.6\&ISO^{N^{N}}PI\&2.16.840.1.113883.4.6&ISO^{N^{N}}PI\&2.16.840.1.11389.4&ISO^{N^{N}}PI^$

 $OBR|1|ORD723222^{N}IST\ EHR^{2}.16.840.1.113883.3.72.5.24^{I}SO|R-783274^{N}IST\ Lab\ Filler^{2}.16.840.1.113883.3.72.5.25^{I}SO|30341-2^{E}rythrocyte\ sedimentation\ rate^{L}N^{5}00^{E}rythrocyte\ sedimentation\ rate^{9}USL^{2}.52^{E}rythrocyte\ sedimentation\ rate||201509251400||||L||||5742200012^{R}adon^{N}icholas^{M}JR^{D}R^{^{N}P}L^{2}.16.840.1.113883.4.6&ISO^{L^{^{N}P}}|||||20150926140551|||F|||10092000194^{H}amlin^{P}afford^{M}Sr.^{D}r.^{^{N}P}L^{2}.16.840.1.113883.4.6&ISO^{L^{^{N}P}}L^{2}129416824^{D}avison^{D}aniel^{D}III^{D}r.^{^{N}P}L^{2}.16.840.1.113883.4.6&ISO^{L^{^{N}P}}L^{2}129416824^{D}avison^{D}aniel^{D}III^{D}r.^{^{N}P}L^{2}.16.840.1.113883.4.6&ISO^{L}C^{^{N}P}L^{2}.16.840.$

NTE|1||Patient is extremely anxious about needles used for drawing blood.\.br\If patient is overly frightened, nervous, or anxious please reschedule blood draw.

NTE|2||Patient is allergic to latex

TO1|1||||||201509251400|201509261400|R^Routine^HL70485^^^^^Routine

 $OBX[1]NM[30341-2^{Erythrocyte}] sedimentation rate^{LN^815117^{ESR^99USL^2.52^{Erythrocyte}} sedimentation rate]^{1^1^1|10}mm^{h^1} $$ ^millimeter per hour^{UCUM^{^^1.9}0 to 17}N[||F|||201509251400|||||20150926130550||||Century Hospital^{^^^}CLIA&2.16.840.1 .113883.4.7&ISO^{XX^{^24D9871327}2070} Test Park^{Los} $$ Angeles^{CA^90067^{USA^B^06037}5432178916^{Knowsalot^Phil^J.^{III^Dr.^^}} $$ ^NPI&2.16.840.1.113883.4.6&ISO^{L^{^NPI}}RSLTUNSP$

 $SPM|1|S-2015-66\&GoodHealthC_EHR\&2.16.840.1.113883.3.72.5.24\&ISO^S-9911-33\&NIST\ Lab\ Filler\&2.16.840.1.113883.3.72.5.25\&I\ SO||119297000^BLD^SCT^BldSpc^Blood^99USL^201509USEd^Blood\ Specimen|||||||||||||||201509251400|||||||COOL^Cool^HL70493^CL^Cool^99USL^^^Cool$